
CAS Spectrometer Driver

SDK

User Manual

4.10.2

1 Introduction	1
2 To-Do List	3
3 What has changed	5
4 Tasks	7
4.1 Interfaces Types and Options	7
4.1.1 Device discovery	7
4.1.2 Letting the user select interface type and option	7
4.1.3 Solutions without a UI	8
4.2 Device Handles (a.k.a. CASIDs) and Interfaces Types	8
4.3 Thread Safety	8
4.3.1 Global settings	9
4.4 Error Handling	9
4.5 Device Parameter (dpid...)	9
4.6 Measurement Parameter (mpid...)	10
4.7 Configuration and Initialization	10
4.8 Working with Calibration files on the device	11
4.8.1 Checking for updated calibrations	11
4.9 Density Filter	11
4.10 Dark Current	12
4.10.1 Dark current array	13
4.11 Performing a Measurement	13
4.12 Getting the Spectrum and Results	14
4.12.1 ColorMetric results	14
4.13 Synchronization	15
4.14 Trigger handshake sequence	16
4.14.1 For this sequence to work	16
4.14.2 Several measurements per DUT	17
4.15 Trigger CAS sequence	17
4.15.1 For this sequence to work	18
4.15.2 Several measurements per DUT	18
4.16 Software trigger sequence	19
4.17 Triggered Measurements	20
4.18 AutoRange measurements	21
4.18.1 AutoRange parameter	21
4.18.2 AutoRange for similar intensities	22
4.18.3 AutoRange for various intensities	22
4.19 Working with Parameter Sets	23
4.20 Using Serial Numbers	24
4.21 MultiTrack measurements	24
4.22 CCD Temperature Monitoring	25

4.23 Transmission measurement	26
5 Namespace Index	29
5.1 Packages	29
6 Class Index	31
6.1 Class List	31
7 Namespace Documentation	33
7.1 InstrumentSystems Namespace Reference	33
7.2 InstrumentSystems.CAS4 Namespace Reference	33
7.2.1 Detailed Description	33
8 Class Documentation	35
8.1 InstrumentSystems.CAS4.CAS4DLL Class Reference	35
8.1.1 Detailed Description	56
8.1.2 Member Function Documentation	56
8.1.2.1 casAssignDeviceEx()	56
8.1.2.2 casCalculateCorrectedData()	57
8.1.2.3 casCalculateCRI()	57
8.1.2.4 casCalculateLambdaDom()	58
8.1.2.5 casCalculateTOPParameter()	58
8.1.2.6 casChangeDevice()	59
8.1.2.7 casClearCalibration()	60
8.1.2.8 casClearDarkCurrent()	60
8.1.2.9 casColorMetric()	61
8.1.2.10 casConvoluteTransmission()	61
8.1.2.11 casCreateDevice()	61
8.1.2.12 casCreateDeviceEx()	62
8.1.2.13 casDeleteParamSet()	62
8.1.2.14 casDoneDevice()	63
8.1.2.15 casFIFOHasData()	63
8.1.2.16 casGetCalibrationFactors()	64
8.1.2.17 casGetCCT()	65
8.1.2.18 casGetCentroid()	66
8.1.2.19 casGetColorCoordinates()	66
8.1.2.20 casGetCRI()	67
8.1.2.21 casGetDarkCurrent()	67
8.1.2.22 casGetData()	68
8.1.2.23 casGetDeviceParameter()	68
8.1.2.24 casGetDeviceParameterString()	69
8.1.2.25 casGetDeviceTypeName()	69
8.1.2.26 casGetDeviceTypeOption()	70
8.1.2.27 casGetDeviceTypeOptionName()	70

8.1.2.28 casGetDeviceTypeOptions()	71
8.1.2.29 casGetDeviceTypes()	71
8.1.2.30 casGetDigitalIn()	71
8.1.2.31 casGetDigitalOut()	72
8.1.2.32 casGetDLLFileName()	72
8.1.2.33 casGetDLLVersionNumber()	73
8.1.2.34 casGetError()	73
8.1.2.35 casGetErrorMessage()	74
8.1.2.36 casGetExtendedColorValues()	74
8.1.2.37 casGetExternalADCValue()	75
8.1.2.38 casGetFIFOData()	77
8.1.2.39 casGetFilterName()	77
8.1.2.40 casGetMeasurementParameter()	78
8.1.2.41 casGetOptions()	78
8.1.2.42 casGetPeak()	79
8.1.2.43 casGetPhotInt()	79
8.1.2.44 casGetRadInt()	80
8.1.2.45 casGetSerialNumberEx()	80
8.1.2.46 casGetShutter()	81
8.1.2.47 casGetTriStimulus()	81
8.1.2.48 casGetWidth()	81
8.1.2.49 casGetWidthEx()	82
8.1.2.50 casGetXArray()	83
8.1.2.51 casInitialize()	83
8.1.2.52 casLoadTestData()	84
8.1.2.53 casMeasure()	84
8.1.2.54 casMeasureDarkCurrent()	85
8.1.2.55 casMultiTrackCopyData()	86
8.1.2.56 casMultiTrackCount()	86
8.1.2.57 casMultiTrackDone()	86
8.1.2.58 casMultiTrackInit()	87
8.1.2.59 casMultiTrackLoadData()	87
8.1.2.60 casMultiTrackSaveData()	88
8.1.2.61 casNmToPixel()	88
8.1.2.62 casPerformAction()	89
8.1.2.63 casPerformActionEx()	89
8.1.2.64 casPixelToNm()	89
8.1.2.65 casSaveCalibration()	90
8.1.2.66 casSaveSpectrum()	90
8.1.2.67 casSetCalibrationFactors()	91
8.1.2.68 casSetDeviceParameter()	92
8.1.2.69 casSetDeviceParameterString()	93

8.1.2.70 casSetDigitalOut()	93
8.1.2.71 casSetMeasurementParameter()	93
8.1.2.72 casSetOptions()	94
8.1.2.73 casSetOptionsOnOff()	94
8.1.2.74 casSetShutter()	95
8.1.2.75 casSetStatusLED()	95
8.1.2.76 casStart()	96
8.1.2.77 casUpdateCalibrations()	96
8.1.2.78 cmXYToDominantWavelength()	97
8.1.3 Member Data Documentation	97
8.1.3.1 mpidACQStateLine	97
8.1.3.2 mpidACQStateLinePolarity	98
8.1.3.3 mpidBusyStateLine	98
8.1.3.4 mpidBusyStateLinePolarity	99
8.1.3.5 mpidCurrentCCDTemperature	99
8.1.3.6 mpidFlashDelayTime	99
8.1.3.7 mpidForceFilter	100
8.1.3.8 mpidLastCCDTemperature	100
8.1.3.9 mpidRelSaturation	100
8.1.3.10 mpidTOPAperture	101
8.1.3.11 mpidTOPDistance	101
8.1.3.12 mpidTOPFieldOfView	101
8.1.3.13 mpidTriggerOptions	101
Index	103

Chapter 1

Introduction

This manual provides information about the functions and procedures of the CAS DLL.

Note

The documentation uses the static .NET class [CAS4DLL](#) to provide a reference for all methods and constants provided and used by the CAS DLL.

As a first step you should use the program CASDLLTEST.EXE to test if your spectrometer works and to get a first impression of the capabilities of the CAS DLL.

The manual is complemented by the examples provided in the subdirectories of the DLL installation directory. Here you find examples for different programming languages:

- A fully functional example written with C++ using *Microsoft Visual C++ 2013* can be found under *VC2013*. You will find the example source code, the CAS4.h header and the lib files for Win32 and Win64.
- In the subdirectory *CSharp* you find the C# solution CAS4DLLTest which uses the class [CAS4DLL](#).
- The subdirectory *Delphi* contains a Delphi sample project that demonstrates basic usage of the CAS DLL.

The DLL itself is installed to the system directory (usually C:\Windows\System32 or C:\Windows\SysWOW64 on Win64 systems).

The lib and header files for the different programming languages are in the directories of their respective examples (see above). If you are missing lib files for other compilers, use their toolchain to generate them from the DLL itself. One example would be C++ Builder, where you use implib.exe for Win32 and mkexp.exe for Win64.

Note

As a special feature the CAS DLL provides a demo mode interface. This interface enables you to test your application without having a spectrometer connected to your PC. The demo mode interface creates artificial noisy sinus spectra.

Chapter 2

To-Do List

This chapter gives an overview about the important topics that should be covered when integration the CAS DLL. It is intended as a check list that can be used to ensure that no vital aspect of the CAS DLL integration are forgotten.

Essential

Check	Chapter in this manual
	Device Handles and Interface Types
	Device Parameter and Measurement Parameter
	Thread Safety
	Error Handling
	Configuration and Initialization
	Density Filter
	Dark Current
	Performing a Measurement
	Getting the Spectrum and Results

Recommended

Check	Chapter in this manual
	Using Serial Numbers
	CCD Temperature Monitoring
	Synchronization

Optional

Check	Chapter in this manual
	Triggered Measurements
	AutoRange Measurements
	Working with Parameter Sets
	MultiTrack Measurements
	Transmission Measurement

Chapter 3

What has changed

This topic lists major changes that have been made to the interface.

casPerformActionEx

A new method [casPerformActionEx](#) has been introduced to support actions which required additional parameter. The old [casPerformAction](#) method has been marked deprecated and should no longer be used for new code and new actions.

If you are porting an application which previously used CAS DLL version 3.4 or earlier there are a few things to consider:

Integer Types

Throughout the DLL interface now only one 32bit integer type is used: int in C, Long in Visual Basic and Integer in Delphi. Note that this also applies to the 64bit version of the DLL which is called CAS4x64.DLL on Windows.

Unified Methods for Accessing Device and Measurement Parameter

Instead of introducing a new set of get- and set-methods for every new measurement or device parameter, generic methods have been introduced. See chapters [Device Parameter](#) and [Measurement Parameter](#).

Enhanced Support for Triggered Measurements

The CAS DLL can now return the triggering capabilities of the device - refer to the [dpidTriggerCapabilities](#) device parameter. If supported by the device, the current state of the acquisition now can also be signaled via DigitalOut ports - refer to the [mpidTriggerOptions](#) measurement parameter.

Chapter 4

Tasks

4.1 Interfaces Types and Options

Two basic aspects of a spectrometer instance within the DLL are the interface type and its option. The interface type defines how the spectrometer is connected to the computer and therefore by which means it is accessed (USB, PCI card etc.). For most interface types, the interface option defines a specific spectrometer that is connected to the PC with this interface type. However, for some interface types like PCI, the interface option only identifies the hardware interface card, not a specific spectrometer.

4.1.1 Device discovery

For the traditional interface types like PCI, PCIe and USB, connected spectrometers are automatically discovered in a synchronous fashion that happens implicitly when the interface options are enumerated. However, some interface types, like Ethernet, can only be discovered asynchronously. This discovery process can be initiated by calling [paSearchForDevices](#). Typically this is called very early during startup of your application and when new devices should be listed, for example whenever the user is presented with a user interface to select interface type and option as described below.

4.1.2 Letting the user select interface type and option

Typically the user is presented with two drop-down boxes where the interface type is selected first and then the available options are listed in the second drop-down.

Retrieve the number of interfaces with [casGetDeviceTypes](#) and their display names with [casGetDeviceTypeName](#). Ignore interface types with an empty display name, as they are unavailable. The number of options for a given interface is returned by [casGetDeviceTypeOptions](#), the value of the actual option with [casGetDeviceTypeOption](#) and the display name with [casGetDeviceTypeOptionName](#).

To get the interface type and option for a given CASID use the [dpidInterfaceType](#) and [dpidInterfaceOption](#) device parameter. If the device associated with a given interface type and option is not available, methods like [casInitialize](#) will return [ErrorAdrControl](#).

Typically the interface type and option properties are stored and restored as part of your applications settings, so that you can start working with the device right away after restarting.

To change the interface type and or option call [casChangeDevice](#).

4.1.3 Solutions without a UI

Some interface type constants are considered constant and it is not very likely that they change in future versions, the most prominent example being [InterfaceUSB](#). So to work with the first available USB spectrometer, you could retrieve the first available interface option for InterfaceUSB create a CAS ID for it with the following code:

[Delphi]

```
//before getting the interface option with index 0, ensure that the
//interface has enough options
intOptionCount = casGetDeviceTypeOptions(InterfaceUSB);
CheckError(intOptionCount);
if intOptionCount = 0 then
    raise Exception.Create('No interface options found for USB');
intOption:= casGetDeviceTypeOption(InterfaceUSB, 0);
...
CASID:= casCreateDeviceEx(InterfaceUSB, intOption);
```

Note

This approach will not work under all circumstances. It also does not detect if a different spectrometer has been connected since the last start (which would need other calibration files etc.). The best solution is still to produce an UI where the user can select interface type and option and then store and restore these values in configuration settings.

Warning

Be aware that interface options are not identical to spectrometers in all cases (e.g. for PCI) and that interfaces might also contain options for spectrometers that are no longer connected to the PC! This can have unexpected consequences when using the approach outlined above.

4.2 Device Handles (a.k.a. CASIDs) and Interfaces Types

When working with the CAS DLL you create objects within the DLL which store information about all aspects of the spectrometer, including hardware information as well as measurement parameters, the measured spectrum etc. These objects are identified via a device handle, which is often called a CASID.

Device Creation and Destruction

A device handle is typically created with [casCreateDeviceEx](#) where interface type and option are passed (see chapter [Interfaces Types and Options](#)).

Once you no longer need access to the spectrometer, you release the device handle with [casDoneDevice](#).

4.3 Thread Safety

Warning

The CAS DLL is not thread safe. Applications using the DLL must ensure that no two threads call a DLL method at the same time! Otherwise device operation might cause unpredicted behavior and/or result in unusual errors!

Typically thread-safety is ensured by the application with protecting every DLL method call by acquiring a global critical section, but any similar synchronization mechanism can be used.

Since every device handle also stores an error code, which is reset at the beginning of each method call, the [error handling](#) must also be done while the critical section is acquired.

4.3.1 Global settings

There are a few settings which are stored globally and therefore affect all devices. This might cause extra trouble in a multi-threaded application which operates several devices. These settings are:

- [mpidCRIMode](#) used when calling [casCalculateCRI](#)
- [mpidObserver](#) used when calling [casColorMetric](#), [casCalculateCRI](#) with [criCIE13_3_95](#), [cmXYToDominantWavelength](#) and [casGetCCT](#)

If these settings really need to be changed while several devices are operated, one possible approach would be to protect the above mentioned calls with a separate critical section to avoid that they are called simultaneously.

4.4 Error Handling

Most methods directly return an error code. For those who do not, the method [casGetError](#) can be called directly afterwards.

This error code can be translated into a error message by calling [casGetErrorMessage](#).

The following code contains a generic error handling method in Delphi which will raise an exception in case of an error. Different usage scenarios are also shown.

[Delphi]

```
procedure CheckError(AErrorCode: Integer);
var
  ach: array[0..255] of AnsiChar;
begin
  if AErrorCode < ErrorNoError then
  begin
    casGetErrorMessage(AErrorCode, ach, 255);
    raise ECASDLLException.CreateFmt('CAS DLL Error (%d): %s', [AErrorCode, StrPas(ach)]);
  end;
end;

...

//Usage with return value:
CheckError(casMeasureDarkCurrent(CasID));

//Usage with casGetError:
factor:= casGetCalibrationFactors(CasID, 25, 36);
CheckError(casGetError(CasID));

//Usage with casGetError too because return value cannot be considered an error code under all
circumstances
triggerOptions:= Round(casGetMeasurementParameter(CasID, mpidTriggerOptions));
CheckError(casGetError(CasID));
```

4.5 Device Parameter (dpid...)

Device parameter are properties or settings of the device and generally things which are not affected by or changed between measurements.

There is a set of four methods to set and retrieve device dependent parameter for a given spectrometer, two for numerical parameter and two for string parameter. The AWhat parameter defines the device parameter which should be retrieved or set. The possible values for AWhat are constants starting with dpid.

- [casGetDeviceParameter](#)
- [casSetDeviceParameter](#)
- [casGetDeviceParameterString](#)
- [casSetDeviceParameterString](#)

[Ansi C]

```
char buf[256];
int ret = casGetDeviceParameterString(CASID, dpidConfigFileName, buf, 255);
if (ret != ErrorNoError)
    goto fail;

...

char buf[256];
buf = ...
int ret = casSetDeviceParameterString(CASID, dpidConfigFileName, buf);
if (ret != ErrorNoError)
    goto fail;
```

4.6 Measurement Parameter (mpid...)

There is a set of two methods to set and retrieve numerical measurement parameter. These are either settings on how the device should measure, but can also document the measurement conditions.

- [casGetMeasurementParameter](#)
- [casSetMeasurementParameter](#)

Warning

The error handling shown in the example below is only possible, because [mpidIntegrationTime](#) is defined as a non-negative integer. With other mpids, it is necessary to call [casGetError](#) separately, to perform error checking, because the return value cannot safely be interpreted as an error code! For example [mpidFlash](#)↔[DelayTime](#) might return -14 as a valid negative delay time of 14ms and does not indicate [ErrorAdrControl](#)!

[Ansi C]

```
double ret = casGetMeasurementParameter(CASID, mpidIntegrationTime);
long intTime = lround(ret);
if (intTime < ErrorNoError)
    goto fail;
```

4.7 Configuration and Initialization

After setting up the device type and option it is essential that the device is configured. At the very least a configuration file (.ini) and a calibration file (.isc) have to be provided with the [dpidConfigFileName](#) and [dpidCalibFileName](#) [device parameter](#). These files typically come in pairs with the same name.

Some devices have configuration and calibration files stored on the device, see chapter [Working with Calibration files on the device](#).

When the configuration and calibration files are set, the device can be initialized with [casInitialize](#). Measurement parameter should only be set after [casInitialize](#) has been called, since only then the configuration has been loaded, which might affect parameter ranges etc.

4.8 Working with Calibration files on the device

You can download the files using [dpidGetFilesFromDevice](#). The path you set this dpid to should be an empty directory where the files will be downloaded to. If the device does not have files or doesn't support this at all, [dpidGetFilesFromDevice](#) will return errors like [ErrorInvalidEEPromType](#) or [ErrorNoFilesOnIdentKey](#).

Each downloaded calibration will consist of a configuration file (.ini), a calibration file (.isc) and possibly additional files. These files need to be used for [Configuration and Initialization](#). The chapter [Using Serial Numbers](#) describes techniques to ensure that the correct calibration is used together with the correct accessory.

4.8.1 Checking for updated calibrations

If a device is recalibrated, it is necessary to redownload the updated calibration files. To check whether the device has newer calibration files, the CAS library supports 2 methods:

- the check is performed during [casInitialize](#), if the option [coCheckConfigUpToDate](#) is enabled. Should an updated calibration be available, [casInitialize](#) would return [ErrorConfigUpToDate](#)
- you can perform the check manually at any time which is suitable, after the device has been initialized, by calling [paCheckConfigUpToDate](#), which will also return [ErrorConfigUpToDate](#), if there's an updated calibration

4.9 Density Filter

The density filter which should be used for the next measurement is defined by [mpidNewDensityFilter](#). Before every measurement, it is essential, that you check the [dpidNeedDensityFilterChange](#) device parameter. If it returns true [mpidDensityFilter](#) should be set to [mpidNewDensityFilter](#) which will turn the filter wheel to the correct position.

[Visual Basic]

```
Public Function MeasureCAS(ByVal ACasID As Integer)

    If casGetDeviceParameter(ACasID, dpidNeedDensityFilterChange) <> 0 Then
        Call casSetMeasurementParameter(ACasID, mpidDensityFilter, casGetMeasurementParameter(ACasID,
            mpidNewDensityFilter))
    End If

    ...

    MeasureCAS = casMeasure(ACasID)
End Function
```

Note

This ensures that all filter options and parameter are taken into account and that the filter wheel is not turned within a measurement cycle, which might interfere with timing restrictions or even destroy the DUT.

If you want to cater for spectrometers without a filter wheel, check whether the configuration supports a manual density filter. This corresponds to [dpidFilterType = 1](#). In this case your application should display a user message when [mpidDensityFilter](#) is changed (see example above). You might want to use [casGetFilterName](#) to display a user-friendly name for the filter.

Monitoring Filterwheel and Shutter operation

To determine if a given device supports shutter and/or filterwheel monitoring, use [casGetOptions](#) and check for the [coGetFilter](#) and [coGetShutter](#) bits.

Should shutter or filterwheel operation fail, the following actions may signal an error either via their return value or by setting error information (see [casGetError](#)).

- reading [mpidCurrentDensityFilter](#)
- setting [mpidDensityFilter](#)
- [casGetShutter](#)
- [casSetShutter](#)
- [casMeasure](#)
- [casMeasureDarkCurrent](#)

If you want to check for shutter and filter wheel errors at any given time, read [mpidCurrentDensityFilter](#) or call [casGetShutter](#).

The [ForceFilter](#) option is automatically enabled for devices which support shutter/filter-monitoring. This causes [dpidNeedDensityFilterChange](#) to always return true, so the filter should be set before each measurement (see section above).

4.10 Dark Current

Before a spectrum can be measured, it is necessary to measure the dark current (or DC) so it can be subtracted from the spectra later. This is done with [casMeasureDarkCurrent](#). Changing the [integration time](#) and or [averages](#) might make it necessary to measure a new dark current. To find out whether this is the case, check [dpidDCRemeasure↔Reasons](#). You can also set up a time after which a DC is marked invalid. Refer to the [coAutoRemeasureDC](#) option (see [casGetOptions](#)) and [mpidRemeasureDCInterval](#).

Use [casSetShutter](#) to close the shutter before calling [casMeasureDarkCurrent](#) and open it afterwards. If you want to display a user message for devices without a shutter, check the [coShutter](#) option with [casGetOptions](#).

Warning

It is absolutely *essential* to perform proper [error handling](#) for any of the above method calls!

[casMeasureDarkCurrent](#) only measures the dark current for the active parameter set. [paMeasureParamSetsDC](#) provides a convenient and efficient way to measure the dark current for some or all parameter sets - after closing the shutter and before opening it as described above. Refer to chapter [Working with Parameter Sets](#) for more information.

Warning

Before the dark current measurement is performed, it is necessary to move the filterwheel to the correct position! This helps avoiding overexposure for spectrometers where the filter wheel also functions as a shutter. For more information refer to the chapter [Density Filter](#).

It is also essential that at least one dark current measurement has been performed, before any spectrum is measured. Otherwise the resulting spectrum would be incorrect, one reason being that the [AmpOffset](#) is determined during a dark current measurement.

An existing dark current can be removed by calling [casClearDarkCurrent](#), but that is usually not necessary. Note that clearing the dark current will not clear the [AmpOffset](#) mentioned above.

4.10.1 Dark current array

A dark current array is supported for all device types. It contains the dark current spectra for various integration times which are measured at once, so that switching to a different integration time does not require a new dark current measurement. This is especially useful for [AutoRange](#) measurements and for spectrometers without a shutter.

Note

With most spectrometers now having a shutter and with introduction of a new [AutoRange](#) variant which requires less DC measurements (see [AutoRange for various intensities](#)), the dark current array has less and less advantages. With the introduction of the new `paMeasureParamSetsDC` (refer to the first part of this chapter), there is now also no longer a need for using the dark current array, even if you have a large number of parameter sets.

The dark current for a given integration time is interpolated from the dark current array. To activate the dark current array, set the `coUseDarkcurrentArray` option (see [casGetOptions](#)). Use [casGetCalibrationFactors](#) to access the actual dark current array. Measuring a dark current array is done the same way as a normal dark current measurement: close the shutter, call [casMeasureDarkCurrent](#) and then open the shutter.

Warning

In rare cases, measuring a dark current array can result in [ErrorDarkArray](#). This indicates that one of the dark current array spectra could not be measured with exactly the desired integration time, resulting in an invalid array. In such cases the dark current array can not be used and `coUseDarkcurrentArray` should be switched off.

Note

If the dark current array has been successfully measured, [dpidDCRemeasureReasons](#) will no longer signal a required dark current measurement when changing the integration time or switching the current parameter set. This might however still happen, if the array hasn't been measured or the range of the dark current array does not cover the integration time to be measured.

Warning

Before the dark current measurement is performed, it is necessary to move the filterwheel to the correct position. This helps avoiding overexposure for spectrometers where the filter wheel also functions as a shutter. For more information refer to the chapter [Density Filter](#).

4.11 Performing a Measurement

A spectrum can be measured by calling [casMeasure](#). This method only returns after the acquisition has been completed.

The following preconditions should be met before `casMeasure` is called:

- The device handle has been created, see [Device Handles \(a.k.a. CASIDs\) and Interfaces Types](#)
- The device has been configured and initialized, see [Configuration and Initialization](#)
- The integration time has been set, see [mpidIntegrationTime](#), or [AutoRange](#) has been enabled
- The density filter wheel has been moved if necessary, i.e. setting [mpidDensityFilter](#) to [mpidNewDensityFilter](#) if [dpidNeedDensityFilterChange](#) returns true.
- The dark current has been measured, if [dpidDCRemeasureReasons](#) is not 0. Refer to chapter [Dark Current](#)
- [casPerformActionEx](#) has been called with `paPrepareMeasurement`

A lot of measurement setups require that the CAS is synchronized with other equipment like a handler or a sourcemeter. For an overview of possible synchronization setups, see chapter [Synchronization](#).

Warning

Validate the measurement by checking that the spectrometer wasn't oversaturated: `mpidMaxADCValue` must not be larger than `dpidADCRange`. If it reaches `dpidADCRange`, this should be considered a hard error.

Another recommended validation is to check `mpidRelSaturation` against `dpidRelSaturationMin` and `dpidRelSaturationMax` to make sure it was sufficiently saturated. If this check fails, it does not have to be considered a hard error, but rather a good/bad indicator.

Note that checking `mpidRelSaturation` alone is not sufficient to check for over-saturation, especially when averaging is used!

Do *not* calculate a relative signal level yourself using `mpidMaxADCValue` as this will be a misleading indicator! `mpidRelSaturation` takes the DC into account and is therefore better suited.

Obviously **error handling** is essential to detect hardware failures and other problems. Without error checking there would be no indication that the measurement failed and the previously measured spectrum would be used.

It is highly recommended to retrieve all measurement parameter after the measurement and document them with your measurement data, because some measurement parameter might have been adjusted. The most prominent example would be `mpidIntegrationTime` which is affected by `AutoRange`, `mpidIntTimeAlignPeriod`, `mpidIntTimeResolution` and for some spectrometers also `dpidIntTimeMax`.

4.12 Getting the Spectrum and Results

After a spectrum has been measured, the intensity values can be retrieved by calling `casGetData` for every visible pixel. The data unit for these values can be retrieved with `dpidCalibrationUnit` once `casInitialize` has been called. `casGetXArray` provides the corresponding wavelength values and `casGetDarkCurrent` returns the dark current per pixel. Make sure that `dpidDeadPixels` and `dpidVisiblePixels` are taken into account, i.e. start with `dpidDeadPixels` and end with `dpidDeadPixels + dpidVisiblePixels - 1` (since the pixel index is 0-based).

4.12.1 ColorMetric results

The CAS SDK can calculate a wide range of colormetric results from a previously measured spectrum. These calculations only use the part of the spectrum defined by the `mpidColormetricStart` and `mpidColormetricStop` wavelength range.

Warning

The default values for `mpidColormetricStart` and `mpidColormetricStop` always have been 380nm and 780nm. It is highly recommended to adjust these values as they might constrain the calculation unexpectedly. Set both of these to zero to use the full range of the spectrum for calculation.

Starting with version 4.10 of the CAS SDK, a new `mpidColormetricType` parameter has been introduced. It allows to optionally switch to the more accurate but also more time-consuming colormetric calculation of SpecWin Pro. `mpidColormetricType` only applies to the given CAS ID but is independent of its **active parameter set**.

`casColorMetric` performs a calculation of several colormetric results from the measured spectrum which can then be retrieved with the following methods.

The photometric and radiometric integrals can be retrieved with `casGetPhotInt` and `casGetRadInt`.

`casGetCentroid` returns the centroid wavelength whereas `casGetPeak` and `casGetWidth` return results related to the peak.

Use `casGetColorCoordinates` to retrieve the CIE color coordinates and `casCalculateLambdaDom` to calculate dominant wavelength and purity.

To get the correlated color temperature (CCT) call `casGetCCT` and to get the tristimulus values call `casGetTriStimulus`.

Further results can be retrieved with `casGetExtendedColorValues`, which include the VIS, UVA, UVB and UVC integrals.

If you need the color rendering indices (CRI), call `casCalculateCRI` to calculate them and `casGetCRI` to retrieve them. Note that the mode of the CRI calculation can be set globally with `mpidCRIMode`.

4.13 Synchronization

In most measurement setups some kind of synchronization is needed to ensure that the CAS measurement takes place at the right point in time. Most likely because of a sourcemeter or handler powering a DUT.

Note

For reliable measurement results, it is essential that the CAS measurement happens only after the DUT is powered on and that the DUT is powered on long enough.

The following table lists possible synchronization methods with their advantages and disadvantages. It starts with the ideal but challenging method and ends with the easiest but most problematic method.

Method	Who triggers whom	Advantages	Disadvantages
<ul style="list-style-type: none"> • trigger handshake recommended! 	<ul style="list-style-type: none"> • CAS triggers sourcemeter start • Sourcemeter triggers CAS start • optional: CAS triggers sourcemeter power-off 	<ul style="list-style-type: none"> • most accurate timings and therefore results • minimal cycle time • minimal power-on time for DUT • trigger delay and settling time possible 	<ul style="list-style-type: none"> • challenging requirements
<ul style="list-style-type: none"> • triggered CAS might be OK 	<ul style="list-style-type: none"> • sourcemeter triggers CAS 	<ul style="list-style-type: none"> • easier to set up 	<ul style="list-style-type: none"> • CAS might not be ready for trigger¹ • longer cycle time • risk of DUT powering off too early
<ul style="list-style-type: none"> • software triggering not recommended 	<ul style="list-style-type: none"> • no triggers 	<ul style="list-style-type: none"> • easiest to set up 	<ul style="list-style-type: none"> • longest cycle time • longest power-on time for DUT • risk of DUT damage² • a varying settling time is unavoidable and might influence results

1. It might be OK to work around this by adding a generous delay before the sourcemeter triggers the CAS. Some CAS types also support the [toAcceptOnlyWhenReady](#) trigger option which would at least help detect a missed trigger.


2. Risk of DUT damage: if [preparation of measurement](#) has not been done properly or if [error handling](#) is missing, the DUT might be powered on for too long, which might cause damage.

4.14 Trigger handshake sequence

The following diagram illustrates a full trigger handshake sequence as described in the [synchronization](#) chapter.

Most steps in the sequence refer to chapters in this documentation with the same name.

See chapter [triggered measurements](#) for details on trigger options and capabilities.



4.14.1 For this sequence to work

- the [dpidTriggerCapabilities](#) has to include tcoCanTrigger

- similarly the toShowACQState bit of [mpidTriggerOptions](#) has to be set, toAcceptOnlyWhenReady is optional but recommended
- the ACQState line settings should be verified: by checking [mpidACQStateLine](#) and [mpidACQStateLine↔Polarity](#). **Not all CAS types support modifying these mpid's!**

4.14.2 Several measurements per DUT

- if several measurements per DUT are needed, use [parameter sets](#)
- for each parameter set, set [dpidCurrentParamSet](#) and then call casMeasure (step 11 to 17)
- after all measurements have been done, again for each parameter set, set [dpidCurrentParamSet](#) and then get spectrum and results (step 19)
- whether you need to start the sourcemeter for each parameter set depends on the sourcemeter and it's capabilities
- note that the parameter sets need to be defined together with setting the measurement parameter (step 4)
- the Dark Current also needs to be measured for all parameter sets

Remarks


- if the sourcemeter does not support being switched off by the CAS (step 15 in the sequence diagram), a possible workaround is to switch off the DUT after a given time. This time - which should be verified by tests - should cover [mpidIntegrationTime](#) * [mpidAverages](#) plus a few milliseconds extra time. The same approach is taken in the [trigger CAS sequence](#).
- in the sequence diagram, the red rectangle on the CAS line shows the time during which the ACQ State line has the polarity defined by [mpidACQStateLinePolarity](#). As soon as the DUT can be switched off, it changes back to the opposite polarity.

4.15 Trigger CAS sequence

The following diagram illustrates the trigger CAS sequence as described in the [synchronization](#) chapter.

Most steps in the sequence refer to chapters in this documentation with the same name.

See chapter [triggered measurements](#) for details on trigger options and capabilities.



4.15.1 For this sequence to work

- the [dpidTriggerCapabilities](#) has to include tcoCanTrigger
- the toAcceptOnlyWhenReady bit of [mpidTriggerOptions](#) should be set, but some CAS types do enforce this bit automatically. It is recommended to verify this bit afterwards by reading mpidTriggerOptions

4.15.2 Several measurements per DUT

- if several measurements per DUT are needed, use [parameter sets](#)

- for each parameter set, set [dpidCurrentParamSet](#) and then call `casMeasure` (step 11 to 15)
- after all measurements have been done, again for each parameter set, set [dpidCurrentParamSet](#) and then get spectrum and results (step 18)
- whether you need to start the sourcemeter for each parameter set depends on the sourcemeter and its capabilities
- note that the parameter sets need to be defined together with setting the measurement parameter (step 4)
- the Dark Current also needs to be measured for all parameter sets


Remarks

- the Sourcemeter delay time (between steps 10 and 11 in the diagram) should be verified by tests. Increase it until you no longer get missed triggers, i.e. [ErrorTriggerTimeout](#).
- the DUT power on time (between steps 13 and 16) should also be verified by tests. Start with [mpid](#)↔[AutoFlowTime](#) as this covers [mpidIntegrationTime](#) * [mpidAverages](#). You might need to add some time until it no longer influences CAS measurement results.

4.16 Software trigger sequence

The following diagram illustrates the software trigger sequence as described in the [synchronization](#) chapter.

Most steps in the sequence refer to chapters in this documentation with the same name.



Remarks

- proper preparation of the measurement (steps 4 to 7 in the diagram) is even more essential here, because an incomplete preparation could extend the unknown time between steps 10 and 11 significantly
- to avoid the risk of DUT damage because of a longer than expected power on time, the DUT could be switched off after a given time, similar to the approach shown in the [trigger CAS sequence](#). But this power on time now also needs to cover the unknown time before measurement starts (between steps 9 and 10), so it must be chosen and verified carefully.

4.17 Triggered Measurements

Triggered measurements using the CAS DLL are done by setting the `mpidTriggerSource` parameter to `trgFlipFlop`. A delay between trigger and acquisition can be set with `mpidTriggerDelayTime`. The trigger timeout should be set appropriately with `mpidTriggerTimeout`.

Some spectrometer types support signaling their state via DigitalOut ports which can help triggering the current source - see [mpidTriggerOptions](#) and the [Synchronization](#) chapter.

Warning

For some interface types a call to [casPerformActionEx](#) with `paPrepareMeasurement` is essential for triggered measurements to work!

Make sure to test whether the spectrometer supports the trigger aspects you want to use! Most importantly check the [dpidTriggerCapabilities](#), but also things like [dpidTriggerDelayTimeMin](#) and [dpidTriggerDelayTimeMax](#).

Be aware that the DigitalOut ports can change their state during the preparation phase! For example, as soon as the [tcoShowACQState](#) trigger option is activated, the [mpidACQStateLine](#) reflects the current state of [mpidACQStateLinePolarity](#). So if the option is activated before the polarity is set, the level will change, which might result in an unwanted trigger. The safest approach is to only arm a sourcemeter triggered by the CAS *after* all the preparation tasks have been completed. This is especially important since dark current measurements could also result in unwanted triggers.

The following example checks if the device supports signaling the ACQ state via DigitalOut ports:

[Visual Basic]

```
Dim Capabilities As Long

Capabilities = Round(casGetDeviceParameter(CASID, dpidTriggerCapabilities))

If (Capabilities And tcoShowACQState) <> 0 Then
    MsgBox "Device supports signaling ACQ state"
Else
    MsgBox "Device does NOT support signaling ACQ state"
End If
```

4.18 AutoRange measurements

For some measurement tasks it might be appropriate to adjust the integration time of the CAS automatically. The CAS Library provides an AutoRange feature that determines the proper integration time and optionally can also change to the appropriate density filter, if necessary. After the AutoRange measurement, these parameter might have been changed and should be retrieved to document the measurement conditions.

Starting with version 4.7 the CAS Library supports two different ways of AutoRange. The one which has been supported for many years now, is aimed at measuring samples with similar intensities. The new AutoRange method is more suitable for various intensities and can significantly reduce the time and number of measurements which have to be made before suitable settings are found.

4.18.1 AutoRange parameter

Use [mpidAutoRangeMinLevel](#), [mpidAutoRangeMaxLevel](#) and [mpidAutoRangeMaxIntTime](#) to specify limits for AutoRange. You can query [dpidAutoRangeFilterMin](#) and [dpidAutoRangeFilterMax](#) to check the density filter range which will possibly be used.

Warning

[mpidAutoRangeMinLevel](#) and [mpidAutoRangeMaxLevel](#) are simple percentages of the complete [dpidADCRange](#). When doing an AutoRange measurement, [mpidMaxADCValue](#) is checked against this range to verify that a sufficient integration time has been found. This means that the DC is not taken into account. So if you specify a low minimal range of e.g. 10%, these 10% - for long integration times - might actually be reached just by the DC itself, so there's no real signal. This is contrary to [mpidRelSaturation](#) where the DC is taken into account, so the percentages of [mpidAutoRangeMinLevel](#) and [mpidAutoRangeMaxLevel](#) cannot be compared with [mpidRelSaturation](#).

4.18.2 AutoRange for similar intensities

To enable AutoRange for similar intensities, use the command `casSetOptionsOnOff` and enable the option `coAutoRangeMeasurement` and, if desired, `coAutoRangeFilter`.

To perform the actual AutoRange measurement, just call `casMeasure` like you would for a normal measurement. This will start a dark-current measurement, if necessary, and then a normal measurement with the current measurement parameter. If the results meet the criteria given via the `AutoRange parameter`, `casMeasure` returns.

So for similar intensities, once suitable parameter have been found during the initial AutoRange measurement, the following calls to `casMeasure` can return after the first acquisition.

If the first measurement does not meet the AutoRange criteria, the integration time and density filter are changed accordingly, a new DC is measured (or calculated from the DC array, see warning below) and then another measurement is made. This will repeat several times until the criteria are met or AutoRange gives up. Note that this procedure will be especially slow if a high number of `averages` are set up.

Warning

An AutoRange measurement for similar intensities might include necessary dark current measurements. For devices without a shutter, it is therefore mandatory to enable the dark current array (see chapter [Dark Current](#)). Otherwise the AutoRange measurement might measure a dark current without the optical path being closed! For the dark current array to actually avoid shutter action, the `mpidAutoRangeMaxIntTime` must also stay within the bounds of the dark array. The integration times of the dark current array can be retrieved with `casGetCalibrationFactors` using `AWhat = gcfDarkArrayValues`.

Note

If you're measuring DUTs which vary greatly you should consider using the new AutoRange method for various intensities. Another approach is to speed up AutoRange measurements by setting `mpidIntegrationTime` and `Density Filter` to low values before measuring. That way you are avoiding a lengthy and over-exposed first acquisition if the previous AutoRange measurement changed these parameter because the DUT had a very low intensity.

Warning

AutoRange is only supported for synchronous measurements using `casMeasure`, but not for asynchronous measurements with `casStart`!

4.18.3 AutoRange for various intensities

For AutoRange with various intensities, simply call `casPerformActionEx` with `paAutoRangeFindParameter` before doing the actual measurement.

`paAutoRangeFindParameter` will start measurements, beginning with `dpidIntTimeMin` and `dpidAutoRangeFilterMin`, but to save time, without averages and without measuring any dark current. Once sufficient integration time and density filter have been found, `paAutoRangeFindParameter` will return (or fail with an error). After that you should do a [dark current measurement](#) if necessary and then perform the measurement like you would normally do, typically by calling `casMeasure`.

To allow `paAutoRangeFindParameter` to change the density filter, enable the option `coAutoRangeFilter` using `casSetOptionsOnOff`. It is NOT recommended to enable `coAutoRangeMeasurement`, because a subsequent measurement would then try to find suiting parameter again.

4.19 Working with Parameter Sets

The CAS DLL provides so called parameter sets. By using parameter sets it is possible to quickly change between different sets of measurement parameter. Since each parameter set holds it's own dark current, spectrum and results, this makes it possible to quickly measure several spectra with varying parameter without the need of dark current measurements in between.

The following parameters can be changed for each defined parameter set:

- [density filter](#)
- [integration time](#)
- [averages](#)
- [trigger delay time](#)
- [flash delay time](#)
- [flash duration](#)
- [CheckStart](#) and [CheckStop](#)
- [ColormetricStart](#) and [ColormetricStop](#)
- [GetWidth Level](#)
- [TOP aperture](#)
- [TOP distance](#)
- [PulseWidth](#)

The following commands and parameter are used to manage parameter sets

- [dpidParamSets](#) - defines the number of parameter sets. When increasing [dpidParamSets](#), new parameter sets are created by copying the currently active parameter set. If the currently active parameter set is deleted because [dpidParamSets](#) was decreased, the first parameter set will become active.
- [dpidCurrentParamSet](#) is the currently active parameter set. This device parameter has to be set prior to every action that is related to measurement parameters or measurement data. It is often used inside a loop to cycle through all parameter sets (see Example below)
- [casDeleteParamSet](#) is used for deleting a specific parameter set, since decreasing the number of parameter sets only deletes the sets at the end

Starting with CAS SDK version 4.10 the new [paMeasureParamSets](#) API can be used to conveniently measure a range of parameter sets. Some spectrometer types even support performing this set of measurements completely autonomously, a feature that is sometimes referred to as recipe mode.

For [dark current measurements](#) use [paMeasureParamSetsDC](#). As with normal dark current measurements, the shutter should be closed before calling the API and opened afterwards.

Note

[paMeasureParamSetsDC](#) avoids measuring multiple DCs with the same settings by copying the DCs into the appropriate parameter sets instead of measuring again.

Both [paMeasureParamSets](#) and [paMeasureParamSetsDC](#) are synchronous, i.e. they only return once the operation has been completed or an error occurred. In case of an error, use [mpidMeasured](#) to determine whether a given parameter set was measured successfully or not.

4.20 Using Serial Numbers

For many years now, all spectrometers by Instrument Systems have their serial number stored internally so that it can be retrieved by software.

Additionally a lot of accessories have their own serial number. If these are connected to the spectrometer correctly, these accessory serial numbers can be retrieved by software as well. To check whether a spectrometer supports reading accessory serial numbers, check the `coGetAccessories` bit with [casGetOptions](#).

The general serial number [dpidSerialNo](#) can consist of up to 3 serial numbers divided by a semicolon. First the serial of the spectrometer, then that of additional accessory and finally the serial number of a TOP, if it is required by the current configuration/calibration. To retrieve the individual serial numbers, use [casGetSerialNumberEx](#).

Warning

You should use the function [casPerformActionEx](#) with `paCheckAccessories` to ensure that the currently used calibration and configuration is intended for this spectrometer and optional accessory. If you want this check to be performed automatically, set the option `coCheckCalibConfigSerials` using the command [casSetOptions↵](#) [OnOff](#). If this option is enabled, [casInitialize](#) may also return error code [ErrorConfigSerialMismatch](#) or [Error↵](#) [CalibSerialMismatch](#) indicating that either the configuration or the calibration file doesn't match the serial number(s) of the device. Additionally the CAS DLL can now verify the serial of a transmission file ([dpid↵](#) [TransmissionFileName](#)). To perform this verification call [casPerformActionEx](#) with `paCheckAccessories↵` `Transmission`. In case of a mismatch it will return [ErrorTransmissionSerialMismatch](#).

Note

All the relevant files from Instrument Systems contain the serial number(s) they are made for. If these files don't contain the intended serial number(s), the `paCheckAccessories` and `paCheckAccessoriesTransmission` mentioned above can never fail.

The ID Key is only read during [casInitialize](#), when validating the serials manually with [casPerformActionEx](#) and `pa↵` `CheckAccessories` or `paCheckAccessoriesTransmission`, when reading [dpidAccessorySerial](#) or when calling [cas↵](#) [GetSerialNumberEx](#) with `casSerialAccessory`.

Warning

If you plug in the ID Key after the device has been initialized, it is not re-read when getting [dpidSerialNo](#) or calling [casGetSerialNumberEx](#) with `casSerialComplete`. Instead use one of the methods mentioned above.

4.21 MultiTrack measurements

The CAS DLL is able to perform so called MultiTrack measurements. These are very fast continuous measurements without any acquisition gaps. This is possible because the acquired raw values are copied into an internal buffer and any correction or calculation can be performed afterwards.

Note

Not all interface types support MultiTrack. Interfaces that do support MultiTrack will include the `coCanMultiTrack` flag in `casGetOptions`.

Starting with CAS DLL 4.5 MultiTrack supports more interface types - even the Demo mode now supports it. Note that the necessary calls have changed (see below). The old way of doing MultiTrack, which involved calling `casMultiTrackReadData` and `casMultiTrackCopySet` is deprecated now and will only work with the PCI interface for the time being.

To perform a MultiTrack measurement these tasks have to be done in the order listed below:

1. All preconditions that are mentioned in [Performing a Measurement](#) are met.
2. Call `casMultiTrackInit` to initialize the internal buffer and specify the number of MultiTrack measurements. This number must not exceed `dpidMultiTrackMaxCount`.
3. You may check the actual size of the buffer with `dpidMultiTrackCount`.
4. Start the MultiTrack measurement by calling `casPerformActionEx` with `paMultiTrackStart`. This will return once all measurements have been made. Triggering the first track is supported.
5. The MultiTrack measurements can be saved to a file with `casMultiTrackSaveData` and loaded using `casMultiTrackLoadData`.
6. For calculating colormetric and other results, call `casMultiTrackCopyData` for each track. This reloads the raw spectrum from the MultiTrack buffer and all other methods use it as if it would have just been measured. After that, proceed as described in the chapter [Getting the Spectrum and Results](#).
7. Finally call `casMultiTrackDone` to release the memory that was allocated for the MultiTrack buffer.

Important notes about MultiTrack files:

- We recommend using the `.swm` file extension
- Before calling `casMultiTrackLoadData` it is essential that the same configuration and calibration files were used to initialize the device! It is perfectly fine to use a demo mode CAS ID, i.e. for loading the spectrometer hardware does not need to be present.
- After loading and working with MultiTrack data, a call to `casMultiTrackDone` is necessary to release the MultiTrack buffer.
- Some MultiTrack specific results are `mpidMultiTrackAcqTime` and `mpidCMTTrackStart`.

4.22 CCD Temperature Monitoring

Various spectrometer types support measuring the CCD temperature, indicated by returning the `coGetTemperature` flag when calling `casGetOptions`.

For these spectrometers, it is highly recommended that you start a temperature measurement from time to time by querying `mpidCurrentCCDTemperature`. Since measuring the temperature might take a few milliseconds, choose a time which is convenient for you, for example after a measurement has been performed and the DUT has been disconnected.

This will ensure that hardware defects are detected (which would cause an `ErrorCCDTemperatureFail`) and also that, if necessary, a new dark current measurement will be requested (via `dpidDCRemeasureReasons`, as described in chapter [Dark Current](#)).

If you want to log the last measured temperature, without actually causing a temperature measurement, query `mpidLastCCDTemperature`. Note that this caching only works if the last temperature measurement is not too old (time frame varies by spectrometer type). Otherwise the temperature is actually measured, which might cause a delay.

Whenever a dark current is measured, the CCD temperature is also acquired. You can retrieve it by querying `mpidDCCCDTemperature`.

4.23 Transmission measurement

This section describes how to perform a transmission measurement where the reference spectrum is measured (as opposed to being loaded from a file etc.).

Warning

This topic should not be confused with taking the transmission of additional optical equipment into account! That can be done using [dpidTransmissionFileName](#) and enabling `coUseTransmission` with [casSetOptions](#)↔[OnOff](#).

Preparations

Since a calibration file is required before calling [casInitialize](#), it is recommended to use a 1-calibration, typically called "one.isc". Alternatively use the original calibration file and after `casInitialize`, call [casClearCalibration](#) with `gcfSensitivityFunction` like in the example provided below. Then set up the measurement parameter and perform a [dark current measurement](#).

Measuring the reference and applying it as a calibration

The reference measurement itself can be started like any normal measurement, either synchronously ([casMeasure](#)) or asynchronously ([casStart](#) with [casFIFOHasData](#) and [casGetFIFOData](#)).

To automatically divide future transmission measurements by the reference spectrum, it is stored as the sensitivity calibration using [casSetCalibrationFactors](#):

[Delphi]

```
//CasID holds the device handle
//the device must have been configured and initialized

//clear the sensitivity to make sure a previous
//reference is not applied
casClearCalibration(CasID, gcfSensitivityFunction);
CheckError(casGetError(CasID));

//set density filter if necessary
if Round(casGetDeviceParameter(CasID, dpidNeedDensityFilterChange)<>0 then
  CheckError(casSetMeasurementParameter(CasID, mpidDensityFilter, casGetMeasurementParameter(CasID,
    mpidNewDensityFilter)));

//dark current measurement if necessary
if Round(casGetDeviceParameter(CasID, dpidDCRemeasureReasons)<>0 then
begin
  casSetShutter(CasID, 1);
  CheckError(casGetError(CasID));

  CheckError(casMeasureDarkCurrent(CasID));

  casSetShutter(CasID, 0);
  CheckError(casGetError(CasID));
end;

//perform the reference measurement
casMeasure(CasID);
CheckError(casGetError(CasID));

//set sensitivity calibration to reference spectrum
//loop through all Pixels
for i:= 0 to Round(casGetDeviceParameter(CasID, dpidPixels))-1 do
begin
  casSetCalibrationFactors(CasID, gcfSensitivityFunction, i, 0, casGetData(CasID, i));
  CheckError(casGetError(CasID));
end;

//make sure the new "calibration" is taken into account
casPerformActionEx(CasID, paUpdateSpectralCalibration, 0, 0);
CheckError(casGetError(CasID));
```


Transmission measurements

After the reference spectrum has been assigned as sensitivity calibration, every subsequent measurement has the reference spectrum automatically applied. To clarify: the CAS DLL divides the spectrum with `gcfSensitivityFunction`, so the spectrum will be divided by the reference, resulting in a transmission spectrum.

Note

A new initialization (see [casInitialize](#)) or [paUpdateCompleteCalibration](#) will restore the original calibration from the calibration file.

Replacing the reference spectrum

To replace a previous reference spectrum, either re-initialize the device (see note above) or call [casClearCalibration](#) with `gcfSensitivityFunction`. This is done at the beginning of the example above, to make sure the old reference spectrum is not taken into account when a new reference is measured.

Chapter 5

Namespace Index

5.1 Packages

Here are the packages with brief descriptions (if available):

InstrumentSystems	33
InstrumentSystems.CAS4 The InstrumentSystems namespace.	33

Chapter 6

Class Index

6.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

InstrumentSystems.CAS4.CAS4DLL	
The CAS DLL interface class	35

Chapter 7

Namespace Documentation

7.1 InstrumentSystems Namespace Reference

Namespaces

- namespace [CAS4](#)
The [InstrumentSystems](#) namespace.

7.2 InstrumentSystems.CAS4 Namespace Reference

The [InstrumentSystems](#) namespace.

Classes

- class [CAS4DLL](#)
The CAS DLL interface class

7.2.1 Detailed Description

The [InstrumentSystems](#) namespace.

Chapter 8

Class Documentation

8.1 InstrumentSystems.CAS4.CAS4DLL Class Reference

The CAS DLL interface class

Static Public Member Functions

- static int [casGetError](#) (int ADevice)
Return error code for a given device/CASID
- static IntPtr [casGetErrorMessage](#) (int AError, StringBuilder ADest, int AMaxLen)
Translates a given error code into a readable error message
- static int [casCreateDevice](#) ()
Deprecated method. Creates a device context aka CASID within the CAS DLL
- static int [casCreateDeviceEx](#) (int AInterfaceType, int AInterfaceOption)
Creates a device context aka CASID within the CAS DLL
- static int [casChangeDevice](#) (int ADevice, int AInterfaceType, int AInterfaceOption)
Change the interface type and/or option of a device / CASID
- static int [casDoneDevice](#) (int ADevice)
Release resources used by the device
- static int [casAssignDeviceEx](#) (int ASourceDevice, int ADestDevice, int AOption)
Assigns properties and or parameters from one device handle to another
- static int [casGetDeviceTypes](#) ()
Retrieves the number of interface types the CAS DLL supports
- static IntPtr [casGetDeviceTypeName](#) (int AInterfaceType, StringBuilder ADest, int AMaxLen)
Retrieves the name of the given interface type
- static int [casGetDeviceTypeOptions](#) (int AInterfaceType)
Retrieves the number of options a given interface type currently supports
- static int [casGetDeviceTypeOption](#) (int AInterfaceType, int AIndex)
Returns the value of the interface option for the given interface type and option index
- static IntPtr [casGetDeviceTypeOptionName](#) (int AInterfaceType, int AInterfaceOptionIndex, StringBuilder ADest, int AMaxLen)
Retrieves the name of the given interface option
- static int [casInitialize](#) (int ADevice, int Perform)
Initializes the hardware of the device after loading the configuration and calibration files
- static double [casGetDeviceParameter](#) (int ADevice, int AWhat)

- Retrieves a float representing a device parameter*

 - static int [casSetDeviceParameter](#) (int ADevice, int AWhat, double AValue)

Sets a numerical device parameter
- static int [casGetDeviceParameterString](#) (int ADevice, int AWhat, StringBuilder ADest, int AMaxLen)

Retrieves a string representing a device parameter
- static int [casSetDeviceParameterString](#) (int ADevice, int AWhat, string AValue)

Sets a string device parameter
- static int [casGetSerialNumberEx](#) (int ADevice, int AWhat, StringBuilder ADest, int AMaxLen)

Retrieves serial numbers of the specified device and/or additional information.
- static int [casGetOptions](#) (int ADevice)

Returns the features and options the device currently supports or performs.
- static void [casSetOptionsOnOff](#) (int ADevice, int AOptions, int AOnOff)

This method can set or clear several options for the device
- static void [casSetOptions](#) (int ADevice, int AOptions)

This method sets and clears all device options for the device
- static int [casMeasure](#) (int ADevice)

Performs a measurement for the given device using the measurement parameter which have been previously set
- static int [casStart](#) (int ADevice)

Starts a measurement for the given device and returns immediately
- static int [casFIFOHasData](#) (int ADevice)

Use this method to check if the spectrometer has data available which can be read
- static int [casGetFIFOData](#) (int ADevice)

This method reads the acquired spectrum from the FIFO and stores it internally
- static int [casMeasureDarkCurrent](#) (int ADevice)

Use casMeasureDarkCurrent to perform a dark current measurement for the given device.
- static int [casPerformActionEx](#) (int ADevice, int AActionID, int AParam1, int AParam2, IntPtr AParam3)

Generic method which performs one of various actions with the specified device. This method replaces the now deprecated method casPerformAction.
- static int [casPerformAction](#) (int ADevice, int AID)

Deprecated. Use casPerformActionEx instead!
- static double [casGetMeasurementParameter](#) (int ADevice, int AWhat)

Returns a numeric measurement parameter for a given spectrometer
- static int [casSetMeasurementParameter](#) (int ADevice, int AWhat, double AValue)

Sets a numeric measurement parameter for a given spectrometer
- static int [casClearDarkCurrent](#) (int ADevice)

Clears a previously measured dark current of the given device
- static int [casDeleteParamSet](#) (int ADevice, int AParamSet)

Deletes a specific parameter set for the given device.
- static int [casGetShutter](#) (int ADevice)

Returns the current shutter position of the given device
- static void [casSetShutter](#) (int ADevice, int OnOff)

Sets the shutter position of the given device.
- static IntPtr [casGetFilterName](#) (int ADevice, int AFilter, StringBuilder Dest, int AMaxLen)

Translates a filter index into a user-readable density filter name.
- static int [casGetDigitalOut](#) (int ADevice, int APort)

Returns the state the specified digital output port was set to by a former call to [casSetDigitalOut](#).
- static void [casSetDigitalOut](#) (int ADevice, int APort, int OnOff)

Sets the state the specified digital output port.
- static int [casGetDigitalIn](#) (int ADevice, int APort)

Returns the state of the specified digital input port.
- static void [casCalculateCorrectedData](#) (int ADevice)

- Applies the spectral correction to the previously acquired raw spectrum.*

 - static void [casConvoluteTransmission](#) (int ADevice)
- Applies the spectral and transmission correction to the previously acquired raw spectrum.*

 - static double [casGetCalibrationFactors](#) (int ADevice, int AWhat, int AIndex, int AExtra)

Returns various details about the configuration/calibration of the given device.
- static void [casSetCalibrationFactors](#) (int ADevice, int AWhat, int AIndex, int AExtra, double AValue)

Changes various details about the configuration/calibration of the given device effectively overriding information that normally comes from configuration/calibration files.
- static void [casUpdateCalibrations](#) (int ADevice)

Updates the calibration information for the given device. Deprecated! Use [paUpdateSpectralCalibration](#).
- static void [casSaveCalibration](#) (int ADevice, string AFileName)

Saves the calibration of the given device to a file.
- static void [casClearCalibration](#) (int ADevice, int AWhat)

Clears the specified calibration part of the given device.
- static double [casGetData](#) (int ADevice, int AIndex)

Returns the intensity of the previously acquired spectrum for the given pixel.
- static double [casGetXArray](#) (int ADevice, int AIndex)

Returns the wavelength of the spectrum for the given pixel.
- static double [casGetDarkCurrent](#) (int ADevice, int AIndex)

Returns the intensity of the previously measured/calculated dark current for the given pixel.
- static void [casGetPhotInt](#) (int ADevice, out double APhotInt, StringBuilder AUnit, int AUnitMaxLen)

Retrieves the previously calculated photometric integral and it's data unit.
- static void [casGetRadInt](#) (int ADevice, out double ARadInt, StringBuilder AUnit, int AUnitMaxLen)

Retrieves the previously calculated radiometric integral and it's data unit.
- static double [casGetCentroid](#) (int ADevice)

Retrieves the previously calculated centroid wavelength.
- static void [casGetPeak](#) (int ADevice, out double x, out double y)

Retrieves the previously calculated peak wavelength and intensity.
- static double [casGetWidth](#) (int ADevice)

Retrieves the peak width in nm.
- static double [casGetWidthEx](#) (int ADevice, int AWhat)

Calculates and retrieves various aspects about the peak width (e.g. full width half maximum aka FWHM aka 50% bandwidth).
- static void [casGetColorCoordinates](#) (int ADevice, ref double x, ref double y, ref double z, ref double u, ref double v1976, ref double v1960)

Retrieves CIE color coordinates of a previously measured spectrum.
- static double [casGetCCT](#) (int ADevice)

Calculates the correlated color temperature CCT of a previously measured spectrum.
- static double [casGetCRI](#) (int ADevice, int Index)

Retrieves one of the previously calculated the color rendering indices.
- static void [casGetTriStimulus](#) (int ADevice, ref double X, ref double Y, ref double Z)

Retrieves the previously calculated X, Y and Z tristimulus values.
- static double [casGetExtendedColorValues](#) (int ADevice, int AWhat)

Retrieves results and conditions for previously performed colormetric calculations.
- static int [casColorMetric](#) (int ADevice)

Calculates colormetric results for the previously measured spectrum.
- static int [casCalculateCRI](#) (int ADevice)

Calculates the color rendering indices (CRI) of a previously measured spectrum.
- static int [cmXYToDominantWavelength](#) (double x, double y, double IllX, double IllY, ref double LambdaDom, ref double Purity)

Calculates dominant wavelength aka *LambdaDom* and purity from a given color coordinate and illuminant reference. This method is deprecated because it does not take the new [mpidColormetricType](#) into account. Use [casCalculateLambdaDom](#) instead.

- static int [casCalculateLambdaDom](#) (int ADevice, double IIX, double IIY, ref double LambdaDom, ref double Purity)
Calculates dominant wavelength aka LambdaDom and purity for the specified illuminant reference and the color coordinates that have been calculated by casColorMetric. This method replaces the now deprecated cmXYToDominantWavelength and does take the new mpidColormetricType into account.
- static IntPtr [casGetDLLFileName](#) (StringBuilder Dest, int AMaxLen)
Retrieves the complete path and file name of CAS library
- static IntPtr [casGetDLLVersionNumber](#) (StringBuilder Dest, int AMaxLen)
Retrieves the version of CAS library
- static int [casSaveSpectrum](#) (int ADevice, string AFileName)
Saves a previously measured spectrum to an .ISD file.
- static double [casGetExternalADCValue](#) (int ADevice, int AIndex)
Obsolete. Use mpidCurrentCCDTemperature. Used to retrieve the CCD temperature.
- static void [casSetStatusLED](#) (int ADevice, int AWhat)
Method to control the status LED of the spectrometer
- static int [casNmToPixel](#) (int ADevice, double nm)
Converts a wavelength into the corresponding CCD pixel index.
- static double [casPixelToNm](#) (int ADevice, int APixel)
Converts a CCD pixel index into the corresponding wavelength.
- static int [casCalculateTOPParameter](#) (int ADevice, int AAperature, double ADistance, ref double ASpotSize, ref double AFieldOfView)
Calculates spot size and field of view of the TOP.
- static int [casMultiTrackInit](#) (int ADevice, int ATracks)
Initializes the MultiTrack buffer
- static int [casMultiTrackDone](#) (int ADevice)
Releases a previously allocated MultiTrack buffer
- static int [casMultiTrackCount](#) (int ADevice)
Deprecated method.
- static int [casMultiTrackCopyData](#) (int ADevice, int ATrack)
Releases a previously allocated MultiTrack buffer
- static int [casMultiTrackSaveData](#) (int ADevice, string AFileName)
Saves a MultiTrack buffer to a .SWM file.
- static int [casMultiTrackLoadData](#) (int ADevice, string AFileName)
Loads a MultiTrack buffer from a .SWM file.
- static void [casLoadTestData](#) (int ADevice, string AFileName)
Loads a spectrum from an .ISD file into the CASID

Public Attributes

- const int [ErrorNoError](#) = 0
No error, i.e. function was successful.
- const int [ErrorUnknown](#) = -1
An unknown or unexpected error. Retrieve the corresponding error message via casGetErrorMessage for additional info.
- const int [ErrorTimeoutRWSNoData](#) = -2
A timeout occurred while waiting for the spectrum data.
- const int [ErrorInvalidDeviceType](#) = -3
Invalid InterfaceType constant passed to casCreateDeviceEx or casChangeDevice

- const int [ErrorAcquisition](#) = -4
Acquisition failed, i.e. an error occurred while the spectrometer was measuring.
- const int [ErrorAccuDataStream](#) = -5
Averaging failed, increase [mpidIntegrationTime](#) or make sure PC is not too busy.
- const int [ErrorPrivilege](#) = -6
no longer used
- const int [ErrorFIFOOverflow](#) = -7
It was detected that the FIFO of the interface card or the spectrometer was overflowing. Spectrum should be discarded.
- const int [ErrorTimeoutEOSScan](#) = -8
The spectrometer FW reported a timeout while waiting for EOS of a spectrum scan.
- const int [ErrorTimeoutEOSDummyScan](#) = -9
The spectrometer FW reported a timeout while waiting for EOS of a dummy scan.
- const int [ErrorFifoFull](#) = -10
It was detected that the FIFO of the interface card or the spectrometer was full. Spectrum should be discarded.
- const int [ErrorPixel1FinalCheck](#) = -11
After reading spectrum data, the FIFO did not contain a pixel marked as first pixel. Spectrum should be discarded.
- const int [ErrorCCDTemperatureFail](#) = -13
The temperature of the CCD is outside of the allowed range. Temperature control failed or ambient temperature too high.
- const int [ErrorAdrControl](#) = -14
The spectrometer hardware associated with the used Device/CASID could not be found or accessed. See chapter [Interfaces Types and Options](#).
- const int [ErrorFloat](#) = -15
Floating point error while calculating calibrated spectrum.
- const int [ErrorTriggerTimeout](#) = -16
Timeout reached while waiting for the trigger. See [mpidTriggerTimeout](#)
- const int [ErrorAbortWaitTrigger](#) = -17
While waiting for the trigger, the operation was aborted using [dpidAbortWaitForTrigger](#)
- const int [ErrorDarkArray](#) = -18
An error occurred while measuring the [Dark current array](#).
- const int [ErrorNoCalibration](#) = -19
Invalid calibration filename specified for [dpidCalibFileName](#)
- const int [ErrorInterfaceVersion](#) = -20
The firmware version of the interface card (PCI with CAS140 CTS) or the driver version is too old to work with this device.
- const int [ErrorCRI](#) = -21
Error calculating CRI (color rendering indices). For more info see [casCalculateCRI](#)
- const int [ErrorNoMultiTrack](#) = -25
Old way of doing MultiTrack not supported for this interface type. See chapter [MultiTrack measurements](#).
- const int [ErrorInvalidTrack](#) = -26
no longer used
- const int [ErrorDetectPixel](#) = -31
An error occurred while trying to detect the number of pixels the spectrometer has.
- const int [ErrorSelectParamSet](#) = -32
Error while activating a parameter set [dpidCurrentParamSet](#). See chapter [Working with Parameter Sets](#).
- const int [ErrorI2CInit](#) = -35
Initializing I2C bus failed.
- const int [ErrorI2CBusy](#) = -36
I2C operation failed because the I2C bus was busy.
- const int [ErrorI2CNotAck](#) = -37

- I2C operation was not acknowledged.*
- const int [ErrorI2CRelease](#) = -38
no longer used
- const int [ErrorI2CTimeOut](#) = -39
I2C operation timed out.
- const int [ErrorI2CTransmission](#) = -40
I2C transmission failed in an unspecified way.
- const int [ErrorI2CController](#) = -41
The I2C controller responded in an unexpected way.
- const int [ErrorDataNotAck](#) = -42
no longer used
- const int [ErrorNoExternalADC](#) = -52
The temperature ADC should be read, but the spectrometer does not have such an ADC.
- const int [ErrorShutterPos](#) = -53
Positioning the shutter or querying it's state failed. See chapter [Density Filter](#).
- const int [ErrorFilterPos](#) = -54
Positioning the filterwheel or querying it's position failed. See chapter [Density Filter](#).
- const int [ErrorConfigSerialMismatch](#) = -55
Configuration file not intended for this spectrometer and possible accessories. The serial number stored in the file [dpidConfigFileName](#) does not match the current [dpidSerialNo](#). Refer to chapter [Using Serial Numbers](#).
- const int [ErrorCalibSerialMismatch](#) = -56
Calibration file not intended for this spectrometer and possible accessories. The serial number stored in the file [dpidCalibFileName](#) does not match the current [dpidSerialNo](#). Refer to chapter [Using Serial Numbers](#).
- const int [ErrorInvalidParameter](#) = -57
Returned if there was an invalid "AWhat" parameter.
- const int [ErrorGetFilterPos](#) = -58
no longer used
- const int [ErrorParamOutOfRange](#) = -59
Returned by many methods indicating that a parameter was out of range; does not apply for ADevice parameter which results in [errCasDeviceNotFound](#)
- const int [ErrorDeviceFileChecksum](#) = -60
Returned by [dpidGetFilesFromDevice](#) if there was a checksum error, so either the device is corrupt or the communication failed.
- const int [ErrorInvalidEEPromType](#) = -61
Returned by [dpidGetFilesFromDevice](#) if the device does not support file storage.
- const int [ErrorDeviceFileTooLarge](#) = -62
Not enough storage on the device.
- const int [ErrorNoCommunication](#) = -63
The communication with the spectrometer ended unexpectedly.
- const int [ErrorNoFilesOnIdentKey](#) = -64
Returned by [dpidGetFilesFromDevice](#) if no files were found on the device.
- const int [ErrorExtraCalibFileInvalid](#) = -66
Additional files of the calibration are either missing or corrupt.
- const int [ErrorFeatureNotSupported](#) = -68
The requested feature is not supported by the device.
- const int [ErrorConfigUpToDate](#) = -70
Config/calib files do not match the ones on the device. Returned by [coCheckConfigUpToDate](#) and [casInitialize](#) if option [coCheckConfigUpToDate](#) is set.
- const int [ErrorCommunicationTimeout](#) = -73
The communication with the spectrometer reached a timeout.
- const int [ErrorTransmissionSerialMismatch](#) = -74

Transmission file not intended for this spectrometer and possible accessories. The serial number stored in the file [dpidTransmissionFileName](#) does not match the current [dpidSerialNo](#). Refer to chapter [Using Serial Numbers](#).

- const int [errCASOK](#) = [ErrorNoError](#)
just an alias for [ErrorNoError](#)
- const int [errCASError](#) = -1000
not returned; base for a few other errors like [errCasDeviceNotFound](#)
- const int [errCasNoConfig](#) = [errCASError](#)-3
Invalid filename specified for [dpidConfigFileName](#)
- const int [errCASDriverMissing](#) = [errCASError](#)-6
no longer used
- const int [errCasDeviceNotFound](#) = [errCASError](#)-10
Invalid ADevice / CASID parameter.
- const int [InterfaceISA](#) = 0
No longer used. Deprecated ISA interface constant.
- const int [InterfacePCI](#) = 1
PCI interface constant. For use with e.g. [casCreateDeviceEx](#). See chapter [Interfaces Types and Options](#).
- const int [InterfaceTest](#) = 3
Demo mode interface constant. For use with e.g. [casCreateDeviceEx](#). See chapter [Interfaces Types and Options](#).
- const int [InterfaceUSB](#) = 5
USB interface constant. For use with e.g. [casCreateDeviceEx](#). See chapter [Interfaces Types and Options](#).
- const int [InterfacePCle](#) = 10
PCle interface constant. For use with e.g. [casCreateDeviceEx](#). See chapter [Interfaces Types and Options](#).
- const int [InterfaceEthernet](#) = 11
Ethernet interface constant. For use with e.g. [casCreateDeviceEx](#). See chapter [Interfaces Types and Options](#).
- const int [aoAssignDevice](#) = 0
Only device specific properties are copied. This is rarely necessary to assign separately. Examples: current filter wheel state, flag whether data is present etc.
- const int [aoAssignParameters](#) = 1
Only parameters are assigned. This includes all measurement parameters but also configuration and calibration as well as dark current. Pretty much everything which is not related to the state of the hardware.
- const int [aoAssignComplete](#) = 2
Assigns everything, i.e. equal to calling [casAssignDeviceEx](#) with [aoAssignDevice](#) and then with [aoAssignParameters](#).
- const int [InitOnce](#) = 0
[casInitialize](#): load config and calib, and only initialize the hardware if not done before
- const int [InitForced](#) = 1
[casInitialize](#): load config and calib, and always initialize the hardware even if it was initialized before
- const int [InitNoHardware](#) = 2
[casInitialize](#): load config and calib, but do not initialize the hardware
- const int [dpidIntTimeMin](#) = 101
float: minimum integration time in ms as supported by the device; check when setting [mpidIntegrationTime](#)
- const int [dpidIntTimeMax](#) = 102
float: maximum integration time in ms as supported by the device; check when setting [mpidIntegrationTime](#)
- const int [dpidDeadPixels](#) = 103
number of dead pixels at the beginning of the pixel array; see [casGetData](#)
- const int [dpidVisiblePixels](#) = 104
number of visible pixels after dead pixels in the pixel array that form the spectrum; see [casGetData](#)
- const int [dpidPixels](#) = 105
the total number of pixels in the pixel array; defined by the configuration file and/or the device - should not be modified; see [casGetData](#)
- const int [dpidParamSets](#) = 106
total number of parameter sets defined for this device - see [Working with Parameter Sets](#); when set to 0, all are cleared but a new one is added immediately

- const int [dpidCurrentParamSet](#) = 107
0-based index of the currently active parameter set - see [Working with Parameter Sets](#)
- const int [dpidADCRange](#) = 108
the ADC level (in Counts) where saturation effects may start, i.e. if [mpidMaxADCValue](#) is larger than [dpidADCRange](#), the measurement was saturated
- const int [dpidADCBits](#) = 109
the bit resolution of the ADC - for information purposes only. To check for saturation, use [dpidADCRange](#)
- const int [dpidSerialNo](#) = 110
the complete serial number of spectrometer and accessories if applicable. See [Using Serial Numbers](#) for an overview
- const int [dpidTOPSerial](#) = 111
deprecated: the serial number of a TOP200 that is required by the current calibration. Rather use [dpidTOPSerialEx](#) and refer to [Using Serial Numbers](#) for an overview
- const int [dpidTransmissionFileName](#) = 112
Returns or sets the path of the transmission correction file (typical extension .isa); see [coUseTransmission](#). Note that an invalid filename will not raise an error. The file is loaded during [casInitialize](#) or when [paUpdateCompleteCalibration](#) is called.
- const int [dpidConfigFileName](#) = 113
Returns or sets the path of the currently used configuration file (extension .ini). This file is required and loaded during [casInitialize](#) or [errCasNoConfig](#) will occur.
- const int [dpidCalibFileName](#) = 114
Returns or sets the path of the currently used calibration file (extension .isc). This file is required and loaded during [casInitialize](#) or [ErrorNoCalibration](#) will occur.
- const int [dpidCalibrationUnit](#) = 115
the calibration unit as defined by the calibration file [dpidCalibFileName](#)
- const int [dpidAccessorySerial](#) = 116
serial number of the ID key connected to the spectrometer. See [Using Serial Numbers](#) for an overview
- const int [dpidTriggerCapabilities](#) = 118
Returns a bit-set describing the trigger capabilities of the spectrometer. [tcoCanTrigger](#) and others. For more info, see task [Triggered Measurements](#).
- const int [dpidAveragesMax](#) = 119
Returns maximum number of averages supported by the device/configuration. Use for validating [mpidAverages](#)
- const int [dpidFilterType](#) = 120
Returns the density filter type: 0 = none, 1 = manual filter, 2,3 = filter wheel. Refer to [Density Filter](#) for an overview.
- const int [dpidRelSaturationMin](#) = 123
Returns the minimum relative saturation in % which is still considered good for this device. Compare against [mpid↔RelSaturation](#)
- const int [dpidRelSaturationMax](#) = 124
Returns the maximum relative saturation in % which is still considered good for this device. Compare against [mpid↔RelSaturation](#)
- const int [dpidInterfaceVersion](#) = 125
Returns an optional version information for the device/configuration. Currently only the PCI-Interface supports an interface version, returning the CPLD version of the used PCI-Card (divide by 100, i.e. 123 for version 1.23), all other interface types return 0.
- const int [dpidTriggerDelayTimeMax](#) = 126
float: maximum trigger delay time in ms supported by the device/configuration. Use to check user input for the [mpidTriggerDelayTime](#) measurement parameter.
- const int [dpidSpectrometerName](#) = 127
Returns the name of the spectrometer type as specified in the configuration file.
- const int [dpidNeedDarkCurrent](#) = 130
Deprecated, rather check [dpidDCRemeasureReasons](#) and refer to chapter [dark current](#) for an overview.
- const int [dpidNeedDensityFilterChange](#) = 131
Returns whether the density filter needs to be moved (all values <> 0). Refer to chapter [density filter](#) for an overview.
- const int [dpidSpectrometerModel](#) = 132

- Returns an ID which identifies the spectrometer type as specified in the configuration file. Rather use [dpid↔SpectrometerName](#)
- const int [dpidLine1FlipFlop](#) = 133

Returns the state of the trigger flipflop on line 1 (CAS140B and CAS140CT only, refer to the hardware manual and the [tcoGetFlipState](#) trigger capability). The value is 0 for low state or 1 for high. When setting this value, the flipflop is reset, regardless of the value which is passed.
 - const int [dpidTimer](#) = 134

Returns the internal timer. Use not recommended. Rather use [mpidTimeSinceScanStart](#)
 - const int [dpidInterfaceType](#) = 135

Returns the interface type for the specified device (or a negative error). See chapter [Interface type and Options](#) for an overview.
 - const int [dpidInterfaceOption](#) = 136

Returns the interface option for the specified device (or a negative error). See chapter [Interface type and Options](#) for an overview.
 - const int [dpidInitialized](#) = 137

Returns whether the device has been correctly initialized (see [casInitialize](#)) or not. A value bigger than 0 indicates it has been initialized, 0 that it hasn't and a negative value indicates an error.
 - const int [dpidDCRemeasureReasons](#) = 138

Returns a set of flags which indicate why a dark current measurement is necessary. Values include [todcrrNeed↔DarkCurrent](#) and [todcrrCCDTemperature](#). See chapter [dark current](#) for an overview.
 - const int [dpidAbortWaitForTrigger](#) = 140

Not recommended - setting [dpidAbortWaitForTrigger](#) to non-zero will abort a triggered acquisition if the trigger has not yet occurred. Only supported by PCI interface.
 - const int [dpidGetFilesFromDevice](#) = 142

Set to the path of an existing directory and immediately a download of the files on the device will start. This might take considerable time, i.e. several seconds. The device does not have to be initialized, but interface type and option have to be set. Might return [ErrorInvalidEEPromType](#) or [ErrorDeviceFileChecksum](#)
 - const int [dpidTOPTType](#) = 143

Returns the type of TOP which is required by the current configuration. Possible values start at [ttNone](#)
 - const int [dpidTOPSerialEx](#) = 144

Returns the serial number of the TOP which is required by the current configuration or an empty string if no TOP is required. See [Using Serial Numbers](#).
 - const int [dpidAutoRangeFilterMin](#) = 145

Identifies the lowest density filter index for [mpidDensityFilter](#) that should be used for [AutoRange measurements](#) with the [coAutorangeFilter](#) option. This constant is part of the configuration file and should not be changed.
 - const int [dpidAutoRangeFilterMax](#) = 146

Identifies the highest density filter index for [mpidDensityFilter](#) that should be used for [AutoRange measurements](#) with the [coAutorangeFilter](#) option. This constant is part of the configuration file and should not be changed.
 - const int [dpidMultiTrackMaxCount](#) = 147

Returns the maximum number of tracks which are supported for [MultiTrack measurements](#) with this device.
 - const int [dpidSLCFileInfo](#) = 148

Returns file path and timestamp for a straylight calibration file, if present. Otherwise an empty string is returned.
 - const int [dpidCheckConfigFileSerial](#) = 149

Write-only string device parameter. Call with [dpidSerialNo](#) from another CASID to check whether this CASID uses the correct configuration file for this accessory.
 - const int [dpidCheckCalibFileSerial](#) = 150

Write-only string device parameter. Call with [dpidSerialNo](#) from another CASID to check whether this CASID uses the correct calibration file for this accessory.
 - const int [dpidExtraTransmissionsFileInfo](#) = 152

Returns file information about additional transmission files the calibration uses or an empty string if it doesn't.
 - const int [dpidMultiTrackCount](#) = 153

Returns the number of MultiTracks which have been initialized/allocated for [MultiTrack measurements](#) with this device.
 - const int [dpidIntTimePossibleResolutions](#) = 154

Returns a string containing the supported values for the [mpidIntTimeResolution](#), separated by semicolons. The string is formatted using the current system locale, times are in microseconds. Examples: "1000;100;25", but might also be "1000" if only 1ms is supported.

- const int [dpidCalibDate](#) = 155
Returns a string in ISO 8601 format containing the calibration date. Only valid after a call to [casInitialize](#).
- const int [dpidTriggerDelayTimeMin](#) = 159
float: minimum trigger delay time in ms supported by the device/configuration. Use to check user input for the [mpidTriggerDelayTime](#) measurement parameter.
- const int [dpidWavelengthCalibrationType](#) = 160
string: returns the type of wavelength calibration for the given CAS ID and it's INI/ISC. Only valid after [casInitialize](#) has been called successfully.
- const int [dpidCheckReferenceSpectrumFile](#) = 162
string, write-only: pass the complete file path of an ACS reference spectrum to check whether that spectrum has matching settings, so it can be compared with spectra of this CAS ID. Returns [ErrorInvalidCalibration](#) if the spectrum file could not be loaded, [ErrorInvalidParameter](#) if the settings of the spectrum do not match that of the CAS ID and [ErrorNoError](#) if they do match.
- const int [dpidFlashDurationMin](#) = 163
float, read-only: minimum supported [mpidFlashDuration](#) in milliseconds. Only applicable if [dpidTriggerCapabilities](#) contains [tcoFlashHardwareDuration](#).
- const int [dpidDebugLogFile](#) = 204
complete filename and path where debug log should be written to. Set to empty path to stop logging. This dpid does not require that you pass valid CASID when getting or setting this parameter (any value will be accepted).
- const int [dpidDebugLogLevel](#) = 205
integer level describing what to log, each higher level includes everything the levels below log, plus the things in this level. See constants starting with [DebugLogLevelErrors](#). This dpid does not require that you pass valid CASID when getting or setting this parameter (any value will be accepted).
- const int [dpidDebugMaxLogSize](#) = 206
integer file size, in bytes, the logfile should not exceed. Once this size is reached, a backup file with the same name but a .bak extension of the log is created and a new logfile starts. An already existing backup file will be deleted. This dpid does not require that you pass valid CASID when getting or setting this parameter (any value will be accepted).
- const int [tcoCanTrigger](#) = 0x00000001
bit for [dpidTriggerCapabilities](#). Spectrometer can be triggered externally, i.e. [mpidTriggerSource](#) supports [trgFlipFlop](#)
- const int [tcoTriggerDelay](#) = 0x00000002
bit for [dpidTriggerCapabilities](#). Spectrometer supports a trigger delay ([mpidTriggerDelayTime](#)).
- const int [tcoTriggerOnlyWhenReady](#) = 0x00000004
bit for [dpidTriggerCapabilities](#). Spectrometer supports changing the [toAcceptOnlyWhenReady](#) trigger option ([mpidTriggerOptions](#))
- const int [tcoAutoRangeTriggering](#) = 0x00000008
bit for [dpidTriggerCapabilities](#). Spectrometer supports the [toForEachAutoRangeTrial](#) trigger option ([mpidTriggerOptions](#))
- const int [tcoShowBusyState](#) = 0x00000010
bit for [dpidTriggerCapabilities](#). Spectrometer supports the [toShowBusyState](#) trigger option ([mpidTriggerOptions](#))
- const int [tcoShowACQState](#) = 0x00000020
bit for [dpidTriggerCapabilities](#). Spectrometer supports the [toShowACQState](#) trigger option ([mpidTriggerOptions](#))
- const int [tcoFlashOutput](#) = 0x00000040
bit for [dpidTriggerCapabilities](#). Spectrometer supports flash output in general, i.e. [mpidFlashType](#). Rather check [tcoFlashHardware](#) for [ftHardware](#) and [tcoFlashSoftware](#) for [ftSoftware](#)
- const int [tcoFlashHardware](#) = 0x00000080
bit for [dpidTriggerCapabilities](#). Spectrometer supports flash hardware, [ftHardware](#) for [mpidFlashType](#)
- const int [tcoFlashHardwareForEachAverage](#) = 0x00000100
bit for [dpidTriggerCapabilities](#). Spectrometer can output a flash signal for each averaged spectrum ([foEveryAverage](#) option of [mpidFlashOptions](#), [ftHardware](#))
- const int [tcoFlashSoftwareDelay](#) = 0x00000200
bit for [dpidTriggerCapabilities](#). Spectrometer supports flash delay ([mpidFlashDelayTime](#)), for [ftSoftware](#)

- const int [tcoFlashSoftwareDelayNegative](#) = 0x00000400
bit for [dpidTriggerCapabilities](#). Spectrometer supports a negative flash delay ([mpidFlashDelayTime](#)) for [ftSoftware](#)
- const int [tcoFlashSoftware](#) = 0x00000800
bit for [dpidTriggerCapabilities](#). Spectrometer supports flash type [ftSoftware](#) ([mpidFlashType](#))
- const int [tcoGetFlipFlopState](#) = 0x00001000
bit for [dpidTriggerCapabilities](#). Spectrometer supports reading the flip flop state using [dpidLine1FlipFlop](#)
- const int [tcoQueryHasData](#) = 0x00002000
bit for [dpidTriggerCapabilities](#). Spectrometer supports querying the FIFO state using [casFIFOHasData](#)
- const int [tcoACQStatePolarity](#) = 0x00004000
bit for [dpidTriggerCapabilities](#). Spectrometer supports changing the polarity of the ACQ state line via [mpidACQ↔StateLinePolarity](#)
- const int [tcoBusyStatePolarity](#) = 0x00008000
bit for [dpidTriggerCapabilities](#). Spectrometer supports changing the polarity of the Busy state line via [mpidBusy↔StateLinePolarity](#)
- const int [tcoFlashHardwareDelay](#) = 0x00010000
bit for [dpidTriggerCapabilities](#). Spectrometer supports flash delay ([mpidFlashDelayTime](#)), for [ftHardware](#)
- const int [tcoFlashHardwareDelayNegative](#) = 0x00020000
bit for [dpidTriggerCapabilities](#). Spectrometer supports a negative flash delay ([mpidFlashDelayTime](#)) for [ftHardware](#)
- const int [tcoFlashHardwareDuration](#) = 0x00040000
bit for [dpidTriggerCapabilities](#). Spectrometer supports an adjustable flash duration ([mpidFlashDuration](#)) for [ftHardware](#)
- const int [todcrrNeedDarkCurrent](#) = 0x0001
bit for [dpidDCRemeasureReasons](#). No valid dark current measurement present. This flag is identical to the now obsolete [dpidNeedDarkCurrent](#).
- const int [todcrrCCDTemperature](#) = 0x0002
bit for [dpidDCRemeasureReasons](#). A new dark current measurement should be done, since the CCD temperature has changed too much (varies depending on device; [mpidDCCCDTemperature](#) and [mpidLastCCDTemperature](#) are compared).
- const int [ttNone](#) = 0
[dpidTOPTType](#) constant for no TOP required
- const int [ttTOP100](#) = 1
[dpidTOPTType](#) constant for TOP 100 required
- const int [ttTOP200](#) = 2
[dpidTOPTType](#) constant for TOP 200 required
- const int [ttTOP150](#) = 3
[dpidTOPTType](#) constant for TOP 150 required
- const int [ttTOPLumiTOP](#) = 4
[dpidTOPTType](#) constant for LumiTOP required
- const int [DebugLogLevelErrors](#) = 1
[dpidDebugLogLevel](#) constant, only log errors
- const int [DebugLogLevelSaturation](#) = 2
[dpidDebugLogLevel](#) constant, only log errors and CCD saturation warnings
- const int [DebugLogLevelHardwareEvents](#) = 3
[dpidDebugLogLevel](#) constant, additionally log any hardware events
- const int [DebugLogLevelParameterChanges](#) = 4
[dpidDebugLogLevel](#) constant, additionally log any [mpid](#) oder [dpid](#) changes
- const int [DebugLogLevelAllMethodCalls](#) = 10
[dpidDebugLogLevel](#) constant, log every API call - warning, this may slow things down and create a big backlog of things which still need to be logged
- const int [casSerialComplete](#) = 0
[casGetSerialNumberEx](#) constant to retrieve complete serial number string, identical to [dpidSerialNo](#)
- const int [casSerialAccessory](#) = 1

- [casGetSerialNumberEx](#) constant to retrieve the serial number of additional equipment, if present, otherwise "N/A".; identical to [dpidAccessorySerial](#)
- const int [casSerialExtInfo](#) = 2

[casGetSerialNumberEx](#) constant to retrieve extended information, i.e. a multiline string with CAS type, ADC bits, serial number, No. of Pixels, ActivePixels and additional information depending on type and version of the firmware.
 - const int [casSerialDevice](#) = 3

[casGetSerialNumberEx](#) constant to retrieve the serial number read from the device; note that this might differ from case [casSerialComplete](#) without accessories or a TOP, since [casSerialComplete](#) falls back to the configuration file or its file name for a serial number (see [Using Serial Numbers](#) for more details)
 - const int [casSerialTOP](#) = 4

[casGetSerialNumberEx](#) constant to retrieve the serial number of a TOP if required by the calibration; identical to [dpidTOPSerialEx](#)
 - const int [coShutter](#) = 0x00000001

[casGetOptions](#) bit: device has a shutter; see [Dark Current](#)
 - const int [coFilter](#) = 0x00000002

[casGetOptions](#) bit: device has a density filter; see [Density Filter](#)
 - const int [coGetShutter](#) = 0x00000004

[casGetOptions](#) bit: device can check the position of the shutter (aka shutter control)
 - const int [coGetFilter](#) = 0x00000008

[casGetOptions](#) bit: device can check the position of the filter wheel (aka filter control)
 - const int [coGetAccessories](#) = 0x00000010

[casGetOptions](#) bit: device can retrieve serial number from additional equipment, see [dpidAccessorySerial](#)
 - const int [coGetTemperature](#) = 0x00000020

[casGetOptions](#) bit: device can measure CCD temperature, see [mpidCurrentCCDTemperature](#)
 - const int [coUseDarkcurrentArray](#) = 0x00000040

[casGetOptions](#) bit: option to enabled dark current array, see [Dark Current](#)
 - const int [coUseTransmission](#) = 0x00000080

[casGetOptions](#) bit: option to automatically apply a transmission correction when measuring as defined by [dpidTransmissionFileName](#)
 - const int [coAutorangeMeasurement](#) = 0x00000100

[casGetOptions](#) bit: option to automatically detect good integration time, i.e. perform a [AutoRange](#) measurements
 - const int [coAutorangeFilter](#) = 0x00000200

[casGetOptions](#) bit: option to use density filter during [AutoRange](#) measurements
 - const int [coCheckCalibConfigSerials](#) = 0x00000400

[casGetOptions](#) bit: option to automatically check if config and calibration files are intended for this device; see [Using Serial Numbers](#)
 - const int [coTOPHasFieldOfViewConfig](#) = 0x00000800

[casGetOptions](#) bit: indication that configuration/calibration includes field of view information for the TOP, so [mpidTOPFieldOfView](#) measurement condition is valid
 - const int [coAutoRemeasureDC](#) = 0x00001000

[casGetOptions](#) bit: option for turning automatic dark current measurements on/off. The interval after which the dark current is marked as invalid is defined by [mpidRemeasureDCInterval](#)
 - const int [coCanMultiTrack](#) = 0x00008000

[casGetOptions](#) bit: device supports [MultiTrack](#) measurements; depends on interface type
 - const int [coCanSwitchLEDOff](#) = 0x00010000

[casGetOptions](#) bit: device can switch off status LEDs during measurement (see option below)
 - const int [coLEDOffWhileMeasuring](#) = 0x00020000

[casGetOptions](#) bit: option for switching off status LEDs during measurements. Only supported if [coCanSwitchLEDOff](#) is present
 - const int [coCheckConfigUpToDate](#) = 0x00040000

[casGetOptions](#) bit: option to check whether the Ident-Key contains newer configuration/calibration files; happens during [casInitialize](#) and might result in [ErrorConfigUpToDate](#)
 - const int [coCanAverageOnDevice](#) = 0x00080000

- [casGetOptions](#) bit: device supports averaging on device, see [mpidAveragingOnDevice](#)
- const int [coCanMultiTrackSensorTemp](#) = 0x00100000

[casGetOptions](#) bit: device supports sensor temperature for MultiTrack, see [mpidCMTSensorTemp](#)
- const int [paPrepareMeasurement](#) = 1

[casPerformActionEx](#) ID: Prepares the spectrometer for the next measurement. This is intended to avoid wasting time within a measurement cycle and should be called before time-critical measurements are started. See chapter [Performing a Measurement](#) for other prerequisites.
- const int [paLoadCalibration](#) = 3

[casPerformActionEx](#) ID: Reloads the calibration file. Might be useful after certain calibration parts have been modified with [casClearCalibration](#) or [casSetCalibrationFactors](#).
- const int [paCheckAccessories](#) = 4

[casPerformActionEx](#) ID: Performs a check whether the current calibration and config files are suitable for the used spectrometer and optional optical probes equipped with an ID Key. The ID key is re-read during the method call. See chapter [Using Serial Numbers](#) for more info.
- const int [paMultiTrackStart](#) = 5

[casPerformActionEx](#) ID: Starts a MultiTrack measurement. The traditional way of starting MultiTrack with [casStart](#) is not supported by most interface types. See chapter [MultiTrack measurements](#) for an overview.
- const int [paAutoRangeFindParameter](#) = 7

[casPerformActionEx](#) ID: performs AutoRange measurements to find suitable measurement parameter. See chapter [AutoRange measurements](#) for details.
- const int [paCheckConfigUpToDate](#) = 12

[casPerformActionEx](#) ID: Performs a check whether the device has more recent configuration/calibration files. Returns [ErrorConfigUpToDate](#) if that is the case and files should be downloaded using [dpidGetFilesFromDevice](#)
- const int [paSearchForDevices](#) = 13

[casPerformActionEx](#) ID: Starts a search for devices (i.e. interface options) for one or all interface types (see below). For some interface types like USB, this search is synchronous, for others like Ethernet, the search is started asynchronously. See chapter [Interfaces Types and Options](#). This action does not require that you pass a valid CASID when calling [casPerformActionEx](#). If you want to search only for a specific interface type, pass the interface type as AParam1 to [casPerformActionEx](#) - which may result in [ErrorInvalidParameter](#) if the interface ID is not valid. If A↔Param2 is 1, some interface types, like USB, perform a full reset and search twice, which may fix problems in console apps where CAS might not appear, which first need a firmware upload.
- const int [paPrepareShutDown](#) = 14

[casPerformActionEx](#) ID: Stops all threads in the CAS library. Depending on how your app loads the library, it might be necessary to call this action from the main thread just before your app shuts down to avoid a dead-lock when the CAS library is unloaded. This action does not require that you pass a valid CASID when calling [casPerformActionEx](#).
- const int [paRecalcSpectrum](#) = 16

[casPerformActionEx](#) ID: Recalculates the spectrum from the raw data of the last measurement. This might be useful after you modified the calibration or something else that might influence the spectrum.
- const int [paUpdateSpectralCalibration](#) = 19

[casPerformActionEx](#) ID: Recalculates the spectral calibration for all parameter sets. A typical reason to call this method, is when you modified the transmission directly via [gcfTransmissionFunction](#). Whenever you change a measurement parameter that influences the spectral calibration (e.g. [mpidNewDensityFilter](#)), this action is called automatically.
- const int [paUpdateCompleteCalibration](#) = 20

[casPerformActionEx](#) ID: Recalculates the X vector for the spectrum and resamples the transmission to the updated wavelength before finally recalculating the spectral calibration. This action should be called whenever the wavelength calibration was modified or when the transmission should be reloaded.
- const int [paCheckAccessoriesTransmission](#) = 21

[casPerformActionEx](#) ID: Performs a check whether the current transmission file is suitable for the used spectrometer and optional optical probes equipped with an ID Key. The ID key is re-read during the method call. See chapter [Using Serial Numbers](#) for more info.
- const int [paMeasureParamSets](#) = 24

[casPerformActionEx](#) ID: Performs measurements for several parameter sets. AParam1 = FirstParamSet, AParam2 = LastParamSet. This action only returns after all measurements have been done or if an error occurred. See [Working with Parameter Sets](#).
- const int [paMeasureParamSetsDC](#) = 25

casPerformActionEx ID: Performs DC measurements for several parameter sets. AParam1 = FirstParamSet, AParam2 = LastParamSet. This action only returns after all DC measurements have been done or if an error occurred. See [Dark Current](#) and [Working with Parameter Sets](#).

- const int **mpidIntegrationTime** = 01

float: the integration time in ms which will be used for the next measurement for the currently active parameter set. Valid values range from [dpidIntTimeMin](#) to [dpidIntTimeMax](#). Note that after changing the integration time, a new dark current measurement might be necessary. See chapter [Dark Current](#) for an overview. Note that both [mpidIntTimeResolution](#) and [mpidIntTimeAlignPeriod](#) will affect the actual integration time that is used when measuring. It is therefore important to verify and document mpidIntegrationTime after the measurement.
- const int **mpidAverages** = 02

integer: the number of averages which will be used for the next measurement for the currently active parameter set. Valid values range from 1 to [dpidAveragesMax](#). Note that after changing the integration time, a new dark current measurement might be necessary. See chapter [Dark Current](#) for an overview.
- const int **mpidTriggerDelayTime** = 03

float: delay time in ms between trigger and start of the spectrum measurement. Valid values range from [dpidTriggerDelayTimeMin](#) to [dpidTriggerDelayTimeMax](#). Only applies to [Triggered Measurements](#). [dpidTriggerCapabilities](#) must include [tcoTriggerDelay](#) for this to work.
- const int **mpidTriggerTimeout** = 04

integer: timeout in ms which must not be exceeded when waiting for the external trigger or [ErrorTriggerTimeout](#) will be returned. Only applies to [Triggered Measurements](#).
- const int **mpidCheckStart** = 05

integer: first pixel which is included for checking the maximum value of the ADC during an acquisition, see [mpidMaxADCValue](#) and [mpidMaxADCPixel](#); leave at 0 to always use the first visible pixel of the calibration; see [dpidDeadPixels](#); use [casPixelToNm](#) and [casNmToPixel](#) for conversion between pixels and wavelength
- const int **mpidCheckStop** = 06

integer: last pixel which is included for checking the maximum value of the ADC during an acquisition, see [mpidMaxADCValue](#) and [mpidMaxADCPixel](#); leave at 0 to always use the last visible pixel of the calibration; see [dpidVisiblePixels](#); use [casPixelToNm](#) and [casNmToPixel](#) for conversion between pixels and wavelength
- const int **mpidColormetricStart** = 07

lower bound in nm of the spectral range which is used for the colormetric calculation. If it is smaller than the wavelength of the first visible pixel, then this wavelength is used (i.e. the colormetric range can never exceed the range of the spectrum). Default value is 380nm. If you want to use the full range, set mpidColormetricStart and mpidColormetricStop both to 0. See also [Getting the Spectrum and Results](#).
- const int **mpidColormetricStop** = 08

upper bound in nm of the spectral range which is used for the colormetric calculation. If it is bigger than the wavelength of the last visible pixel, then this wavelength is used (i.e. the colormetric range can never exceed the range of the spectrum). Default value is 780nm. If you want to use the full range, set mpidColormetricStart and mpidColormetricStop both to 0. See also [Getting the Spectrum and Results](#).
- const int **mpidACQTime** = 10

integer: returns the time in ms the last acquisition took. This includes integration time for all averages as well as the time to read out the FIFO. The accuracy of this timing depends on the interface type, with PCI being the most accurate.
- const int **mpidMaxADCValue** = 11

integer: returns the maximum value of the ADC during the last successful acquisition. Check this mpid for ADC-overflow by comparing it with the [dpidADCRange](#) device parameter. Do not calculate a relative signal level using mpidMaxADCValue yourself, but query [mpidRelSaturation](#) instead! [mpidMaxADCPixel](#) returns the pixel at which this maximum occurred. The pixel range which is taken into account can be customized by [mpidCheckStart](#) and [mpidCheckStop](#). See also the section about validating a measurement in [Performing a Measurement](#).
- const int **mpidMaxADCPixel** = 12

integer: returns the pixel which had the maximum ADC value [mpidMaxADCValue](#) during the last successful acquisition. The pixel range which is taken into account can be customized by [mpidCheckStart](#) and [mpidCheckStop](#). See also the section about validating a measurement in [Performing a Measurement](#).
- const int **mpidTriggerSource** = 14

integer: the current trigger source. Either [trgSoftware](#) for software triggering or [trgFlipFlop](#) when a hardware trigger is used. See chapter [Triggered Measurements](#) for an overview.
- const int **mpidAmpOffset** = 15

- integer*: returns the calculated amplifier offset of the last dark current measurement which has been performed with this device
- const int [mpidSkipLevel](#) = 16

intensity below which spectral intensities are not taken into account during colormetric calculations (the so called SkipLevel). This maybe useful to suppress noise from affecting the colormetric calculation. Note that skip-leveling needs to be enabled with [mpidSkipLevelEnabled](#).
 - const int [mpidSkipLevelEnabled](#) = 17

integer: determines whether skip-leveling is enabled (values <> 0). Set the actual intensity level using [mpidSkipLevel](#)
 - const int [mpidScanStartTime](#) = 18

integer: returns the value the internal timer had, when the last measurement was started. For [Triggered Measurements](#) the scan start is just before the spectrometer started accepting triggers, but not when the trigger actually arrived. See also [mpidTimeSinceScanStart](#) which might be useful for determining the timing at high-accuracy.
 - const int [mpidAutoRangeMaxIntTime](#) = 19

float: maximum integration time in milliseconds which must not be exceeded during an [AutoRange](#) measurements.
 - const int [mpidAutoRangeLevel](#) = 20

deprecated; use [mpidAutoRangeMinLevel](#) below
 - const int [mpidAutoRangeMinLevel](#) = 20

minimum relative maxADCValue between 0 and 100 percent which must be reached during an [AutoRange](#) measurements. See [mpidAutoRangeMaxLevel](#) for an upper limit.
 - const int [mpidDensityFilter](#) = 21

integer: returns previously set density filter or -1 if it never has been set. Use [mpidCurrentDensityFilter](#) to check the current physical filter position. When set, the filter wheel is immediately positioned. Valid values range from 0 to 7, otherwise [ErrorParamOutOfRange](#) will be returned. See chapter [Density Filter](#) for an overview and [Performing a Measurement](#) for prerequisites. [mpidDensityFilter](#) is not part of the [parameter sets](#), but [mpidNewDensityFilter](#) is. Use [casGetFilterName](#) to translate a filter index into a filter name.
 - const int [mpidCurrentDensityFilter](#) = 22

integer: returns the current physical position of the density filter wheel, if the spectrometer supports reading the filter position (see [coGetFilter](#) option). If the device has no filter (i.e. [coFilter](#) option missing), -1 will be returned, otherwise the filter index between 0 and 7.
 - const int [mpidNewDensityFilter](#) = 23

integer: the filter wheel position between 0 and 7 which should be used for the next measurement of the [active parameter set](#). See chapter [Density Filter](#) for an overview and [Performing a Measurement](#) for prerequisites. Reading [mpidNewDensityFilter](#) does not check the physical filter wheel position; use [mpidCurrentDensityFilter](#) for that.
 - const int [mpidLastDCAge](#) = 24

integer: returns the number of milliseconds that have passed, since the last [Dark Current](#) was measured for the [active parameter set](#) or -1 if no DC has been measured.
 - const int [mpidRelSaturation](#) = 25

returns the relative saturation (between 0 and 100%) of the previous successful measurement for the [active parameter set](#). [mpidRelSaturation](#) can be checked against the [dpidRelSaturationMin](#) and [dpidRelSaturationMax](#) range, to ensure that the measurement has a sufficient signal level. However, to check for oversaturated measurements use [mpidMaxADCValue](#) and check it against [dpidADCRange](#). Because with averaging, [mpidRelSaturation](#) might be < 100% even though there were saturated spectra!
 - const int [mpidPulseWidth](#) = 27

If this parameter is non-zero, it will cause the spectrum to be corrected by the ratio of [mpidPulseWidth](#) divided by [mpidIntegrationTime](#). This might be useful to correct spectra where the DUT is on for only a fraction of the integration time.
 - const int [mpidRemeasureDCInterval](#) = 28

integer: the time in ms, after which an automatic dark current measurement should be invalidated, so [dpidDCRemeasureReasons](#) becomes non-zero. This setting only has an effect, if the [coAutoRemeasureDC](#) option has been activated. The interval is checked against [mpidLastDCAge](#). See chapter [Dark Current](#) for an overview.
 - const int [mpidFlashDelayTime](#) = 29

float: the delay time of the flash output signal in milliseconds. Applies only when flash is activated, i.e. [mpidFlashType](#) <> [fitNone](#). If the delay time is negative, the flash signal and the delay happen before starting the measurement, otherwise the delay and the flash signal occur after the measurement started.
 - const int [mpidTOPAperture](#) = 30

integer: the TOP aperture which should be taken into account when the calibration is applied to the spectrum. The TOP aperture parameter is a 0-based index ranging from 0 to 6 corresponding to the TOP apertures 1 to 7. This parameter doesn't actually adjust the aperture of the TOP, but controls which aperture calibration factor will be applied.

- const int `mpidTOPDistance` = 31

the distance in mm from the DUT to the reference plane of the TOP. This parameter doesn't actually adjust the distance of the TOP, but controls which calibration factor will be applied.

- const int `mpidTOPSpotSize` = 32

returns the size of the measurement spot for the current `mpidTOPAperture` and `mpidTOPDistance`.

- const int `mpidTriggerOptions` = 33

integer: bit-set describing current trigger options. The set consists of the various `to<XXX>` bits, starting with `toAcceptOnlyWhenReady`.

- const int `mpidForceFilter` = 34

integer: flag which controls whether the filter wheel is moved for every measurement, even if it had been set to the same position previously. Internally the flag only affects `dpidNeedDensityFilterChange`, it will always return True if the ForceFilter flag is set.

- const int `mpidFlashType` = 35

integer: the type of the flash signal which should be emitted during or before a measurement. This signal is a short pulse during the integration time which is typically used to trigger flash lamps. To check whether the spectrometer supports a flash signal in general, check the `tcoFlashOutput` flag in `dpidTriggerCapabilities`. Use `mpidFlashOptions` to control further options related to the flash. Possible values for `mpidFlashType` are all `ft<XXX>` constants, starting at `ftNone`.

- const int `mpidFlashOptions` = 36

integer: bit-set describing the flash options of the device. Only applies, if flash is enabled, i.e. `mpidFlashType` is different from `ftNone`. The bits start with `foEveryAverage`.

- const int `mpidACQStateLine` = 37

integer: number of the digital out port which should be used for the `toShowACQState` option in `mpidTriggerOptions`. The port specified by this parameter is set to the level given by `mpidACQStateLinePolarity` after the spectrometer has been started and is ready to accept a trigger. The level is restored after the acquisition, i.e. after integration and read-out time. For possible port values refer to `casSetDigitalOut` and the hardware manual of the spectrometer.

- const int `mpidACQStateLinePolarity` = 38

integer: level of the ACQ state line. Only applies if the `toShowACQState` option is enabled in `mpidTriggerOptions`. The value is 0 for low level, all other values indicate high level. The port specified by `mpidACQStateLine` is set to this level after the spectrometer has been started and is ready to accept a trigger. The level is restored after the acquisition, i.e. after integration and read-out time.

- const int `mpidBusyStateLine` = 39

integer: number of the digital out port which should be used for the `toShowBusyState` option in `mpidTriggerOptions`. The port specified by this parameter is set to the level given by `mpidBusyStateLinePolarity` after the spectrometer has been started and is ready to accept a trigger. The level is restored after the trigger has been received. For possible port values refer to `casSetDigitalOut` and the hardware manual of the spectrometer.

- const int `mpidBusyStateLinePolarity` = 40

integer: level of the Busy state line. Only applies if the `toShowBusyState` option is enabled in `mpidTriggerOptions`. The value is 0 for low level, all other values indicate high level. The port specified by `mpidACQStateLine` is set to this level after the spectrometer has been started and is ready to accept a trigger. The level is restored after the acquisition, i.e. after integration and read-out time.

- const int `mpidAutoFlowTime` = 41

integer: returns the minimum flow time (in ms) the current source should switch on the DUT for, when it is triggered by either the busy or ACQ state line (`mpidTriggerOptions`). This time includes the `mpidIntegrationTime` for all `mpidAverages`, a `mpidTriggerDelayTime` as well as a flash delay, if it is negative (see `mpidFlashDelayTime`). The read-out time of the CCD is not included. See chapter [Synchronization](#), namely the section [Trigger CAS sequence](#)

- const int `mpidCRIMode` = 42

integer: controls the way the CRI is calculated (see `casCalculateCRI`). This is a global setting which affects all devices! Possible values are `criDIN6169`, which is the default, and `criCIE13_3_95`.

- const int `mpidObserver` = 43

integer: determines, which observer is used for `casColorMetric`. This is a global setting which affects all devices! Possible values are `2°cieObserver1931`, which is the default, and `10°cieObserver1964`.

- const int `mpidTOPFieldOfView` = 44

- returns the field of view for the current [mpidTOPAperture](#) and [mpidTOPDistance](#).
- const int [mpidCurrentCCDTemperature](#) = 46
 reading this [mpid](#) performs a temperature measurement and returns the CCD temperature in degrees Celsius. If the device does not support temperature measurements the return value will be NAN. To verify that the device supports temperature measurements, check the [coGetTemperature](#) flag in [casGetOptions](#).
 - const int [mpidLastCCDTemperature](#) = 47
 returns the previously measured CCD temperature in degrees Celsius. If the last temperature measurement is too old, a new temperature measurement will be performed. The interval which triggers a new temperature measurement is device dependent. If the device does not support temperature measurements the return value will be NAN. To verify that the device supports temperature measurements, check the [coGetTemperature](#) flag in [casGetOptions](#).
 - const int [mpidDCCCDTemperature](#) = 48
 returns the CCD temperature in degrees Celsius which was previously measured during the dark current measurement. If the dark current hasn't been measured or if the device does not support temperature measurements, the return value will be NAN. To verify that the device supports temperature measurements, check the [coGetTemperature](#) flag in [casGetOptions](#).
 - const int [mpidAutoRangeMaxLevel](#) = 49
 maximum relative [maxADCValue](#) between 0 and 100 percent which must not be exceeded during an [AutoRange](#) measurements. See [mpidAutoRangeMinLevel](#) for the lower limit.
 - const int [mpidMultiTrackAcqTime](#) = 50
 integer: returns the acquisition time in milliseconds of the complete [MultiTrack](#) measurement series. The accuracy of the timing varies depending on the spectrometer interface type.
 - const int [mpidTimeSinceScanStart](#) = 51
 integer: returns the time in milliseconds since the spectrometer was started, i.e. [mpidScanStartTime](#). This [mpid](#) can replace calls to [casStopTime](#), which has been deprecated.
 - const int [mpidCMTTrackStart](#) = 52
 float: returns the time in milliseconds between the start of the [MultiTrack](#) measurement series and the start of the current track (which was previously set with [casMultiTrackCopyData](#)). For the first track this is typically 0, the second it would be [mpidIntegrationTime](#), and so on.
 - const int [mpidColormetricWidthLevel](#) = 54
 integer: level in % which should be used for the next calculation of [casGetWidthEx](#). Default is 50% to calculate FWHM. Independent for each parameter set, see [Working with Parameter Sets](#)
 - const int [mpidIntTimeResolution](#) = 55
 float in microseconds, resolution for [mpidIntegrationTime](#). Supported values can be queried using [dpidIntTimePossibleResolutions](#) after the device has been initialized. Might return [ErrorParamOutOfRange](#) when set to an unsupported value. Attention: changing the [mpidIntTimeResolution](#) will affect [dpidIntTimeMax](#), so [mpidIntegrationTime](#) needs to be checked for being inside the allowed range.
 - const int [mpidIntTimeAlignPeriod](#) = 56
 float in ms: if non-zero positive, every measurement (also [AutoRange](#), but not DC array measurements) uses an integration time which is a multiple of this period. Check [mpidIntegrationTime](#) after the measurement for the actual integration time that was measured with.
 - const int [mpidColormetricType](#) = 57
 integer: determines the type of colormetric calculations that should be performed. Default is [cmtDefaultColormetric](#). Affects a wide range of calculation methods, e.g. [casColorMetric](#), [casGetCCT](#), [casCalculateCRI](#) and [casCalculateLambdaDom](#). Note: the wavelength range covered by the photopic luminosity function $V(\lambda)$ and color matching functions $x(\lambda)$, $y(\lambda)$, $z(\lambda)$ depend on [mpidColormetricType](#). For [cmtDefaultColormetric](#) the functions range from 380nm to 780nm and for [cmtSWPColormetric](#) they have a broader range of 360nm to 830nm. When using [cmtSWPColormetric](#) it is especially important to adjust [mpidColormetricStart](#) and [mpidColormetricStop](#), because their historic default values of 380nm and 780nm would unexpectedly constrain the calculation of the affected results. Set both [mpidColormetricStart](#) and [mpidColormetricStop](#) to zero in order to use the full range of the spectrum for calculation.
 - const int [mpidAveragingOnDevice](#) = 58
 integer: determines whether averaging happens on the device (values $\neq 0$) or in software, i.e. the CAS SDK. Only supported by some devices; check the [coCanAverageOnDevice](#) capability. If the device does not support averaging on device, getting and setting this [mpid](#) will return [ErrorFeatureNotSupported](#). Averaging on device might be helpful for a large number of averages with very short integration times, as only the averaged spectrum is transferred. It should however be taken into account that with this option enabled, the ADC non-linearity correction is only performed on the averaged spectrum and not on each measured spectrum before they are averaged.

- const int `mpidMeasured` = 59
integer (read-only): returns True (value <> 0) if the current parameter set was successfully measured. This flag is reset to False when a new measurement or DC measurement is started. When measuring single parameter sets, evaluating this flag is typically not necessary, as error checking the calls to `casMeasure` etc. is sufficient. When measuring multiple parameter sets with `paMeasureParamSets` or `paMeasureParamSetsDC`, all affected parameter sets get their measured flag reset immediately. This allows determining which spectra were correctly measured, even if an error occurred later on in the call.
- const int `mpidCMTSensorTemp` = 61
float (read-only): returns the sensor temperature in degree Celsius for the `MultiTrack` measurement series of the current track (which was previously set with `casMultiTrackCopyData`). Only devices returning `coCanMultiTrackSensorTemp` support this, others will just return 0.
- const int `mpidFlashDuration` = 63
float: the duration of the flash signal in milliseconds. Minimum supported value is `dpidFlashDurationMin`, but only if `dpidTriggerCapabilities` contains `tcoFlashHardwareDuration`
- const int `mpidColormetricPeakDiffWidth` = 67
float: the width in nm that is used for differentiation when calculating the peak wavelength `casGetPeak`. If negative, a width of 7nm or at least 6 pixels are used - like in previous versions. Default value is -1. Applies to the current CAS ID and the currently active ParamSet only!
- const int `cmtDefaultColormetric` = 0
`mpidColormetricType` value constant: default CAS SDK colormetric like in CAS SDK version 4.9 and before
- const int `cmtSWPColormetric` = 100
`mpidColormetricType` value constant: SWP colormetric like in SWP 3.6 and beyond. Note that SWP colormetric calculation might be taking more time than `cmtDefaultColormetric` due to the increased accuracy of the results.
- const int `toAcceptOnlyWhenReady` = 1
`mpidTriggerOptions` bit: the hardware trigger is only accepted when the spectrometer is ready for acquisition, i.e. the `dpidLine1FlipFlop` is reset before the spectrometer starts waiting for a trigger. This option flag can only be modified if `tcoTriggerOnlyWhenReady` is included in `dpidTriggerCapabilities`. Otherwise it is either enforced or cleared, depending on the spectrometer type.
- const int `toForEachAutoRangeTrial` = 2
`mpidTriggerOptions` bit: for hardware triggered `AutoRange` measurement this flag controls, whether each acquisition during `AutoRange` requires a hardware trigger or only the first one. Since the number of acquisitions during `AutoRange` can't be determined in advance, the spectrometer should be triggered whenever either the busy state or the ACQ state (see below) indicate that it is waiting for a trigger. Only supported if `tcoAutoRangeTriggering` flag is set in `dpidTriggerCapabilities`! Has no effect if `mpidTriggerSource` is not `trgFlipFlop`.
- const int `toShowBusyState` = 4
`mpidTriggerOptions` bit: enables the busy signal. As soon as the spectrometer is ready to accept a trigger, the signal on the digital out port identified by `mpidBusyStateLine` changes to the level specified by `mpidBusyStateLinePolarity`. The signal is reset as soon as the trigger was received. Only supported if `tcoShowBusyState` flag is set in `dpidTriggerCapabilities`!
- const int `toShowACQState` = 8
`mpidTriggerOptions` bit: enables the acquisition signal. The digital out port identified by `mpidACQStateLine` changes to the level given by `mpidACQStateLinePolarity` as soon as the spectrometer is waiting for a trigger and is reset only after the acquisition has been finished, i.e. after integration and read-out time. Only supported if `tcoShowACQState` flag is set in `dpidTriggerCapabilities`!
- const int `ftNone` = 0
`mpidFlashType` value constant: No flash signal is emitted. Default value and "supported" by all spectrometers.
- const int `ftHardware` = 1
`mpidFlashType` value constant: The flash signal is controlled by the spectrometer. Refer to the hardware manual for details. Depending on the flag `tcoFlashHardwareDelay` and related flags, this might enable the `mpidFlashDelayTime`. Only supported if the `tcoFlashHardware` flag is set in `dpidTriggerCapabilities`.
- const int `ftSoftware` = 2
`mpidFlashType` value constant: The flash signal is controlled by the CAS-DLL. Depending on the flag `tcoFlashSoftwareDelay` and related flags, this might enable the `mpidFlashDelayTime`. Software flash in general is only supported if the `tcoFlashSoftware` flag is set in `dpidTriggerCapabilities`.
- const int `foEveryAverage` = 1

- mpidFlashOptions* bit: if set, the flash signal is emitted for each acquisition when *mpidAverages* is bigger than 1. Only applies to hardware flash and only supported by some devices (check *tcoFlashHardwareForEachAverage* in *dpidTriggerCapabilities*).
- const int *trgSoftware* = 0
mpidTriggerSource value constant: no external trigger, *casMeasure* immediately starts the acquisition
 - const int *trgFlipFlop* = 3
mpidTriggerSource value constant: external trigger at PIN 1 (DB9) is used to start the acquisition. See chapter *Triggered Measurements* for a detailed overview on this topic.
 - const int *criDIN6169* = 0
mpidCRIMode value constant: CRI calculation according to DIN 6169
 - const int *criCIE13_3_95* = 1
mpidCRIMode value constant: CRI calculation according to CIE 13.3-95
 - const int *cieObserver1931* = 0
mpidObserver value constant: 2° observer according to CIE 1931
 - const int *cieObserver1964* = 1
mpidObserver value constant: 10° observer according to CIE 1964
 - const int *casShutterInvalid* = -1
return value of *casGetShutter* if current shutter state is invalid (device failure) or device has no shutter (check using *casGetOptions*)
 - const int *casShutterOpen* = 0
return value of *casGetShutter* if current shutter state is open, i.e. normal spectra will be measured
 - const int *casShutterClose* = 1
return value of *casGetShutter* if current shutter state is closed, i.e. dark current will be measured
 - const int *gcfDensityFunction* = 0
casGetCalibrationFactors AWhat: spectral calibration of the density filters
 - const int *gcfSensitivityFunction* = 1
casGetCalibrationFactors AWhat: spectral calibration
 - const int *gcfTransmissionFunction* = 2
casGetCalibrationFactors AWhat: transmission of additional optical hardware
 - const int *gcfDensityFactor* = 3
casGetCalibrationFactors AWhat: density filter factors
 - const int *gcfTOPApertureFactor* = 4
casGetCalibrationFactors AWhat: absolute calibration for TOP
 - const int *gcfTOPDistanceFactor* = 5
casGetCalibrationFactors AWhat: distance calibration factors for TOP
 - const int *gcfTDDistanceMin* = -3
casGetCalibrationFactors AIndex with AWhat=*gcfTOPDistanceFactor* to retrieve the minimum calibrated TOP distance
 - const int *gcfTDDistanceMax* = -2
casGetCalibrationFactors AIndex with AWhat=*gcfTOPDistanceFactor* to retrieve the maximum calibrated TOP distance
 - const int *gcfTDCount* = -1
casGetCalibrationFactors AIndex with AWhat=*gcfTOPDistanceFactor* to retrieve the number of TOP distance factors
 - const int *gcfTDExtraDistance* = 1
casGetCalibrationFactors AExtra with AWhat=*gcfTOPDistanceFactor* to retrieve a given TOP distance at AIndex
 - const int *gcfTDExtraFactor* = 2
casGetCalibrationFactors AExtra with AWhat=*gcfTOPDistanceFactor* to retrieve a given TOP distance factor at AIndex
 - const int *gcfWLCalibrationChannel* = 6
casGetCalibrationFactors AWhat: wavelength calibration channels, get/set with 0-based AIndex to get/set channels, i.e. Pixels with calibration points
 - const int *gcfWLExtraFirstCoefficient* = -6

- [*casGetCalibrationFactors*](#) AExtra with AWhat=[*gcfWLCalibrationChannel*](#) to retrieve the first coefficient C0 of a polynom wavelength calibration; check for non-zero to find out whether polynom calibration is used.
- const int [*gcfWLExtraLastCoefficient*](#) = -2

[*casGetCalibrationFactors*](#) AExtra with AWhat=[*gcfWLCalibrationChannel*](#) to retrieve the last coefficient C4 of a polynom wavelength calibration
- const int [*gcfWLCalibPointCount*](#) = -1

[*casGetCalibrationFactors*](#) AExtra with AWhat=[*gcfWLCalibrationChannel*](#) to retrieve the number of wavelength calibration points
- const int [*gcfWLExtraCalibrationDelete*](#) = 1

[*casSetCalibrationFactors*](#) AExtra with AWhat=[*gcfWLCalibrationChannel*](#): set to delete WL calibration point at AIndex
- const int [*gcfWLExtraCalibrationDeleteAll*](#) = 2

[*casSetCalibrationFactors*](#) AExtra with AWhat=[*gcfWLCalibrationChannel*](#): set to delete ALL calibration points
- const int [*gcfWLCalibrationAlias*](#) = 7

[*casGetCalibrationFactors*](#) AWhat: wavelength calibration wavelengths, get/set with 0-based AIndex to get/set wavelengths for corresponding channels/pixels
- const int [*gcfWLCalibrationSave*](#) = 8

[*casGetCalibrationFactors*](#) AWhat: set to any value, to update configuration file with current wavelength calibration
- const int [*gcfDarkArrayValues*](#) = 9

[*casGetCalibrationFactors*](#) AWhat: get/set integration times and actual dark current spectra of dark current array
- const int [*gcfDarkArrayDepth*](#) = -1

[*casGetCalibrationFactors*](#) AExtra with AWhat=[*gcfDarkArrayValues*](#): get/set depth/length of dark current array
- const int [*gcfDarkArrayIntTime*](#) = -2

[*casGetCalibrationFactors*](#) AExtra with AWhat=[*gcfDarkArrayValues*](#): get/set integration time of dark current array at AIndex
- const int [*gcfTOPParameter*](#) = 11

[*casGetCalibrationFactors*](#) AWhat: get various TOP parameter
- const int [*gcfTOPApertureSize*](#) = 0

[*casGetCalibrationFactors*](#) AExtra with AWhat=[*gcfTOPParameter*](#): get aperture size of TOP aperture AIndex
- const int [*gcfTOPSpotSizeDenominator*](#) = 1

[*casGetCalibrationFactors*](#) AExtra with AWhat=[*gcfTOPParameter*](#): get spot size distance denominator to calculate spot size from distance
- const int [*gcfTOPSpotSizeOffset*](#) = 2

[*casGetCalibrationFactors*](#) AExtra with AWhat=[*gcfTOPParameter*](#): get spot size distance offset to calculate spot size from distance
- const int [*gcfLinearityFunction*](#) = 12

[*casGetCalibrationFactors*](#) AWhat: linearity calibration of the ADC
- const int [*gcfLinearityCounts*](#) = 0

[*casGetCalibrationFactors*](#) AExtra with AWhat=[*casGetCalibrationFactors*](#): get counts of linearity calibration array at AIndex
- const int [*gcfLinearityFactor*](#) = 1

[*casGetCalibrationFactors*](#) AExtra with AWhat=[*casGetCalibrationFactors*](#): get factor of linearity calibration array at AIndex
- const int [*gcfRawData*](#) = 14

[*casGetCalibrationFactors*](#) AWhat: retrieve raw data of previously acquired spectrum, AIndex = channel/pixel
- const int [*cLambdaWidth*](#) = 0

AWhat constant for [*casGetWidthEx*](#): returns the peak width in nm, identical to [*casGetWidth*](#)
- const int [*cLambdaLow*](#) = 1

AWhat constant for [*casGetWidthEx*](#): start wavelength of the width.
- const int [*cLambdaMiddle*](#) = 2

AWhat constant for [*casGetWidthEx*](#): center wavelength of the width, i.e. Center - Start = Stop - Center.
- const int [*cLambdaHigh*](#) = 3

AWhat constant for [*casGetWidthEx*](#): stop wavelength of the width.
- const int [*cLambdaOuterWidth*](#) = 4

- AWhat constant for [casGetWidthEx](#): identical to cLambdaWidth, but using old algorithm.*

 - const int [cLambdaOuterLow](#) = 5

AWhat constant for [casGetWidthEx](#): identical to cLambdaLow, but using old algorithm.

 - const int [cLambdaOuterMiddle](#) = 6

AWhat constant for [casGetWidthEx](#): identical to cLambdaMiddle, but using old algorithm.

 - const int [cLambdaOuterHigh](#) = 7

AWhat constant for [casGetWidthEx](#): identical to cLambdaHigh, but using old algorithm.

 - const int [ecvVisualEffect](#) = 2

AWhat constant for [casGetExtendedColorValues](#): Visual Effect in %. Tristimulus Y divided by radiometric VIS integral (380..780nm)

 - const int [ecvUVA](#) = 3

AWhat constant for [casGetExtendedColorValues](#): UVA radiometric integral 315..400 nm.

 - const int [ecvUVB](#) = 4

AWhat constant for [casGetExtendedColorValues](#): UVB radiometric integral 280..315 nm.

 - const int [ecvUVC](#) = 5

AWhat constant for [casGetExtendedColorValues](#): UVC radiometric integral 200..280 nm.

 - const int [ecvVIS](#) = 6

AWhat constant for [casGetExtendedColorValues](#): VIS radiometric integral 380..780 nm.

 - const int [ecvCRICCT](#) = 7

AWhat constant for [casGetExtendedColorValues](#): correlated color temperature in K used for CRI calculation. Not identical to [casGetCCT](#)!

 - const int [ecvCDI](#) = 8

AWhat constant for [casGetExtendedColorValues](#): color discrimination index.

 - const int [ecvDistance](#) = 9

AWhat constant for [casGetExtendedColorValues](#): distance between the color locus (u, v) and the Planck curve.

 - const int [ecvCalibMin](#) = 10

AWhat constant for [casGetExtendedColorValues](#): minimum wavelength of the currently used calibration.

 - const int [ecvCalibMax](#) = 11

AWhat constant for [casGetExtendedColorValues](#): maximum wavelength of the currently used calibration.

 - const int [ecvScotopicInt](#) = 12

AWhat constant for [casGetExtendedColorValues](#): scotopic integral, unit identical to photometric integral retrieved by [casGetPhotInt](#)

 - const int [ecvCRIFirst](#) = 100

AWhat constant for [casGetExtendedColorValues](#): retrieves first (common) CRI, increase up to [ecvCRILast](#) to retrieve others.

 - const int [ecvCRILast](#) = 116

AWhat constant for [casGetExtendedColorValues](#): retrieves last CRI, see [ecvCRIFirst](#).

 - const int [ecvCRITriKXFirst](#) = 120

AWhat constant for [casGetExtendedColorValues](#): intermediate result of CRI calculation. See overview table at [cas↔GetExtendedColorValues](#).

 - const int [ecvCRITriKXLast](#) = 136

AWhat constant for [casGetExtendedColorValues](#): intermediate result of CRI calculation. See overview table at [cas↔GetExtendedColorValues](#).

 - const int [ecvCRITriKYFirst](#) = 140

AWhat constant for [casGetExtendedColorValues](#): intermediate result of CRI calculation. See overview table at [cas↔GetExtendedColorValues](#).

 - const int [ecvCRITriKYLast](#) = 156

AWhat constant for [casGetExtendedColorValues](#): intermediate result of CRI calculation. See overview table at [cas↔GetExtendedColorValues](#).

 - const int [ecvCRITriKZFirst](#) = 160

AWhat constant for [casGetExtendedColorValues](#): intermediate result of CRI calculation. See overview table at [cas↔GetExtendedColorValues](#).

 - const int [ecvCRITriKZLast](#) = 176

- AWhat constant for [casGetExtendedColorValues](#): intermediate result of CRI calculation. See overview table at [cas↔GetExtendedColorValues](#).*

 - const int [ecvCRITriRXordUFirst](#) = 180

AWhat constant for [casGetExtendedColorValues](#): intermediate result of CRI calculation. See overview table at [cas↔GetExtendedColorValues](#).

 - const int [ecvCRITriRXordULast](#) = 196

AWhat constant for [casGetExtendedColorValues](#): intermediate result of CRI calculation. See overview table at [cas↔GetExtendedColorValues](#).

 - const int [ecvCRITriRYordVFirst](#) = 200

AWhat constant for [casGetExtendedColorValues](#): intermediate result of CRI calculation. See overview table at [cas↔GetExtendedColorValues](#).

 - const int [ecvCRITriRYordVLast](#) = 216

AWhat constant for [casGetExtendedColorValues](#): intermediate result of CRI calculation. See overview table at [cas↔GetExtendedColorValues](#).

 - const int [ecvCRITriRZordWFirst](#) = 220

AWhat constant for [casGetExtendedColorValues](#): intermediate result of CRI calculation. See overview table at [cas↔GetExtendedColorValues](#).

 - const int [ecvCRITriRZordWLast](#) = 236

AWhat constant for [casGetExtendedColorValues](#): intermediate result of CRI calculation. See overview table at [cas↔GetExtendedColorValues](#).

 - const int [extNoError](#) = 0

AWhat constant for [casSetStatusLED](#)

 - const int [extExternalError](#) = 1

AWhat constant for [casSetStatusLED](#)

 - const int [extFilterBlink](#) = 2

AWhat constant for [casSetStatusLED](#)

 - const int [extShutterBlink](#) = 4

AWhat constant for [casSetStatusLED](#)

8.1.1 Detailed Description

The CAS DLL interface class

This static class provides access to all methods the CAS DLL exports as well as defining all the constants that are used for calling the methods or interpreting the results. Even if you're accessing the CAS DLL directly and not via this interface class, you can still use this documentation as a reference for the available methods and constants.

8.1.2 Member Function Documentation

8.1.2.1 [casAssignDeviceEx\(\)](#)

```
static int InstrumentSystems.CAS4.CAS4DLL.casAssignDeviceEx (
    int ASourceDevice,
    int ADestDevice,
    int AOption ) [static]
```

Assigns properties and or parameters from one device handle to another

Parameters

<i>ASourceDevice</i>	Source device
<i>ADestDevice</i>	Destination device
<i>AOption</i>	Controls which aspects get assigned from ASourceDevice to ADestDevice

Returns

This method returns 0 if the assignment was successful. Otherwise it returns an Error Code

Assigning devices might be useful if you want to store different measurement setups for a given device. By keeping copies of all parameters with a second device handle it is possible to store measurement setups which may even differ by calibration or dark current, something which cannot be achieved using parameter sets.

8.1.2.2 casCalculateCorrectedData()

```
static void InstrumentSystems.CAS4.CAS4DLL.casCalculateCorrectedData (
    int ADevice ) [static]
```

Applies the spectral correction to the previously acquired raw spectrum.

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

This method causes the spectral correction to be applied to the previously acquired raw spectrum but without a transmission correction defined by [dpidTransmissionFileName](#)! This is done automatically after a measurement (but taking the transmission correction into account if applicable), so normally there is no need to call this method explicitly.

If you want to apply the transmission correction as well, use [casConvoluteTransmission](#) instead.

Note

This method call is useful if for one spectrum, you want the same spectrum with and without the transmission correction applied.

8.1.2.3 casCalculateCRI()

```
static int InstrumentSystems.CAS4.CAS4DLL.casCalculateCRI (
    int ADevice ) [static]
```

Calculates the color rendering indices (CRI) of a previously measured spectrum.

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

The return value is 0, if the CRI calculation was successful or a negative error code. Most likely is [ErrorCRI](#), which can indicate that there's not enough signal to calculate CRI or that the CCT is outside a supported range.

Since different standards for calculating CRI are supported, use the global [mpidCRIMode](#) to select the standard before calling `casCalculateCRI`.

Note

[casGetCCT](#) must have been called before calling [casCalculateCRI](#), because the CRI calculation uses a special CCT ([ecvCRICCT](#)).

To retrieve the actual color rendering indices, use [casGetCRI](#).

[Delphi]

```
...
casColorMetric(CASID);
CheckError(casGetError(CASID));
CCT:= casGetCCT(CASID); //also calculates CRICCT and stores it internally
CheckError(casGetError(CASID));
CheckError(casCalculateCRI(CASID));
CRIRa:= casGetCRI(CASID, 0); //Index = 0 retrieves common CRI
CheckError(casGetError(CASID));
...
```

8.1.2.4 casCalculateLambdaDom()

```
static int InstrumentSystems.CAS4.CAS4DLL.casCalculateLambdaDom (
    int ADevice,
    double IIIX,
    double IIY,
    ref double LambdaDom,
    ref double Purity ) [static]
```

Calculates dominant wavelength aka LambdaDom and purity for the specified illuminant reference and the color coordinates that have been calculated by [casColorMetric](#). This method replaces the now deprecated [cmXYToDominantWavelength](#) and does take the new [mpidColormetricType](#) into account.

Parameters

<i>ADevice</i>	The device / CASID
<i>IIIX</i>	The x color coordinate of the illuminant reference
<i>IIY</i>	The y color coordinate of the illuminant reference
<i>LambdaDom</i>	The variable which will receive the calculated dominant wavelength
<i>Purity</i>	The variable which will receive the calculated purity

Returns

The return value is 0, if the calculation was successful, or a negative error code (see chapter [Error Handling](#)).

The device must have a valid spectrum and [casColorMetric](#) must have been called successfully. IIIX and IIY are the color coordinates of the illuminant reference. One of the typical references is Illuminant E where x and y are 0.3333. Purity is 0 if x = IIIX and y = IIY. Purity will never exceed 1.

8.1.2.5 casCalculateTOPParameter()

```
static int InstrumentSystems.CAS4.CAS4DLL.casCalculateTOPParameter (
    int ADevice,
    int AAperture,
```



```
double ADistance,
ref double ASpotSize,
ref double AFieldOfView ) [static]
```

Calculates spot size and field of view of the TOP.

Parameters

<i>ADevice</i>	The device / CASID
<i>AAperture</i>	The 0-based index of the TOP aperture, e.g. mpidTOPAperture
<i>ADistance</i>	The TOP distance in mm, e.g. mpidTOPDistance
<i>ASpotSize</i>	Variable which will receive the calculated spot size in mm.
<i>AFieldOfView</i>	Variable which will receive the calculated field of view in °.

Returns

The return value is 0 if the method was successful (even if [coTOPHasFieldOfViewConfig](#) is not present!). Otherwise the return value is an error code.

This method calculates the spot size and field of view of the TOP for the given aperture and distance and returns them in the *ASpotSize* and *AFieldOfView* parameter.

```
//this demo displays the FieldOfView in ° and SpotSize in mm for the current TOP aperture
...
var
    lSpotSize, lFieldOfView: Double;
...
if casGetOptions(CASID) and coTOPHasFieldOfViewConfig <> 0 then
begin
    CheckError(casCalculateTOPParameter(CASID, CurrentTOPAperture, CurrentTOPDistance, lSpotSize,
        lFieldOfView));
    lblFieldOfView.Caption:= FormatFloat('0.## °', lFieldOfView);
    lblSpotSize.Caption:= FormatFloat('0.## mm', lSpotSize);
end else
begin
    lblFieldOfView.Caption:= 'No Field Of View Info';
    lblSpotSize.Caption:= 'N/A';
end;
```

8.1.2.6 casChangeDevice()

```
static int InstrumentSystems.CAS4.CAS4DLL.casChangeDevice (
    int ADevice,
    int AInterfaceType,
    int AInterfaceOption ) [static]
```

Change the interface type and/or option of a device / CASID

Parameters

<i>ADevice</i>	The device / CASID
<i>AInterfaceType</i>	New interface type
<i>AInterfaceOption</i>	New interface option

Returns

The CASID of the changed device or a negative error like [ErrorInvalidDeviceType](#) or [errCasDeviceNotFound](#)

Using this method it is possible to change the interface type and option of a device created with [casCreateDeviceEx](#). Refer to the chapter [Interfaces Types and Options](#) for an overview about interface types and options.

8.1.2.7 casClearCalibration()

```
static void InstrumentSystems.CAS4.CAS4DLL.casClearCalibration (
    int ADevice,
    int AWhat ) [static]
```

Clears the specified calibration part of the given device.

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	One of the gcf<XXX> constants, specifying which calibration part to delete. See table below for possible values.

The following table lists the possible values for AWhat

AWhat	Description
gcfDensityFunction	Spectral calibration of all density filters
gcfSensitivityFunction	Spectral calibration of the spectrometer
gcfTransmissionFunction	Transmission correction of additional optical accessories
gcfTOPApertureFactor	Absolute calibration of the TOP apertures
gcfTOPDistanceFactor	TOP distance factors
gcfWLCalibrationChannel	Clears the wavelength calibration
gcfLinearityFunction	Linearity calibration of the ADC

There is no need to call [casUpdateCalibrations](#) after [casClearCalibration](#), as this is done implicitly.

8.1.2.8 casClearDarkCurrent()

```
static int InstrumentSystems.CAS4.CAS4DLL.casClearDarkCurrent (
    int ADevice ) [static]
```

Clears a previously measured dark current of the given device

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

The return value is 0 if no error occurred, otherwise it is one of the error codes.

Normally there is no need to call this method since a new dark current measurement (see [casMeasureDarkCurrent](#)) automatically replaces the previous one. This method only affects the dark current of the current parameter set (see [Working with Parameter Sets](#)). Refer to chapter [Dark Current](#) for an overview over related methods and dpid's.

8.1.2.9 casColorMetric()

```
static int InstrumentSystems.CAS4.CAS4DLL.casColorMetric (
    int ADevice ) [static]
```

Calculates colormetric results for the previously measured spectrum.

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

The return value is 0, if the method was successful, or a negative error code (see chapter [Error Handling](#)).

The spectral range taken into account for the calculation is defined by [mpidColormetricStart](#) and [mpidColormetricStop](#). The calculation is also influenced by [mpidObserver](#). Refer to the chapter [Getting the Spectrum and Results](#) for an overview of the methods (and the order they have to be called) to retrieve the results of this calculation.

8.1.2.10 casConvoluteTransmission()

```
static void InstrumentSystems.CAS4.CAS4DLL.casConvoluteTransmission (
    int ADevice ) [static]
```

Applies the spectral and transmission correction to the previously acquired raw spectrum.

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

This method causes the spectral and transmission correction to be applied to the previously acquired raw spectrum! This is done automatically after a measurement if the [coUseTransmission](#) option is enabled, so normally there is no need to call this method explicitly. The file containing the transmission correction must have been set using [dpidTransmissionFileName](#) before the device was initialized with [casInitialize](#)!

If you don't want to apply the transmission correction as well, use [casCalculateCorrectedData](#) instead.

Note

This method call is useful if for one spectrum, you want the same spectrum with and without the transmission correction applied.

8.1.2.11 casCreateDevice()

```
static int InstrumentSystems.CAS4.CAS4DLL.casCreateDevice ( ) [static]
```

Deprecated method. Creates a device context aka CASID within the CAS DLL

Returns

If the function was successful, the return value is a device handle (≥ 0). A negative value indicates an error

Use [casCreateDeviceEx](#) instead

8.1.2.12 casCreateDeviceEx()

```
static int InstrumentSystems.CAS4.CAS4DLL.casCreateDeviceEx (
    int AInterfaceType,
    int AInterfaceOption ) [static]
```

Creates a device context aka CASID within the CAS DLL

Parameters

<i>AInterfaceType</i>	Type of the interface, i.e. the way the spectrometer is connected to the PC
<i>AInterfaceOption</i>	The interface option, identifying either the spectrometer itself or the interface card

Returns

If the function was successful, the return value is a device handle (≥ 0). A negative value indicates an error, e.g. -3 = `ErrorInvalidDeviceType`

This method must be called to work with a spectrometer. It returns a device handle which is later used with subsequent function calls to identify a specific device (typically called the `ADevice` parameter, also referred to as CASID).

For every device handle created, you should call [casDoneDevice](#) once it is no longer needed.

For interface type and option, there are some predefined interface constants (like [InterfaceUSB](#)), but it is generally recommended to enumerate interface types and options and let the user choose from a set of two lists. Refer to chapter [Interfaces Types and Options](#) for a detailed overview.

Note

`casCreateDeviceEx` does not perform hardware initialization. Use [casInitialize](#) for that.

8.1.2.13 casDeleteParamSet()

```
static int InstrumentSystems.CAS4.CAS4DLL.casDeleteParamSet (
    int ADevice,
    int AParamSet ) [static]
```

Deletes a specific parameter set for the given device.

Parameters

<i>ADevice</i>	The device / CASID
<i>AParamSet</i>	The 0-based index of the parameter set which should be deleted

Returns

The return value is 0 if no error occurred, otherwise it is one of the error codes, e.g. [ErrorSelectParamSet](#)

AParamSet may range from 0 to [dpidParamSets](#) - 1. If the [dpidCurrentParamSet](#) is deleted, the first parameter set becomes active. There must be at least one parameter set, so [casDeleteParamSet](#) cannot delete the last remaining one. Refer to chapter [Working with Parameter Sets](#) for an overview.

8.1.2.14 casDoneDevice()

```
static int InstrumentSystems.CAS4.CAS4DLL.casDoneDevice (
    int ADevice ) [static]
```

Release resources used by the device

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

The method returns 0 if finalization of the device was successful, otherwise an error code

Call this method for every device which was created with [casCreateDeviceEx](#). Note that after a successful call to this method, the CASID ADevice is no longer valid!

8.1.2.15 casFIFOHasData()

```
static int InstrumentSystems.CAS4.CAS4DLL.casFIFOHasData (
    int ADevice ) [static]
```

Use this method to check if the spectrometer has data available which can be read

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

The return value is 0 if the spectrometer is still performing the acquisition, and 1 if FIFO data is ready to be read. Negative return values indicate an error which can be translated into an error message using [casGet↔ErrorMessage](#)

Calling [casFIFOHasData](#) is only necessary if a measurement was started with [casStart](#). It allows to check if the integration time of the CCD is over.

Note

Not all spectrometer types do support this - check [dpidTriggerCapabilities](#), namely the [tcoQueryHasData](#) bit.

8.1.2.16 `casGetCalibrationFactors()`

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetCalibrationFactors (
    int ADevice,
    int AWhat,
    int AIndex,
    int AExtra ) [static]
```

Returns various details about the configuration/calibration of the given device.

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	What calibration part to retrieve: one of the <code>gcf<XXX></code> constants, e.g. gcfDensityFunction .
<i>AINdex</i>	For calibration parts which consist of lists or spectra, this is the index or pixel for which to retrieve the calibration information. See the documentation for the <code>gcf<XXX></code> constant used as <i>AWhat</i> .
<i>AExtra</i>	Additional constant to indicate what calibration info to retrieve. Often another <code>gcf<XXX></code> constant, e.g. gcfDarkArrayDepth .

Returns

The return value is the calibration information identified by the parameter which were passed.

The following table lists the possible values for *AINdex* and *AExtra* for the different *AWhat* constants.

AWhat	AINdex	AExtra	Description
gcfDensityFunction	0..Pixel-1	Density Filter (0..7)	Spectral calibration of the density filters
gcfSensitivityFunction	0..Pixel-1	Unused	Spectral calibration of the CAS
gcfTransmissionFunction	0..Pixel-1	Unused	Transmission of additional optical hardware
gcfDensityFactor	0..7	Unused	Absolute calibration for the density filters
gcfTOPApertureFactor	0..6	Unused	Absolute calibration factor for TOP, aperture given by Index. Check against 0 to find out, which TOP apertures are calibrated.
gcfTOPDistanceFactor	gcfTDDistanceMin	Unused	The minimum calibrated TOP distance
	gcfTDDistanceMax	Unused	The maximum calibrated TOP distance
	gcfTDCount	Unused	Number of distance and factor pairs
	0..n	gcfTDExtraDistance	TOP distance in mm
		gcfTDExtraFactor	TOP calibration factor for this distance
gcfWLCalibration↔ Channel	gcfWLCalibPointCount	Unused	Number of calibration points of the wavelength calibration
	0..CalibPointCount-1	Unused	Wavelength calibration↔ : pixel of the CalibPoint specified by Index

gcfWLCalibrationAlias	0..CalibPointCount-1	Unused	Wavelength calibration↔ : wavelength of the CalibPoint specified by Index
gcfDarkArrayValues	unused	gcfDarkArrayDepth	Depth of the dark array, i.e. the number of dark currents in the array
	0..DarkArrayDepth-1	gcfDarkArrayIntTime	The integration time of the dark current measurement identified by Index.
	0..Pixel-1	0..DarkArrayDepth-1	The intensity of the dark current at the pixel given by Index and the dark current measurement identified by Extra.
gcfTOPParameter	0..6	gcfTOPApertureSize	Size of the TOP aperture in mm
	Unused	gcfTOPSpotSize↔Denominator	Denominator for calculation of the spot size of the TOP, rather use cas↔CalculateTOPParameter
	Unused	gcfTOPSpotSizeOffset	Offset for calculation of the spot size of the T↔OP, rather use cas↔CalculateTOPParameter
gcfLinearityFunction	0..n	gcfLinearityCounts	ADC counts for linearity correction
		gcfLinearityFactor	Correction Factor for A↔DC range
gcfRawData	0..Pixel - 1	Unused	Spectral raw data, normalized to 1 ms and 1 average

Some of these gcf<XXX> constants are also used by [casClearCalibration](#).

Warning

Since this method returns the actual calibration information, error handling has to be done using [casGet↔Error](#) calls. Typical error codes are [ErrorInvalidParameter](#) and [ErrorParamOutOfRange](#). When enumerating some calibration information lists, there might be no other way than increasing the AIndex parameter until [ErrorParamOutOfRange](#) occurs.

8.1.2.17 casGetCCT()

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetCCT (
    int ADevice ) [static]
```

Calculates the correlated color temperature CCT of a previously measured spectrum.

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

The return value is the CCT in K. Use [casGetError](#) to check for errors.

Before calling `casGetCCT`, a call to [casColorMetric](#) is necessary, to calculate all results from the current spectrum.

As with all colormetric calculations, the range of the spectrum which is used for calculating the CCT is defined by [mpidColormetricStart](#) and [mpidColormetricStop](#). The calculation is also affected by [mpidObserver](#).

Use [casGetError](#) to check for errors.

8.1.2.18 `casGetCentroid()`

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetCentroid (
    int ADevice ) [static]
```

Retrieves the previously calculated centroid wavelength.

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

The return value is the centroid wavelength in nm. Use [casGetError](#) to check for errors.

Before calling `casGetCentroid`, a call to [casColorMetric](#) is necessary, to calculate all results from the current spectrum.

8.1.2.19 `casGetColorCoordinates()`

```
static void InstrumentSystems.CAS4.CAS4DLL.casGetColorCoordinates (
    int ADevice,
    ref double x,
    ref double y,
    ref double z,
    ref double u,
    ref double v1976,
    ref double v1960 ) [static]
```

Retrieves CIE color coordinates of a previously measured spectrum.

Parameters

<i>ADevice</i>	The device / CASID
<i>x</i>	Receives the x color coordinate according to CIE 1931
<i>y</i>	Receives the y color coordinate according to CIE 1931
<i>z</i>	Receives the z color coordinate according to CIE 1931
<i>u</i>	Receives the u / u' color coordinate according to CIE 1960 / CIE 1976
<i>v1976</i>	Receives the v' color coordinate according to CIE 1976
<i>v1960</i>	Receives the v color coordinate according to CIE 1960

Before calling `casGetColorCoordinates`, a call to `casColorMetric` is necessary, to calculate all results from the current spectrum.

As with all colormetric calculations, the range of the spectrum which is used for calculating the color coordinates is defined by `mpidColormetricStart` and `mpidColormetricStop`. The calculation is also affected by `mpidObserver`.

Use `casGetError` to check for errors.

8.1.2.20 `casGetCRI()`

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetCRI (
    int ADevice,
    int Index ) [static]
```

Retrieves one of the previously calculated the color rendering indices.

Parameters

<i>ADevice</i>	The device / CASID
<i>Index</i>	Index ranging from 0 to 16, identifying which index should be returned. 0 returns the common index, which is the average of the indices 1..16.

Returns

The return value is the CRI. Use `casGetError` to check for errors.

The following preconditions must have been met, before calling `casGetCRI`:

- the device must have a valid spectrum
- set CRI calculation mode globally using `mpidCRIMode`
- colormetric calculation using `casColorMetric`
- calculation of the CCT using (`casGetCCT`)
- calculation of the CRI using (`casCalculateCRI`)

8.1.2.21 `casGetDarkCurrent()`

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetDarkCurrent (
    int ADevice,
    int AIndex ) [static]
```

Returns the intensity of the previously measured/calculated dark current for the given pixel.

Parameters

<i>ADevice</i>	The device / CASID
<i>AIndex</i>	The 0-based index of the pixel in the complete pixel array for which the dark current should be returned

Returns

The return value is the dark current for the given pixel in negative calibration units. Use [casGetError](#) to check for [ErrorInvalidParameter](#) and other errors.

AIndex can range from 0 to [dpidPixels](#) - 1, but the valid dark current ranges from [dpidDeadPixels](#) to [dpidDeadPixels](#) + [dpidVisiblePixels](#) - 1.

Use [casGetXArray](#) to retrieve the corresponding wavelength and [casGetData](#) to retrieve the intensity of a measured spectrum.

8.1.2.22 casGetData()

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetData (
    int ADevice,
    int AIndex ) [static]
```

Returns the intensity of the previously acquired spectrum for the given pixel.

Parameters

<i>ADevice</i>	The device / CASID
<i>AIndex</i>	The 0-based index of the pixel in the complete pixel array for which the intensity should be returned

Returns

The return value is the intensity for the given pixel in the calibration unit [dpidCalibrationUnit](#). Use [casGetError](#) to check for [ErrorInvalidParameter](#) and other errors.

AIndex can range from 0 to [dpidPixels](#) - 1, but the valid calibrated spectrum ranges from [dpidDeadPixels](#) to [dpidDeadPixels](#) + [dpidVisiblePixels](#) - 1.

Use [casGetXArray](#) to retrieve the corresponding wavelength. A previously measured dark current can be retrieved using [casGetDarkCurrent](#).

[Visual Basic]

```
DeadPixels = casGetDeviceParameter(CasID, dpidDeadPixels)
VisiblePixels = casGetDeviceParameter(CasID, dpidVisiblePixels)

'get Lambda and Intensity, skip dead pixels!
For i = 0 To VisiblePixels - 1
    Spectrum(0, i) = casGetXArray(CasID, DeadPixels + i)
    Spectrum(1, i) = casGetData(CasID, DeadPixels + i)
Next i
```

8.1.2.23 casGetDeviceParameter()

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetDeviceParameter (
    int ADevice,
    int AWhat ) [static]
```

Retrieves a float representing a device parameter

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	dpid constant identifying which device parameter to retrieve

Returns

Either an error code or a float whose meaning depends on the *AWhat* parameter. Refer to the documentation of the dpid constant used for *AWhat*.

8.1.2.24 `casGetDeviceParameterString()`

```
static int InstrumentSystems.CAS4.CAS4DLL.casGetDeviceParameterString (
    int ADevice,
    int AWhat,
    StringBuilder ADest,
    int AMaxLen ) [static]
```

Retrieves a string representing a device parameter

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	dpid constant identifying which device parameter to retrieve
<i>ADest</i>	StringBuilder/PAnsiChar buffer to retrieve the string
<i>AMaxLen</i>	number of characters ADest can hold - including a terminating zero

Returns

0 if successful or an error code

8.1.2.25 `casGetDeviceTypeName()`

```
static IntPtr InstrumentSystems.CAS4.CAS4DLL.casGetDeviceTypeName (
    int AInterfaceType,
    StringBuilder ADest,
    int AMaxLen ) [static]
```

Retrieves the name of the given interface type

Parameters

<i>AInterfaceType</i>	The interface type for which to retrieve the name
<i>ADest</i>	Destination string for the interface type name
<i>AMaxLen</i>	The maximum number of characters ADest can hold

Returns

The return value should not be used!

Use this method when iterating over all interface names. For `AInterfaceType` use constants from 0 to [casGetDeviceTypes](#) - 1.

Note

If an empty string is returned for an interface name, this interface type should not be used nor presented to the user!

For more details on interface types and options, refer to chapter [Interfaces Types and Options](#)

8.1.2.26 casGetDeviceTypeOption()

```
static int InstrumentSystems.CAS4.CAS4DLL.casGetDeviceTypeOption (
    int AInterfaceType,
    int AIndex ) [static]
```

Returns the value of the interface option for the given interface type and option index

Parameters

<i>AInterfaceType</i>	The interface type for which to retrieve the option
<i>AIndex</i>	0-based index of the option

Returns

Returns an interface option. Do not interpret as an error code!

Use this method when iterating over all interface options. `AIndex` can range from 0 to [casGetDeviceTypeOptions](#) - 1. The returned interface option can be used when calling [casCreateDeviceEx](#) or [casChangeDevice](#).

For more details on interface types and options, refer to chapter [Interfaces Types and Options](#)

8.1.2.27 casGetDeviceTypeOptionName()

```
static IntPtr InstrumentSystems.CAS4.CAS4DLL.casGetDeviceTypeOptionName (
    int AInterfaceType,
    int AInterfaceOptionIndex,
    StringBuilder ADest,
    int AMaxLen ) [static]
```

Retrieves the name of the given interface option

Parameters

<i>AInterfaceType</i>	The interface type for which to retrieve the interface option name
<i>AInterfaceOptionIndex</i>	The 0-based index of the interface option for which to retrieve the name
<i>ADest</i>	Destination string for the interface option name
<i>AMaxLen</i>	The maximum number of characters <code>ADest</code> can hold

Returns

The return value should not be used!

Use this method when iterating over all interface options. For `AInterfaceOptionIndex` use constants from 0 to [cas↔GetDeviceTypeOptions](#) - 1.

For more details on interface types and options, refer to chapter [Interfaces Types and Options](#)

8.1.2.28 casGetDeviceTypeOptions()

```
static int InstrumentSystems.CAS4.CAS4DLL.casGetDeviceTypeOptions (
    int AInterfaceType ) [static]
```

Retrieves the number of options a given interface type currently supports

Parameters

<i>AInterfaceType</i>	The interface type for which to retrieve the number of options
-----------------------	--

Returns

Returns a negative error code like [ErrorInvalidParameter](#) or a number indicating how many interface options the given interface type supports

Use this method when iterating over interface options using [casGetDeviceTypeOption](#)

For more details on interface types and options, refer to chapter [Interfaces Types and Options](#)

8.1.2.29 casGetDeviceTypes()

```
static int InstrumentSystems.CAS4.CAS4DLL.casGetDeviceTypes ( ) [static]
```

Retrieves the number of interface types the CAS DLL supports

Returns

Returns a positive number indicating how many interface types are supported

Use this method when iterating over all interface names using [casGetDeviceTypeName](#)

For more details on interface types and options, refer to chapter [Interfaces Types and Options](#)

8.1.2.30 casGetDigitalIn()

```
static int InstrumentSystems.CAS4.CAS4DLL.casGetDigitalIn (
    int ADevice,
    int APort ) [static]
```

Returns the state of the specified digital input port.

Parameters

<i>ADevice</i>	The device / CASID
<i>APort</i>	Identifies which DigitalIn port state to retrieve. The possible values differ with the spectrometer type and interface - refer to the hardware manual for more information.

Returns

The return value is 0 for low and 1 for high state. If APort is invalid for this device, [ErrorParamOutOfRange](#) will be returned. Other error codes are possible

8.1.2.31 `casGetDigitalOut()`

```
static int InstrumentSystems.CAS4.CAS4DLL.casGetDigitalOut (
    int ADevice,
    int APort ) [static]
```

Returns the state the specified digital output port was set to by a former call to [casSetDigitalOut](#).

Parameters

<i>ADevice</i>	The device / CASID
<i>APort</i>	Identifies which DigitalOut port state to retrieve. The possible values differ with the spectrometer type and interface - refer to the hardware manual for more information.

Returns

The return value is 0 for low and 1 for high state. If APort is invalid for this device, [ErrorParamOutOfRange](#) will be returned. Other error codes are possible

APort=1 is pin 8 and APort=2 for pin 7 of the DB9 connector of the PCI/ISA card or the trigger connector on the rear panel of the CAS140CT with USB interface.

8.1.2.32 `casGetDLLFileName()`

```
static IntPtr InstrumentSystems.CAS4.CAS4DLL.casGetDLLFileName (
    StringBuilder Dest,
    int AMaxLen ) [static]
```

Retrieves the complete path and file name of CAS library

Parameters

<i>Dest</i>	A buffer which can hold at least the number of characters specified in AMaxLen
<i>AMaxLen</i>	Number of characters including a trailing zero that Dest can hold

Returns

The return value is the pointer passed in Dest and can be ignored

There are A and W overloads for ANSI and Unicode versions of the method. The returned string is null-terminated.

8.1.2.33 casGetDLLVersionNumber()

```
static IntPtr InstrumentSystems.CAS4.CAS4DLL.casGetDLLVersionNumber (
    StringBuilder Dest,
    int AMaxLen ) [static]
```

Retrieves the version of CAS library

Parameters

<i>Dest</i>	A buffer which can hold at least the number of characters specified in AMaxLen
<i>AMaxLen</i>	Number of characters including a trailing zero that Dest can hold

Returns

The return value is the pointer passed in Dest and can be ignored

There are A and W overloads for ANSI and Unicode versions of the method. The returned string is null-terminated.

Note

some platforms might not support version info and will return "N/A".

8.1.2.34 casGetError()

```
static int InstrumentSystems.CAS4.CAS4DLL.casGetError (
    int ADevice ) [static]
```

Return error code for a given device/CASID

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Some methods provided by the CAS DLL do not return error codes. Call casGetError afterwards to check for errors which might have occurred during these method calls. A negative value indicates an error. [ErrorNoError](#) indicates that the previous action was successful. See chapter [Error Handling](#) for more details.

Note

Every subsequent call into the CAS DLL will clear the previous error for the given device! [casGetErrorMessage](#) can be used to translate the returned error into an error message.

8.1.2.35 casGetErrorMessage()

```
static IntPtr InstrumentSystems.CAS4.CAS4DLL.casGetErrorMessage (
    int AError,
    StringBuilder ADest,
    int AMaxLen ) [static]
```

Translates a given error code into a readable error message

Parameters

<i>AError</i>	The error code which should be translated
<i>ADest</i>	Destination for the error message
<i>AMaxLen</i>	The maximum number of characters ADest can hold

Returns

The return value should not be used

Use this method to translate an error constant (AError) into a user-readable error message. The error constant can either be the return value of a method of CAS DLL or can be explicitly retrieved using [casGetError](#). The message is copied into the buffer ADest is pointing to. This buffer must be able to hold at least the number of characters specified by AMaxLen plus a trailing zero.

8.1.2.36 casGetExtendedColorValues()

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetExtendedColorValues (
    int ADevice,
    int AWhat ) [static]
```

Retrieves results and conditions for previously performed colormetric calculations.

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	Determines which colormetric value to return. One of the ecv constants starting with ecvVisualEffect

Returns

The return value depends on the ecv constant passed as What parameter. Use [casGetError](#) to check for errors.

Before calling [casGetExtendedColorValues](#) the values have to be calculated by calling [casColorMetric](#), [casGetCCT](#) and [casCalculateCRI](#). Refer to these methods to find out which other properties influence their results, like colormetric range or [mpidCRIMode](#).

The ecv constants starting with [ecvCRIFirst](#) allow access to intermediate- and end-results of the CRI calculation. The following table lists their meaning when [mpidCRIMode](#) is [criDIN6169](#) and [mpidColormetricType](#) is [cmtDefaultColormetric](#). For [cmtSWPColormetric](#) these results are not available and will return 0.

AWhat range	Description
ecvCRIFirst..ecvCRILast	Returns the CRI value, identical to casGetCRI , i.e. AWhat = ecvCRIFirst returns the common CRI, ecvCRIFirst + 1 .. ecvCRILast return CRIs 1 to 16.
ecvCRITriKXFirst..ecvCRITriKXLast	Returns the Tristimulus X of the spectrum without a test color (ecvCRITriKXFirst) or with test color 1..16 (ecvCRITriKXFirst + 1 .. ecvCRITriKXLast)
ecvCRITriKYFirst..ecvCRITriKYLlast	Returns the Tristimulus Y of the spectrum without a test color (ecvCRITriKYFirst) or with test color 1..16 (ecvCRITriKYFirst + 1 .. ecvCRITriKYLlast)
ecvCRITriKZFirst..ecvCRITriKZLast	Returns the Tristimulus Z of the spectrum without a test color (ecvCRITriKZFirst) or with test color 1..16 (ecvCRITriKZFirst + 1 .. ecvCRITriKZLast)
ecvCRITriRXordUFirst..ecvCRITriRXordULast	Returns the Tristimulus X of the calculated reference spectrum without a test color (ecvCRITriRXordUFirst) or with test color 1..16 (ecvCRITriRXordUFirst + 1 .. ecvCRITriRXordULast)
ecvCRITriRYordVFirst..ecvCRITriRYordVLast	Returns the Tristimulus Y of the calculated reference spectrum without a test color (ecvCRITriRYordVFirst) or with test color 1..16 (ecvCRITriRYordVFirst + 1 .. ecvCRITriRYordVLast)
ecvCRITriRZordWFirst..ecvCRITriRZordWLast	Returns the Tristimulus Z of the calculated reference spectrum without a test color (ecvCRITriRZordWFirst) or with test color 1..16 (ecvCRITriRZordWFirst + 1 .. ecvCRITriRZordWLast)

The following table lists their meaning when [mpidCRIMode](#) is [criCIE13_3_95](#), regardless of [mpidColormetricType](#).

AWhat range	Description
ecvCRIFirst..ecvCRILast	Returns the CRI value, identical to casGetCRI , i.e. AWhat = ecvCRIFirst returns the common CRI, ecvCRIFirst + 1 .. ecvCRILast return CRIs 1 to 16.
ecvCRITriKXFirst..ecvCRITriKXLast	Returns the Tristimulus X of the spectrum without a test color (ecvCRITriKXFirst) or with test color 1..16 (ecvCRITriKXFirst + 1 .. ecvCRITriKXLast)
ecvCRITriKYFirst..ecvCRITriKYLlast	Returns the Tristimulus Y of the spectrum without a test color (ecvCRITriKYFirst) or with test color 1..16 (ecvCRITriKYFirst + 1 .. ecvCRITriKYLlast)
ecvCRITriKZFirst..ecvCRITriKZLast	Returns the Tristimulus Z of the spectrum without a test color (ecvCRITriKZFirst) or with test color 1..16 (ecvCRITriKZFirst + 1 .. ecvCRITriKZLast)
ecvCRITriRXordUFirst + 1..ecvCRITriRXordULast	Returns dU for CRI 1..16 (ecvCRITriRXordUFirst + 1 .. ecvCRITriRXordULast)
ecvCRITriRYordVFirst + 1..ecvCRITriRYordVLast	Returns dV for CRI 1..16 (ecvCRITriRYordVFirst + 1 .. ecvCRITriRYordVLast)
ecvCRITriRZordWFirst + 1..ecvCRITriRZordWLast	Returns dW for CRI 1..16 (ecvCRITriRZordWFirst + 1 .. ecvCRITriRZordWLast)

8.1.2.37 [casGetExternalADCValue\(\)](#)

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetExternalADCValue (
    int ADevice,
    int AIndex ) [static]
```

Obsolete. Use [mpidCurrentCCDTemperature](#). Used to retrieve the CCD temperature.

Parameters

<i>ADevice</i>	The device / CASID
<i>AIndex</i>	The ADC channel to retrieve

Returns

The return value is the CCD temperature

This method should no longer be used. It might not work with some CAS models! Refer to chapter [CCD Temperature Monitoring](#) for a complete overview of the topic.

8.1.2.38 casGetFIFOData()

```
static int InstrumentSystems.CAS4.CAS4DLL.casGetFIFOData (
    int ADevice ) [static]
```

This method reads the acquired spectrum from the FIFO and stores it internally

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

Negative return values indicate an error which can be translated into an error message using [casGetError↔Message](#). All other return values indicate that the spectrum was read successfully

The spectrum is sensitivity corrected and subsequent averages are performed if [mpidAverages](#) is > 1.

Calling [casGetFIFOData](#) should only be done after a measurement was started with [casStart](#).

To access the acquired spectrum use [casGetData](#)

8.1.2.39 casGetFilterName()

```
static IntPtr InstrumentSystems.CAS4.CAS4DLL.casGetFilterName (
    int ADevice,
    int AFilter,
    StringBuilder Dest,
    int AMaxLen ) [static]
```

Translates a filter index into a user-readable density filter name.

Parameters

<i>ADevice</i>	The device / CASID
<i>AFilter</i>	The 0-based filter index ranging from 0 to 7
<i>Dest</i>	Destination buffer for the filter name
<i>AMaxLen</i>	The maximum number of characters Dest can hold

Returns

The return value is identical to Dest or nil if there was an error

The density filter names are defined in the configuration file of the device, so the device must have been initialized (see [casInitialize](#)). A default name "Filter x" will be returned when the configuration doesn't define a name for a given filter.

8.1.2.40 casGetMeasurementParameter()

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetMeasurementParameter (
    int ADevice,
    int AWhat ) [static]
```

Returns a numeric measurement parameter for a given spectrometer

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	Integer constant, defining which parameter to return. One of the mpid<XXX> constants starting at mpidIntegrationTime

Returns

The return value is the measurement parameter itself and value depends on AWhat. Call [casGetError](#) immediately afterwards to do proper See [Error Handling](#)

After performing a measurement, it is often recommended to query the various measurement parameter to build a list of measurement conditions for documentation.

8.1.2.41 casGetOptions()

```
static int InstrumentSystems.CAS4.CAS4DLL.casGetOptions (
    int ADevice ) [static]
```

Returns the features and options the device currently supports or performs.

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

An integer that has the bits starting with [coShutter](#) set - each bit standing for an option or feature.

For historic reasons the options is a mixture of capabilities - which are defined by the device and cannot be changed - and actual options that can be enabled or disabled.

Use the method [casSetOptionsOnOff](#) to turn specific options on or off.

The following example checks if the device supports measuring the CCD temperature:

[Delphi]

```
if (casGetOptions(CASID) and coGetTemperature) <> 0 then
  ShowMessage('Device supports measuring CCD temperature') else
  ShowMessage('Device does NOT support measuring CCD temperature');
```

8.1.2.42 casGetPeak()

```
static void InstrumentSystems.CAS4.CAS4DLL.casGetPeak (
    int ADevice,
    out double x,
    out double y ) [static]
```

Retrieves the previously calculated peak wavelength and intensity.

Parameters

<i>ADevice</i>	The device / CASID
<i>x</i>	Receives the wavelength in nm of the interpolated peak
<i>y</i>	Receives the intensity in calibration units of the interpolated peak

Before calling `casGetPeak`, a call to `casColorMetric` is necessary, to calculate all results from the current spectrum.

Use `casGetError` to check for errors.

8.1.2.43 casGetPhotInt()

```
static void InstrumentSystems.CAS4.CAS4DLL.casGetPhotInt (
    int ADevice,
    out double APhotInt,
    StringBuilder AUnit,
    int AUnitMaxLen ) [static]
```

Retrieves the previously calculated photometric integral and it's data unit.

Parameters

<i>ADevice</i>	The device / CASID
<i>APhotInt</i>	Double that will contain the photometric integral after a successful method call.
<i>AUnit</i>	Destination buffer for string representing the photometric data unit
<i>AUnitMaxLen</i>	The maximum number of characters AUnit can hold

Before calling `casGetPhotInt`, a call to `casColorMetric` is necessary, to calculate all results from the current spectrum.

The value and data unit will always be returned using the basic unit, e.g. `APhotInt = 1.23E-5` and `AUnit = "lx"` and not in "mlx" or "klx".

The integral is calculated from 380nm to 780nm for `cmtDefaultColormetric` or from 360nm to 830nm for `cmtSW↔PColormetric` unless further constrained using `mpidColormetricStart` and `mpidColormetricStop`. For integration, the spectrum is multiplied with $V(\lambda)$ and the photometric equivalent Km. It is independent of `mpidObserver`.

Use `casGetError` to check for errors.

`casGetRadInt` retrieves the radiometric integral.

8.1.2.44 casGetRadInt()

```
static void InstrumentSystems.CAS4.CAS4DLL.casGetRadInt (
    int ADevice,
    out double ARadInt,
    StringBuilder AUnit,
    int AUnitMaxLen ) [static]
```

Retrieves the previously calculated radiometric integral and it's data unit.

Parameters

<i>ADevice</i>	The device / CASID
<i>ARadInt</i>	Double that will contain the radiometric integral after a successful method call.
<i>AUnit</i>	Destination buffer for string representing the radiometric data unit
<i>AUnitMaxLen</i>	The maximum number of characters AUnit can hold

Before calling `casGetRadInt`, a call to [casColorMetric](#) is necessary, to calculate all results from the current spectrum.

The value and data unit will always be returned using the basic unit, e.g. `ARadInt = 1.23E-5` and `AUnit = "W"` and not in "mW" or "kW".

The integral is calculated from 380nm to 780nm unless further constrained using [mpidColormetricStart](#) and [mpidColormetricStop](#).

Use [casGetError](#) to check for errors.

[casGetPhotInt](#) retrieves the photometric integral.

8.1.2.45 casGetSerialNumberEx()

```
static int InstrumentSystems.CAS4.CAS4DLL.casGetSerialNumberEx (
    int ADevice,
    int AWhat,
    StringBuilder ADest,
    int AMaxLen ) [static]
```

Retrieves serial numbers of the specified device and/or additional information.

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	constant starting with "casSerial" identifying which serial number to retrieve
<i>ADest</i>	StringBuilder/PAnsiChar buffer to retrieve the string
<i>AMaxLen</i>	number of characters ADest can hold - including a terminating zero

Returns

0 if successful or an error code

As opposed to [dpidSerialNo](#), this method can retrieve the serial number of the device and additional equipment separately. See chapter [Using Serial Numbers](#) for an overview. Possible values for `AWhat` start with [casSerialComplete](#)

8.1.2.46 casGetShutter()

```
static int InstrumentSystems.CAS4.CAS4DLL.casGetShutter (
    int ADevice ) [static]
```

Returns the current shutter position of the given device

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

The return value is one of the casShutter constants starting with [casShutterInvalid](#) or one of the error codes, e.g. [errCasDeviceNotFound](#)

Note that the actual shutter position of the device is only checked if the device supports the [coGetShutter](#). Otherwise the most recently set shutter state is returned, which may not reflect the actual shutter position.

8.1.2.47 casGetTriStimulus()

```
static void InstrumentSystems.CAS4.CAS4DLL.casGetTriStimulus (
    int ADevice,
    ref double X,
    ref double Y,
    ref double Z ) [static]
```

Retrieves the previously calculated X, Y and Z tristimulus values.

Parameters

<i>ADevice</i>	The device / CASID
<i>X</i>	Variable to receive the tristimulus X
<i>Y</i>	Variable to receive the tristimulus Y
<i>Z</i>	Variable to receive the tristimulus Z

Before calling `casGetTriStimulus`, [casColorMetric](#) must have been called for an already measured spectrum. Calculation is influenced by [mpidObserver](#). The tristimulus is calculated between 380 and 780 nm, but the range can be further restricted by [mpidColormetricStart](#) and [mpidColormetricStop](#).

8.1.2.48 casGetWidth()

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetWidth (
    int ADevice ) [static]
```

Retrieves the peak width in nm.

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

The return value is the peak width in nm. Use [casGetError](#) to check for errors.

`casGetWidth` calculates the width of the peak and returns just the width in nm. It is recommended to use [casGetWidthEx](#) instead, as it also returns other aspects of the peak width. The level of the peak width is 50% by default so the FWHM is calculated. It can be adjusted with [mpidColormetricWidthLevel](#) before calling `casGetWidth` or `casGetWidthEx`.

As with all colormetric calculations, the range of the spectrum which is used for determining the width is defined by [mpidColormetricStart](#) and [mpidColormetricStop](#).

8.1.2.49 casGetWidthEx()

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetWidthEx (
    int ADevice,
    int AWhat ) [static]
```

Calculates and retrieves various aspects about the peak width (e.g. full width half maximum aka FWHM aka 50% bandwidth).

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	One of the <code>cLambda</code> constants specifying what to retrieve; see table below

Returns

The return value is the detail about the width that was specified by *AWhat*. Use [casGetError](#) to check for errors.

By default the FWHM is calculated, i.e. the peak level used for the calculation is 50%. You can modify this level independently for each parameter set using [mpidColormetricWidthLevel](#) before calling `casGetWidthEx`.

The following table lists possible values for the *AWhat* parameter and a description which aspect of the peak width is returned. The meaning of the algorithm and calculates columns are explained below.

AWhat	Calculates?*	Algorithm**	Description
cLambdaWidth	Yes	Inner width	The peak width in nm
cLambdaLow	No	Inner width	The start wavelength of the peak width
cLambdaMiddle	No	Inner width	The center wavelength of the peak width, in the middle between start and stop wavelength
cLambdaHigh	No	Inner width	The stop wavelength of the peak width
cLambdaOuterWidth	Yes	Outer width	Returns the peak width
cLambdaOuterLow	Yes	Outer width	The start wavelength of the peak width
cLambdaOuterMiddle	Yes	Outer width	The center wavelength of the peak width, in the middle between start and stop wavelength
cLambdaOuterHigh	Yes	Outer width	The stop wavelength of the peak width

As with all colormetric calculations, the range of the spectrum which is used for determining the peak width is defined by [mpidColormetricStart](#) and [mpidColormetricStop](#).

* Some AWhat constants only retrieve a previously calculated value and some perform the actual calculation. For example, before you can retrieve cLambdaLow, you must retrieve cLambdaWidth, as only the latter performs the actual calculation using the Inner width algorithm

** Normally the peak width is determined by "walking" the spectrum up and down from the intensity maximum to check when the intensity falls below the level defined by [mpidColormetricWidthLevel](#). That is the so called "Inner width" algorithm which should normally be used. The "Outer width" algorithm is mainly supported because very old CAS DLL implementations (version 3 and before) used this. This algorithm actually starts walking from the ends of the colormetric range towards the center and stops when the width level of the intensity maximum is reached.

8.1.2.50 casGetXArray()

```
static double InstrumentSystems.CAS4.CAS4DLL.casGetXArray (
    int ADevice,
    int AIndex ) [static]
```

Returns the wavelength of the spectrum for the given pixel.

Parameters

<i>ADevice</i>	The device / CASID
<i>AINdex</i>	The 0-based index of the pixel in the complete pixel array for which the wavelength should be returned

Returns

The return value is the wavelength for the given pixel in nm. Use [casGetError](#) to check for [ErrorInvalidParameter](#) and other errors.

AINdex can range from 0 to [dpidPixels](#) - 1, but the valid calibrated spectrum ranges from [dpidDeadPixels](#) to [dpidDeadPixels](#) + [dpidVisiblePixels](#) - 1.

Use [casGetData](#) to retrieve the corresponding intensity at this wavelength.

Note

Contrary to [casGetData](#), [casGetXArray](#) does not require that a spectrum has been measured. A successful initialization with [casInitialize](#) is sufficient.

[Visual Basic]

```
DeadPixels = casGetDeviceParameter(CasID, dpidDeadPixels)
VisiblePixels = casGetDeviceParameter(CasID, dpidVisiblePixels)

'get Lambda and Intensity, skip dead pixels!
For i = 0 To VisiblePixels - 1
    Spectrum(0, i) = casGetXArray(CasID, DeadPixels + i)
    Spectrum(1, i) = casGetData(CasID, DeadPixels + i)
Next i
```

8.1.2.51 casInitialize()

```
static int InstrumentSystems.CAS4.CAS4DLL.casInitialize (
    int ADevice,
    int Perform ) [static]
```

Initializes the hardware of the device after loading the configuration and calibration files

Parameters

<i>ADevice</i>	The device / CASID
<i>Perform</i>	One of the following constants: InitOnce , InitForced and InitNoHardware

Returns

Returns 0 if the device was successfully initialized or a negative error constant like [errCasDeviceNotFound](#), [errCasNoConfig](#), [ErrorNoCalibration](#) or [ErrorAdrControl](#)

8.1.2.52 `casLoadTestData()`

```
static void InstrumentSystems.CAS4.CAS4DLL.casLoadTestData (
    int ADevice,
    string AFileName ) [static]
```

Loads a spectrum from an .ISD file into the CASID

Parameters

<i>ADevice</i>	The device / CASID
<i>AFileName</i>	Full null-terminated path of the .ISD file which should be loaded

Use `casLoadTestData` to load a spectrum from the ISD-file specified by `AFileName`. The spectrum is stored internally for the given device almost as if it would have just been measured, so a subsequent call to `casColorMetric` etc. will perform calculations for this spectrum.

Warning

`casLoadTestData` will only work for a CASID of a real spectrometer, which is connected to the PC and has been initialized. Otherwise it will return [ErrorFeatureNotSupported!](#)

Note

Loaded spectra will not work with something like [paRecalcSpectrum](#) since this recalculates the spectrum from the previously acquired raw data - which is not affected by `casLoadTestData`.

The loaded spectrum will be resampled to the wavelength calibration of the spectrometer, so it is only safe to use `casLoadTestData` with a CASID that has been initialized with the same calibration files that were used to measure the saved spectrum.

Call [casGetError](#) after `casLoadTestData` for error handling.

There are A and W overloads for ANSI and Unicode versions of the method.

8.1.2.53 `casMeasure()`

```
static int InstrumentSystems.CAS4.CAS4DLL.casMeasure (
    int ADevice ) [static]
```

Performs a measurement for the given device using the measurement parameter which have been previously set

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

Negative return values indicate an error which can be translated into an error message using [casGetError↵Message](#). All other return values indicate that the measurement was successful

The spectrum is stored internally and can be accessed using [casGetData](#).

Before starting a measurement, you should check whether the [density filter](#) needs to be updated and whether a [dark current](#) measurement is necessary. Calling [casPerformActionEx](#) with `paPrepareMeasurement` is also recommended.

See chapter [performing a measurement](#) for all prerequisites and more details. If the `AutoRange` option is enabled (see [casGetOptions](#)):

- the integration time is automatically determined (within the bounds of the [mpidAutoRangeMinLevel](#), [mpid↵AutoRangeMaxLevel](#) and [mpidAutoRangeMaxIntTime](#)). Refer to the chapter [AutoRange measurements](#) for more details.
- using a dark array may be helpful (see chapter [Dark current](#)), because `casMeasure` may automatically perform dark current measurements and even change the density filter if necessary and enabled ([coAutorangeFilter](#))

8.1.2.54 casMeasureDarkCurrent()

```
static int InstrumentSystems.CAS4.CAS4DLL.casMeasureDarkCurrent (
    int ADevice ) [static]
```

Use `casMeasureDarkCurrent` to perform a dark current measurement for the given device.

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

Negative return values indicate an error which can be translated into an error message using [casGetError↵Message](#). All other return values indicate that the dark current was measured successfully

Note

This method does not close the shutter! Call [casSetShutter](#) to close the shutter before calling `casMeasure↵DarkCurrent` and again to open it afterwards

Warning

Before measuring dark current, check whether the [density filter](#) needs to be moved!

Refer to chapter [Dark Current](#) for more details and sample code.

To access the measured dark current use [casGetDarkCurrent](#)

8.1.2.55 casMultiTrackCopyData()

```
static int InstrumentSystems.CAS4.CAS4DLL.casMultiTrackCopyData (
    int ADevice,
    int ATrack ) [static]
```

Releases a previously allocated MultiTrack buffer

Parameters

<i>ADevice</i>	The device / CASID
<i>ATrack</i>	The 0-based track index which should be copied/activated

Returns

The return value is 0 if the method was successful or a negative error code

casMultiTrackCopyData retrieves the spectrum specified by ATrack from the MultiTrack measurement buffer (A↔Track is 0-based, so it may range from 0 to [dpidMultiTrackCount](#) - 1). The MultiTrack measurement may have been performed ([paMultiTrackStart](#)) or loaded ([casMultiTrackLoadData](#)). The spectrum is stored internally as if it would have just been measured, so a subsequent call to [casColorMetric](#) calculates the colormetric results for it etc. For an overview, refer to the chapter [MultiTrack measurements](#).

8.1.2.56 casMultiTrackCount()

```
static int InstrumentSystems.CAS4.CAS4DLL.casMultiTrackCount (
    int ADevice ) [static]
```

Deprecated method.

Returns

Returns the number of MultiTracks already allocated. A negative value indicates an error

Use [dpidMultiTrackCount](#) instead

8.1.2.57 casMultiTrackDone()

```
static int InstrumentSystems.CAS4.CAS4DLL.casMultiTrackDone (
    int ADevice ) [static]
```

Releases a previously allocated MultiTrack buffer

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

The return value is 0 if the method was successful or a negative error code

`casMultiTrackDone` releases the MultiTrack buffer after it has been allocated with `casMultiTrackInit` or [casMultiTrackLoadData](#). Obviously, after that, the raw spectra in the MultiTrack buffer are no longer available. For an overview, refer to the chapter [MultiTrack measurements](#).

8.1.2.58 casMultiTrackInit()

```
static int InstrumentSystems.CAS4.CAS4DLL.casMultiTrackInit (
    int ADevice,
    int ATracks ) [static]
```

Initializes the MultiTrack buffer

Parameters

<i>ADevice</i>	The device / CASID
<i>ATracks</i>	The number of MultiTracks that should later be measured

Returns

The return value is 0 if the method was successful or a negative error code

`casMultiTrackInit` prepares a [MultiTrack measurements](#) for the given device. `ATracks` defines the number of tracks (i.e. raw spectra) which will be acquired during the MultiTrack measurement and should not exceed [dpidMultiTrackMaxCount](#). Since the memory used for saving MultiTrack measurements needs to be allocated beforehand, it is not possible to change the number of tracks on the fly. You may use [dpidMultiTrackCount](#) to check how many tracks have been allocated.

Note

It is mandatory to call [casMultiTrackDone](#) to discard the allocated MultiTrack memory when the MultiTrack data is no longer needed. `casMultiTrackInit` and [casMultiTrackLoadData](#) implicitly call [casMultiTrackDone](#) before initializing/loading new MultiTrack data.

8.1.2.59 casMultiTrackLoadData()

```
static int InstrumentSystems.CAS4.CAS4DLL.casMultiTrackLoadData (
    int ADevice,
    string AFileName ) [static]
```

Loads a MultiTrack buffer from a .SWM file.

Parameters

<i>ADevice</i>	The device / CASID
<i>AFileName</i>	Full null-terminated path of the .ISD file which should be loaded

Returns

The return value is the number of tracks loaded if the method was successful or a negative error code

Loads MultiTrack data for the specified device from the file specified by *AFileName* and discards any MultiTrack data which may have previously been loaded or acquired ([casMultiTrackDone](#) is called implicitly). The file must have been saved with [casMultiTrackSaveData](#). Use [dpidMultiTrackCount](#) to check how many tracks have been loaded. If you want to retrieve colormetric results or the spectra, call [casMultiTrackCopyData](#), but make sure the same configuration and calibration files are used for the device!

Note

It is mandatory to call [casMultiTrackDone](#) to discard the allocated MultiTrack memory, when the MultiTrack data is no longer needed. For an overview, refer to the chapter [MultiTrack measurements](#). There are A and W overloads for ANSI and Unicode versions of the method.

8.1.2.60 casMultiTrackSaveData()

```
static int InstrumentSystems.CAS4.CAS4DLL.casMultiTrackSaveData (
    int ADevice,
    string AFileName ) [static]
```

Saves a MultiTrack buffer to a .SWM file.

Parameters

<i>ADevice</i>	The device / CASID
<i>AFileName</i>	Full null-terminated path of the .SWM file which should be created

Returns

The return value is 0 if the method was successful or a negative error code

Saves the MultiTrack buffer to *AFileName*. Depending on the number of tracks, the file can get rather big and saving it may take a while. The resulting .SWM file can be loaded by SpecWin Pro or by calling [casMultiTrackLoadData](#). An existing file at *AFileName* will be overwritten. For an overview, refer to the chapter [MultiTrack measurements](#). There are A and W overloads for ANSI and Unicode versions of the method.

8.1.2.61 casNmToPixel()

```
static int InstrumentSystems.CAS4.CAS4DLL.casNmToPixel (
    int ADevice,
    double nm ) [static]
```

Converts a wavelength into the corresponding CCD pixel index.

Parameters

<i>ADevice</i>	The device / CASID
<i>nm</i>	The wavelength which should be converted into a pixel

Returns

The return value is the 0-based CCD pixel index or a negative [error code](#).

The conversion is done using the wavelength calibration of the specified device. Therefore the device must have been initialized using [casInitialize](#) or the wavelength calibration updated with [paUpdateCompleteCalibration](#). The return value is the first 0-based pixel whose wavelength is greater or equal than the nm parameter. Only visible pixel are taken into account. A negative return value indicates an error. [casPixelToNm](#) does the conversion the other way around. The returned pixel index can be used with methods like [casGetXArray](#), [casGetData](#) or [casGetDarkCurrent](#).

8.1.2.62 casPerformAction()

```
static int InstrumentSystems.CAS4.CAS4DLL.casPerformAction (
    int ADevice,
    int AID ) [static]
```

Deprecated. Use [casPerformActionEx](#) instead!

8.1.2.63 casPerformActionEx()

```
static int InstrumentSystems.CAS4.CAS4DLL.casPerformActionEx (
    int ADevice,
    int AActionID,
    int AParam1,
    int AParam2,
    IntPtr AParam3 ) [static]
```

Generic method which performs one of various actions with the specified device. This method replaces the now deprecated method `casPerformAction`.

Parameters

<i>ADevice</i>	The device / CASID. For some AActionID constants, a valid CASID is not required, so any value will be accepted.
<i>AActionID</i>	Integer defining the action to perform. One of the constants starting at paPrepareMeasurement
<i>AParam1</i>	First integer parameter for the action to perform. Meaning depends on the action to be performed
<i>AParam2</i>	Second integer parameter for the action to perform. Meaning depends on the action to be performed
<i>AParam3</i>	Third pointer parameter for the action to perform. Meaning depends on the action to be performed

Returns

0 for success or a negative error code.

8.1.2.64 casPixelToNm()

```
static double InstrumentSystems.CAS4.CAS4DLL.casPixelToNm (
    int ADevice,
    int APixel ) [static]
```

Converts a CCD pixel index into the corresponding wavelength.

Parameters

<i>ADevice</i>	The device / CASID
<i>APixel</i>	The 0-based CCD pixel for which to retrieve the corresponding wavelength

Returns

The return value is the wavelength. Use [casGetError](#) for error checking.

The conversion is done using the wavelength calibration of the specified device. Therefore the device must have been initialized using [casInitialize](#) or the wavelength calibration updated with [paUpdateCompleteCalibration](#). [casNmToPixel](#) does the conversion the other way around.

8.1.2.65 casSaveCalibration()

```
static void InstrumentSystems.CAS4.CAS4DLL.casSaveCalibration (
    int ADevice,
    string AFileName ) [static]
```

Saves the calibration of the given device to a file.

Parameters

<i>ADevice</i>	The device / CASID
<i>AFileName</i>	The full path where the calibration file should be stored. Typical file extension is .ISC

This calibration file does not contain all details of the calibration, e.g. the wavelength calibration - which is saved in the configuration file. Use [casSetCalibrationFactors](#) with `gfcWLCalibrationSave` to save the wavelength calibration in the current configuration file.

8.1.2.66 casSaveSpectrum()

```
static int InstrumentSystems.CAS4.CAS4DLL.casSaveSpectrum (
    int ADevice,
    string AFileName ) [static]
```

Saves a previously measured spectrum to an .ISD file.

Parameters

<i>ADevice</i>	The device / CASID
<i>AFileName</i>	Full null-terminated path of the .ISD file which should be created.

Returns

The return value is 0 if the method was successful or a negative error code

The ISD file format is an ANSI/UTF-8 text file which can be loaded by SpecWin Pro and other programs. An existing file at AFileName will be overwritten. There are A and W overloads for ANSI and Unicode versions of the method.

8.1.2.67 casSetCalibrationFactors()

```
static void InstrumentSystems.CAS4.CAS4DLL.casSetCalibrationFactors (
    int ADevice,
    int AWhat,
    int AIndex,
    int AExtra,
    double AValue ) [static]
```

Changes various details about the configuration/calibration of the given device effectively overriding information that normally comes from configuration/calibration files.

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	What calibration part to modify: one of the gcf<XXX> constants, e.g. gcfDensityFunction .
<i>AINdex</i>	For calibration parts which consist of lists or spectra, this is the index or pixel for which to modify the calibration information. See the documentation for the gcf<XXX> constant used as AWhat.
<i>AExtra</i>	Additional constant to indicate what calibration info to change. Often another gcf<XXX> constant, e.g. gcfDarkArrayDepth .
<i>AValue</i>	New value for the calibration information that is identified by the other parameter.

The following table lists the possible values for AIndex and AExtra for the different AWhat constants.

AWhat	AINdex	AExtra	Description
gcfDensityFunction	0..Pixel-1	Density Filter (0..7)	Spectral calibration of the density filters
gcfSensitivityFunction	0..Pixel-1	Unused	Spectral calibration of the CAS
gcfTransmissionFunction	0..Pixel-1	Unused	Transmission of additional optical hardware
gcfDensityFactor	0..7	Unused	Absolute calibration for the density filters
gcfTOPApertureFactor	0..6	Unused	Absolute calibration factor for TOP, aperture given by Index. Check against 0 to find out, which TOP apertures are calibrated.
gcfTOPDistanceFactor	gcfTDCount	Unused	Number of distance and factor pairs
	0..n	gcfTDExtraDistance	TOP distance in mm
		gcfTDExtraFactor	TOP calibration factor for this distance
gcfWLCalibration↔Channel	Channel/Pixel to delete	gcfWLExtraCalibration↔Delete	Deletes the calibration point for the channel/pixel value passed in Index. No error if no calibration point at this pixel.
	Unused	gcfWLExtraCalibration↔DeleteAll	Deletes all calibration points. Identical to calling casClearCalibration with gcfWLCalibration↔Channel

	Use gcfWLCalibrationAlias to define calibration points in one go.		
gcfWLCalibrationAlias	0..Pixels-1	Unused	Add/modify calibration point for pixel/channel given by Index and set it to the wavelength passed in AValue.
gcfWLCalibrationSave	Unused	Unused	Immediately saves the modified wavelength calibration in the current configuration file specified by dpidConfigFileName .
gcfDarkArrayValues	unused	gcfDarkArrayDepth	Depth of the dark array, i.e. the number of dark currents in the array
	0..DarkArrayDepth-1	gcfDarkArrayIntTime	Sets the integration time of the DarkArray at Index to the integration time passed in value. Note↔ : the integration times in the dark array must be in ascending order and start with dpidIntTimeMin !
gcfLinearityFunction	0..n	ADCCounts	Adds/modifies the linearity correction array to the correction factor passed in AValue for the ADC counts passed in AIndex. The array must be sorted in ascending order by A↔ DCCounts.

Warning

After changing the calibration [paUpdateSpectralCalibration](#) should be called before the next measurement is performed, so the new calibration factors are used! Some actions (like changing the active parameter set or the density filter) do this implicitly so it might not always be necessary. The wavelength calibration is updated automatically every time it is modified.

Use [casSaveCalibration](#) if you want to save the modified calibration as a .ISC file. The wavelength calibration can be saved by calling [casSetCalibrationFactors](#) with [gcfWLCalibrationSave](#) as mentioned above. All the other information that is stored in the configuration file cannot be saved by the CAS-DLL.

8.1.2.68 [casSetDeviceParameter\(\)](#)

```
static int InstrumentSystems.CAS4.CAS4DLL.casSetDeviceParameter (
    int ADevice,
    int AWhat,
    double AValue ) [static]
```

Sets a numerical device parameter

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	dpid constant identifying which device parameter to change
<i>AValue</i>	the numerical value the device parameter should be set to

Returns

0 if successful or an error code

8.1.2.69 casSetDeviceParameterString()

```
static int InstrumentSystems.CAS4.CAS4DLL.casSetDeviceParameterString (
    int ADevice,
    int AWhat,
    string AValue ) [static]
```

Sets a string device parameter

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	dpid constant identifying which device parameter to change
<i>AValue</i>	the null terminated string value the device parameter should be set to

Returns

0 if successful or an error code

8.1.2.70 casSetDigitalOut()

```
static void InstrumentSystems.CAS4.CAS4DLL.casSetDigitalOut (
    int ADevice,
    int APort,
    int OnOff ) [static]
```

Sets the state the specified digital output port.

Parameters

<i>ADevice</i>	The device / CASID
<i>APort</i>	Identifies which DigitalOut port state to change. The possible values differ with the spectrometer type and interface - refer to the hardware manual for more information.
<i>OnOff</i>	The desired state: 0 for low, all other values for high state

See [casGetDigitalOut](#) for more information.

Warning

Proper error handling using [casGetError](#) is essential! If APort is not supported by the device, [ErrorParam↵ OutOfRange](#) will be returned!

8.1.2.71 casSetMeasurementParameter()

```
static int InstrumentSystems.CAS4.CAS4DLL.casSetMeasurementParameter (
    int ADevice,
```

```
int AWhat,
double AValue ) [static]
```

Sets a numeric measurement parameter for a given spectrometer

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	Integer constant, defining which parameter to modify. One of the mpid<XXX> constants starting at mpidIntegrationTime
<i>AValue</i>	Double holding the new value the measurement parameter should have

Returns

The return value is 0 if successful. Negative values indicate an error code. See [Error Handling](#)

Typical error codes include [ErrorParamOutOfRange](#) and [ErrorInvalidParameter](#).

8.1.2.72 casSetOptions()

```
static void InstrumentSystems.CAS4.CAS4DLL.casSetOptions (
    int ADevice,
    int AOptions ) [static]
```

This method sets and clears all device options for the device

Parameters

<i>ADevice</i>	The device / CASID
<i>AOptions</i>	Integer which has all bits set for the corresponding options which should be enabled and all bits cleared whose corresponding options should be disabled. Bits start with coShutter

Since some options are defined by the hardware and should not be changed manually, it is recommend to use the method [casSetOptionsOnOff](#) instead as it allows setting and clearing only the bits you actually want to modify.

8.1.2.73 casSetOptionsOnOff()

```
static void InstrumentSystems.CAS4.CAS4DLL.casSetOptionsOnOff (
    int ADevice,
    int AOptions,
    int AOnOff ) [static]
```

This method can set or clear several options for the device

Parameters

<i>ADevice</i>	The device / CASID
<i>AOptions</i>	Integer which has all bits set for the corresponding options which should modified. Bits start with coShutter
<i>AOnOff</i>	0 if the options should be disabled, all other values will enable them

Some options are defined by the device and should not be changed manually, see [casGetOptions](#).

[Delphi]

```
//turn AutoRange and AutoRange filter on
casSetOptionsOnOff(CASID, coAutorangeMeasurement or coAutorangeFilter, 1);

//turn UseTransmission off
casSetOptionsOnOff(CASID, coUseTransmission, 0);
```

8.1.2.74 casSetShutter()

```
static void InstrumentSystems.CAS4.CAS4DLL.casSetShutter (
    int ADevice,
    int OnOff ) [static]
```

Sets the shutter position of the given device.

Parameters

<i>ADevice</i>	The device / CASID
<i>OnOff</i>	Desired shutter state: either casShutterOpen or casShutterClose

If the device does not have a shutter ([coShutter](#) option missing), this method will not raise an error!

Warning

Error handling for this method using [casGetError](#) is essential!

8.1.2.75 casSetStatusLED()

```
static void InstrumentSystems.CAS4.CAS4DLL.casSetStatusLED (
    int ADevice,
    int AWhat ) [static]
```

Method to control the status LED of the spectrometer

Parameters

<i>ADevice</i>	The device / CASID
<i>AWhat</i>	One of the ext constants starting with extNoError describing the state the status LED should change to

In general, the status LED is controlled by the firmware of the spectrometer and an error of filter or shutter operation is signalled via the status-LED to the operator. The `casSetStatusLED` command can be used to override this and to generate a user defined error (e.g. if the sensor temperature rises above -8 °C). If the LED signals a malfunction of the shutter or filterwheel, a successful operation of the same component brings the LED back to the "normal" state. For example when you manually signal a shutter error using `casSetStatusLED`, a successful shutter operation causes the LED to stop blinking (the LED is green again). If you operate the filter wheel successfully this has no influence on the status LED. In general, using this method is not recommended nor necessary. Refer to the hardware manual of the spectrometer for more details about the status LED and it's possible states.

The following table list the states of the CAS 140CT

AWhat	LED State	Description
N/A	red	after power-on, not initialized
N/A	orange	spectrometer is busy, e.g. moving filter wheel
extNoError	green	initialized, no error
extFilterBlink	red, blinking	filter error
extShutterBlink	green, blinking	shutter error
extExternalError	orange, blinking	user defined error

Note

Use [casGetError](#) for error handling!

8.1.2.76 [casStart\(\)](#)

```
static int InstrumentSystems.CAS4.CAS4DLL.casStart (
    int ADevice ) [static]
```

Starts a measurement for the given device and returns immediately

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

Returns

The return value depends on the interface type. Use [casGetError](#) to retrieve error information

Unlike [casMeasure](#), this method does not wait until the measurement has been performed and it also doesn't store the measured spectrum internally. Use the method [casFIFOHasData](#) to wait until the measurement has finished and [casGetFIFOData](#) to store the spectrum internally.

Note

It is not guaranteed that all interface types actually return immediately. Some may return only after the integration time is over.

Before starting a measurement, you should check whether the [density filter](#) needs to be updated and whether a [dark current](#) measurement is necessary. Calling [casPerformActionEx](#) with `paPrepareMeasurement` is also recommended.

8.1.2.77 [casUpdateCalibrations\(\)](#)

```
static void InstrumentSystems.CAS4.CAS4DLL.casUpdateCalibrations (
    int ADevice ) [static]
```

Updates the calibration information for the given device. Deprecated! Use [paUpdateSpectralCalibration](#).

Parameters

<i>ADevice</i>	The device / CASID
----------------	--------------------

This method recalculates the spectral calibration factors after changing them via [casSetCalibrationFactors](#). If you don't edit the spectral calibration manually, there is no need to call this method, since it is called automatically whenever necessary (for example when changing the density filter or the currently active parameter set).

8.1.2.78 cmXYToDominantWavelength()

```
static int InstrumentSystems.CAS4.CAS4DLL.cmXYToDominantWavelength (
    double x,
    double y,
    double IllX,
    double IllY,
    ref double LambdaDom,
    ref double Purity ) [static]
```

Calculates dominant wavelength aka LambdaDom and purity from a given color coordinate and illuminant reference. This method is deprecated because it does not take the new [mpidColormetricType](#) into account. Use [casCalculateLambdaDom](#) instead.

Parameters

<i>x</i>	The x color coordinate of the spectrum
<i>y</i>	The y color coordinate of the spectrum
<i>IllX</i>	The x color coordinate of the illuminant reference
<i>IllY</i>	The y color coordinate of the illuminant reference
<i>LambdaDom</i>	The variable which will receive the calculated dominant wavelength
<i>Purity</i>	The variable which will receive the calculated purity

Returns

The return value is 0 and should be ignored

IllX and IllY are the color coordinates of the illuminant reference. One of the typical references is Illuminant E where x and y are 0.3333. To use the color coordinates of the previously measured spectrum, retrieve them with [casGetColorCoordinates](#) and pass them as x and y. Purity is 0 if x = IllX and y = IllY. Purity will never exceed 1.

8.1.3 Member Data Documentation

8.1.3.1 mpidACQStateLine

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidACQStateLine = 37
```

integer: number of the digital out port which should be used for the [toShowACQState](#) option in [mpidTriggerOptions](#)). The port specified by this parameter is set to the level given by [mpidACQStateLinePolarity](#) after the spectrometer has been started and is ready to accept a trigger. The level is restored after the acquisition, i.e. after integration and read-out time. For possible port values refer to [casSetDigitalOut](#) and the hardware manual of the spectrometer.

Note

When setting this parameter the level of the given line will be changed immediately if the device has been initialized and the [toShowACQState](#) trigger option is enabled.

Warning

Not all device types support changing this parameter which also causes varying default values depending on the device type. Therefore always verify this parameter, if you're relying on a specific line.

See chapter [Triggered Measurements](#) for an overview of all related options and parameters.

8.1.3.2 mpidACQStateLinePolarity

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidACQStateLinePolarity = 38
```

integer: level of the ACQ state line. Only applies if the [toShowACQState](#) option is enabled in [mpidTriggerOptions](#). The value is 0 for low level, all other values indicate high level. The port specified by [mpidACQStateLine](#) is set to this level after the spectrometer has been started and is ready to accept a trigger. The level is restored after the acquisition, i.e. after integration and read-out time.

Note

When setting this parameter the level of the digital port will be changed immediately, if the device has been initialized and the [toShowACQState](#) trigger option is enabled.

Warning

Not all device types support changing this parameter, which also causes varying default values depending on the device type. Therefore always verify this parameter, if you're relying on a specific polarity. We recommend using the low level for ACQ. You might also check the [tcoBusyStatePolarity](#) capability to see if the spectrometer supports a custom polarity.

See chapter [Triggered Measurements](#) for an overview of all related options and parameters.

8.1.3.3 mpidBusyStateLine

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidBusyStateLine = 39
```

integer: number of the digital out port which should be used for the [toShowBusyState](#) option in [mpidTriggerOptions](#). The port specified by this parameter is set to the level given by [mpidBusyStateLinePolarity](#) after the spectrometer has been started and is ready to accept a trigger. The level is restored after the trigger has been received. For possible port values refer to [casSetDigitalOut](#) and the hardware manual of the spectrometer.

Note

When setting this parameter the level of the given line will be changed immediately if the device has been initialized and the [toShowBusyState](#) trigger option is enabled.

Warning

Not all device types support changing this parameter, which also causes varying default values depending on the device type. Therefore always verify this parameter, if you're relying on a specific line. You might also check the [tcoBusyStatePolarity](#) capability to see if the spectrometer supports a custom polarity.

See chapter [Triggered Measurements](#) for an overview of all related options and parameters.

8.1.3.4 mpidBusyStateLinePolarity

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidBusyStateLinePolarity = 40
```

integer: level of the Busy state line. Only applies if the [toShowBusyState](#) option is enabled in [mpidTriggerOptions](#). The value is 0 for low level, all other values indicate high level. The port specified by [mpidACQStateLine](#) is set to this level after the spectrometer has been started and is ready to accept a trigger. The level is restored after the acquisition, i.e. after integration and read-out time.

Note

When setting this parameter the level of the digital port will be changed immediately, if the device has been initialized and the [toShowACQState](#) trigger option is enabled.

Warning

Not all device types support changing this parameter, which also causes varying default values depending on the device type. Therefore always verify this parameter, if you're relying on a specific polarity. We recommend using the low level for ACQ.

See chapter [Triggered Measurements](#) for an overview of all related options and parameters.

8.1.3.5 mpidCurrentCCDTemperature

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidCurrentCCDTemperature = 46
```

reading this mpid performs a temperature measurement and returns the CCD temperature in degrees Celsius. If the device does not support temperature measurements the return value will be NAN. To verify that the device supports temperature measurements, check the [coGetTemperature](#) flag in [casGetOptions](#).

Note

This method call might cause an error if the CCD temperature is outside of the allowed range ([ErrorCCDTemperatureFail](#)). Therefore it is especially important to check [casGetError](#) after retrieving the current CCD temperature.

8.1.3.6 mpidFlashDelayTime

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidFlashDelayTime = 29
```

float: the delay time of the flash output signal in milliseconds. Applies only when flash is activated, i.e. [mpidFlashType](#) <> [ftNone](#). If the delay time is negative, the flash signal and the delay happen before starting the measurement, otherwise the delay and the flash signal occur after the measurement started.

Note

Not all spectrometers support flash delay! Check the relevant bits in [dpidTriggerCapabilities](#) to ensure that the desired behaviour is actually supported by the spectrometer.

8.1.3.7 mpidForceFilter

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidForceFilter = 34
```

integer: flag which controls whether the filter wheel is moved for every measurement, even if it had been set to the same position previously. Internally the flag only affects [dpidNeedDensityFilterChange](#), it will always return True if the ForceFilter flag is set.

Note

Some device types might enable ForceFilter during [casInitialize](#). This is typically done for devices which can read the current filter state from the device ([coGetFilter](#) in [casGetOptions](#)). If only one application and device handle access the spectrometer, setting mpidForceFilter back to 0 can help minimize the measurement cycle time, but note that some devices might even enforce that ForceFilter is enabled. This is typically done for devices where monitoring the filter is a very fast operation (<1ms). Refer to the chapter [Density Filter](#) for an overview of all methods and parameter related to the density filter.

8.1.3.8 mpidLastCCDTemperature

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidLastCCDTemperature = 47
```

returns the previously measured CCD temperature in degrees Celsius. If the last temperature measurement is too old, a new temperature measurement will be performed. The interval which triggers a new temperature measurement is device dependent. If the device does not support temperature measurements the return value will be NAN. To verify that the device supports temperature measurements, check the [coGetTemperature](#) flag in [casGetOptions](#).

Note

If a new temperature measurement had to be performed, this method call might cause an error if the CCD temperature is outside of the allowed range ([ErrorCCDTemperatureFail](#)), so it is especially important to check [casGetError](#).

8.1.3.9 mpidRelSaturation

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidRelSaturation = 25
```

returns the relative saturation (between 0 and 100%) of the previous successful measurement for the [active parameter set](#). mpidRelSaturation can be checked against the [dpidRelSaturationMin](#) and [dpidRelSaturationMax](#) range, to ensure that the measurement has a sufficient signal level. However, to check for oversaturated measurements use [mpidMaxADCValue](#) and check it against [dpidADCRange](#). Because with averaging, mpidRelSaturation might be < 100% even though there were saturated spectra!

Note

Contrary to [mpidMaxADCValue](#), the relative saturation is also affected by the [Dark Current](#). It returns the percentage of the signal of the actual ADC range available for light, i.e. [dpidADCRange](#) minus the dark current at the [mpidMaxADCPixel](#)

8.1.3.10 mpidTOPAperture

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidTOPAperture = 30
```

integer: the TOP aperture which should be taken into account when the calibration is applied to the spectrum. The TOP aperture parameter is a 0-based index ranging from 0 to 6 corresponding to the TOP apertures 1 to 7. This parameter doesn't actually adjust the aperture of the TOP, but controls which aperture calibration factor will be applied.

Note

If no TOP is used, leave [mpidTOPAperture](#) at 0. To find out whether a specific TOP aperture is calibrated, use [casGetCalibrationFactors](#) with [gcfTOPApertureFactor](#) to check whether the factor for the aperture is not 0.

8.1.3.11 mpidTOPDistance

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidTOPDistance = 31
```

the distance in mm from the DUT to the reference plane of the TOP. This parameter doesn't actually adjust the distance of the TOP, but controls which calibration factor will be applied.

Note

To check which distance range is covered by the calibration, use [casGetCalibrationFactors](#) with [gcfTOP↔DistanceFactor](#). If the calibration does not contain TOP distance factors, this parameter has no effect.

8.1.3.12 mpidTOPFieldOfView

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidTOPFieldOfView = 44
```

returns the field of view for the current [mpidTOPAperture](#) and [mpidTOPDistance](#).

Note

Not all configurations support the calculation of field of view - in these cases [mpidTOPFieldOfView](#) will always return 0. It is recommended to only use this measurement condition when the [coTOPHasFieldOfViewConfig](#) option is returned by [casGetOptions](#).

8.1.3.13 mpidTriggerOptions

```
const int InstrumentSystems.CAS4.CAS4DLL.mpidTriggerOptions = 33
```

integer: bit-set describing current trigger options. The set consists of the various to<XXX> bits, starting with [toAcceptOnlyWhenReady](#).

Note

Some options apply to hardware and software triggers! See chapter [Triggered Measurements](#) for an overview. Additionally a device type might enforce specific trigger options, for example [toShowACQState](#) is always enabled for some CAS USB devices with hardware triggered measurements ([mpidTriggerSource](#) = [trgFlip↔Flop](#)).

The documentation for this class was generated from the following file:

- InstrumentSystems.CAS4.CAS4DLL.cs

Index

casAssignDeviceEx
InstrumentSystems::CAS4::CAS4DLL, 56

casCalculateCRI
InstrumentSystems::CAS4::CAS4DLL, 57

casCalculateCorrectedData
InstrumentSystems::CAS4::CAS4DLL, 57

casCalculateLambdaDom
InstrumentSystems::CAS4::CAS4DLL, 58

casCalculateTOPParameter
InstrumentSystems::CAS4::CAS4DLL, 58

casChangeDevice
InstrumentSystems::CAS4::CAS4DLL, 59

casClearCalibration
InstrumentSystems::CAS4::CAS4DLL, 60

casClearDarkCurrent
InstrumentSystems::CAS4::CAS4DLL, 60

casColorMetric
InstrumentSystems::CAS4::CAS4DLL, 61

casConvoluteTransmission
InstrumentSystems::CAS4::CAS4DLL, 61

casCreateDevice
InstrumentSystems::CAS4::CAS4DLL, 61

casCreateDeviceEx
InstrumentSystems::CAS4::CAS4DLL, 62

casDeleteParamSet
InstrumentSystems::CAS4::CAS4DLL, 62

casDoneDevice
InstrumentSystems::CAS4::CAS4DLL, 63

casFIFOHasData
InstrumentSystems::CAS4::CAS4DLL, 63

casGetCCT
InstrumentSystems::CAS4::CAS4DLL, 65

casGetCRI
InstrumentSystems::CAS4::CAS4DLL, 67

casGetCalibrationFactors
InstrumentSystems::CAS4::CAS4DLL, 63

casGetCentroid
InstrumentSystems::CAS4::CAS4DLL, 66

casGetColorCoordinates
InstrumentSystems::CAS4::CAS4DLL, 66

casGetDLLFileName
InstrumentSystems::CAS4::CAS4DLL, 72

casGetDLLVersionNumber
InstrumentSystems::CAS4::CAS4DLL, 73

casGetDarkCurrent
InstrumentSystems::CAS4::CAS4DLL, 67

casGetData
InstrumentSystems::CAS4::CAS4DLL, 68

casGetDeviceParameter
InstrumentSystems::CAS4::CAS4DLL, 68

casGetDeviceParameterString
InstrumentSystems::CAS4::CAS4DLL, 69

casGetDeviceTypeName
InstrumentSystems::CAS4::CAS4DLL, 69

casGetDeviceTypeOption
InstrumentSystems::CAS4::CAS4DLL, 70

casGetDeviceTypeOptionName
InstrumentSystems::CAS4::CAS4DLL, 70

casGetDeviceTypeOptions
InstrumentSystems::CAS4::CAS4DLL, 71

casGetDeviceTypes
InstrumentSystems::CAS4::CAS4DLL, 71

casGetDigitalIn
InstrumentSystems::CAS4::CAS4DLL, 71

casGetDigitalOut
InstrumentSystems::CAS4::CAS4DLL, 72

casGetError
InstrumentSystems::CAS4::CAS4DLL, 73

casGetErrorMessage
InstrumentSystems::CAS4::CAS4DLL, 73

casGetExtendedColorValues
InstrumentSystems::CAS4::CAS4DLL, 74

casGetExternalADCValue
InstrumentSystems::CAS4::CAS4DLL, 75

casGetFIFOData
InstrumentSystems::CAS4::CAS4DLL, 77

casGetFilterName
InstrumentSystems::CAS4::CAS4DLL, 77

casGetMeasurementParameter
InstrumentSystems::CAS4::CAS4DLL, 78

casGetOptions
InstrumentSystems::CAS4::CAS4DLL, 78

casGetPeak
InstrumentSystems::CAS4::CAS4DLL, 79

casGetPhotInt
InstrumentSystems::CAS4::CAS4DLL, 79

casGetRadInt
InstrumentSystems::CAS4::CAS4DLL, 79

casGetSerialNumberEx
InstrumentSystems::CAS4::CAS4DLL, 80

casGetShutter
InstrumentSystems::CAS4::CAS4DLL, 80

casGetTriStimulus
InstrumentSystems::CAS4::CAS4DLL, 81

casGetWidth
InstrumentSystems::CAS4::CAS4DLL, 81

casGetWidthEx
InstrumentSystems::CAS4::CAS4DLL, 82

- casGetXArray
 - InstrumentSystems::CAS4::CAS4DLL, [83](#)
- casInitialize
 - InstrumentSystems::CAS4::CAS4DLL, [83](#)
- casLoadTestData
 - InstrumentSystems::CAS4::CAS4DLL, [84](#)
- casMeasure
 - InstrumentSystems::CAS4::CAS4DLL, [84](#)
- casMeasureDarkCurrent
 - InstrumentSystems::CAS4::CAS4DLL, [85](#)
- casMultiTrackCopyData
 - InstrumentSystems::CAS4::CAS4DLL, [85](#)
- casMultiTrackCount
 - InstrumentSystems::CAS4::CAS4DLL, [86](#)
- casMultiTrackDone
 - InstrumentSystems::CAS4::CAS4DLL, [86](#)
- casMultiTrackInit
 - InstrumentSystems::CAS4::CAS4DLL, [87](#)
- casMultiTrackLoadData
 - InstrumentSystems::CAS4::CAS4DLL, [87](#)
- casMultiTrackSaveData
 - InstrumentSystems::CAS4::CAS4DLL, [88](#)
- casNmToPixel
 - InstrumentSystems::CAS4::CAS4DLL, [88](#)
- casPerformAction
 - InstrumentSystems::CAS4::CAS4DLL, [89](#)
- casPerformActionEx
 - InstrumentSystems::CAS4::CAS4DLL, [89](#)
- casPixelToNm
 - InstrumentSystems::CAS4::CAS4DLL, [89](#)
- casSaveCalibration
 - InstrumentSystems::CAS4::CAS4DLL, [90](#)
- casSaveSpectrum
 - InstrumentSystems::CAS4::CAS4DLL, [90](#)
- casSetCalibrationFactors
 - InstrumentSystems::CAS4::CAS4DLL, [90](#)
- casSetDeviceParameter
 - InstrumentSystems::CAS4::CAS4DLL, [92](#)
- casSetDeviceParameterString
 - InstrumentSystems::CAS4::CAS4DLL, [93](#)
- casSetDigitalOut
 - InstrumentSystems::CAS4::CAS4DLL, [93](#)
- casSetMeasurementParameter
 - InstrumentSystems::CAS4::CAS4DLL, [93](#)
- casSetOptions
 - InstrumentSystems::CAS4::CAS4DLL, [94](#)
- casSetOptionsOnOff
 - InstrumentSystems::CAS4::CAS4DLL, [94](#)
- casSetShutter
 - InstrumentSystems::CAS4::CAS4DLL, [95](#)
- casSetStatusLED
 - InstrumentSystems::CAS4::CAS4DLL, [95](#)
- casStart
 - InstrumentSystems::CAS4::CAS4DLL, [96](#)
- casUpdateCalibrations
 - InstrumentSystems::CAS4::CAS4DLL, [96](#)
- cmXYToDominantWavelength
 - InstrumentSystems::CAS4::CAS4DLL, [97](#)
- InstrumentSystems, [33](#)
- InstrumentSystems.CAS4, [33](#)
- InstrumentSystems.CAS4.CAS4DLL, [35](#)
- InstrumentSystems::CAS4::CAS4DLL
 - casAssignDeviceEx, [56](#)
 - casCalculateCRI, [57](#)
 - casCalculateCorrectedData, [57](#)
 - casCalculateLambdaDom, [58](#)
 - casCalculateTOPParameter, [58](#)
 - casChangeDevice, [59](#)
 - casClearCalibration, [60](#)
 - casClearDarkCurrent, [60](#)
 - casColorMetric, [61](#)
 - casConvoluteTransmission, [61](#)
 - casCreateDevice, [61](#)
 - casCreateDeviceEx, [62](#)
 - casDeleteParamSet, [62](#)
 - casDoneDevice, [63](#)
 - casFIFOHasData, [63](#)
 - casGetCCT, [65](#)
 - casGetCRI, [67](#)
 - casGetCalibrationFactors, [63](#)
 - casGetCentroid, [66](#)
 - casGetColorCoordinates, [66](#)
 - casGetDLLFileName, [72](#)
 - casGetDLLVersionNumber, [73](#)
 - casGetDarkCurrent, [67](#)
 - casGetData, [68](#)
 - casGetDeviceParameter, [68](#)
 - casGetDeviceParameterString, [69](#)
 - casGetDeviceTypeName, [69](#)
 - casGetDeviceTypeOption, [70](#)
 - casGetDeviceTypeOptionName, [70](#)
 - casGetDeviceTypeOptions, [71](#)
 - casGetDeviceTypes, [71](#)
 - casGetDigitalIn, [71](#)
 - casGetDigitalOut, [72](#)
 - casGetError, [73](#)
 - casGetErrorMessage, [73](#)
 - casGetExtendedColorValues, [74](#)
 - casGetExternalADCValue, [75](#)
 - casGetFIFOData, [77](#)
 - casGetFilterName, [77](#)
 - casGetMeasurementParameter, [78](#)
 - casGetOptions, [78](#)
 - casGetPeak, [79](#)
 - casGetPhotInt, [79](#)
 - casGetRadInt, [79](#)
 - casGetSerialNumberEx, [80](#)
 - casGetShutter, [80](#)
 - casGetTriStimulus, [81](#)
 - casGetWidth, [81](#)
 - casGetWidthEx, [82](#)
 - casGetXArray, [83](#)
 - casInitialize, [83](#)
 - casLoadTestData, [84](#)
 - casMeasure, [84](#)
 - casMeasureDarkCurrent, [85](#)

- casMultiTrackCopyData, [85](#)
 - casMultiTrackCount, [86](#)
 - casMultiTrackDone, [86](#)
 - casMultiTrackInit, [87](#)
 - casMultiTrackLoadData, [87](#)
 - casMultiTrackSaveData, [88](#)
 - casNmToPixel, [88](#)
 - casPerformAction, [89](#)
 - casPerformActionEx, [89](#)
 - casPixelToNm, [89](#)
 - casSaveCalibration, [90](#)
 - casSaveSpectrum, [90](#)
 - casSetCalibrationFactors, [90](#)
 - casSetDeviceParameter, [92](#)
 - casSetDeviceParameterString, [93](#)
 - casSetDigitalOut, [93](#)
 - casSetMeasurementParameter, [93](#)
 - casSetOptions, [94](#)
 - casSetOptionsOnOff, [94](#)
 - casSetShutter, [95](#)
 - casSetStatusLED, [95](#)
 - casStart, [96](#)
 - casUpdateCalibrations, [96](#)
 - cmXYToDominantWavelength, [97](#)
 - mpidACQStateLine, [97](#)
 - mpidACQStateLinePolarity, [98](#)
 - mpidBusyStateLine, [98](#)
 - mpidBusyStateLinePolarity, [98](#)
 - mpidCurrentCCDTemperature, [99](#)
 - mpidFlashDelayTime, [99](#)
 - mpidForceFilter, [99](#)
 - mpidLastCCDTemperature, [100](#)
 - mpidRelSaturation, [100](#)
 - mpidTOPAperture, [100](#)
 - mpidTOPDistance, [101](#)
 - mpidTOPFieldOfView, [101](#)
 - mpidTriggerOptions, [101](#)
 - InstrumentSystems::CAS4::CAS4DLL, [101](#)
 - mpidTOPFieldOfView
 - InstrumentSystems::CAS4::CAS4DLL, [101](#)
 - mpidTriggerOptions
 - InstrumentSystems::CAS4::CAS4DLL, [101](#)
-
- mpidACQStateLine
 - InstrumentSystems::CAS4::CAS4DLL, [97](#)
 - mpidACQStateLinePolarity
 - InstrumentSystems::CAS4::CAS4DLL, [98](#)
 - mpidBusyStateLine
 - InstrumentSystems::CAS4::CAS4DLL, [98](#)
 - mpidBusyStateLinePolarity
 - InstrumentSystems::CAS4::CAS4DLL, [98](#)
 - mpidCurrentCCDTemperature
 - InstrumentSystems::CAS4::CAS4DLL, [99](#)
 - mpidFlashDelayTime
 - InstrumentSystems::CAS4::CAS4DLL, [99](#)
 - mpidForceFilter
 - InstrumentSystems::CAS4::CAS4DLL, [99](#)
 - mpidLastCCDTemperature
 - InstrumentSystems::CAS4::CAS4DLL, [100](#)
 - mpidRelSaturation
 - InstrumentSystems::CAS4::CAS4DLL, [100](#)
 - mpidTOPAperture
 - InstrumentSystems::CAS4::CAS4DLL, [100](#)
 - mpidTOPDistance
 - InstrumentSystems::CAS4::CAS4DLL, [101](#)